

# MBA

# 2026.4

## Dental Handpiece Replacement Parts

Make	Pages	Make	Pages
ABU®	3	KMD	115
ACETON SATELEC®	4	Kool Torque®	116
ADEC®	5	Lares®	117-118
ALPHA AIR®	6	Marathon®	119-120
B.A. International®	7-11	Micro Mega®	121
BANNER®	12	Midwest®	122-125
BEAVER®	13	Miscellaneous	126
Bien Air®	14-22	MK-Dent®	127-136
Borden®	23	Morita®	137
Brasseler®	24	MTC	138
Bravo Bedetrice®	25	MTI	139
Buffalo®	26	Muss®	140
C. De Giorgi®	27	NSK®	141-166
Castellini®	28-29	Ritter®	167
Chirana®	30	Sable®	168-169
COXO®	31	Saeshin®	170-171
Dabi Atlante®	32-33	Schick®	172-173
Densim®	34-36	Siemens®	174-182
DME	37	Sirona®	183-190
Dynadent®	38	SS White®	191
Dynasurge®	39	Star®	192-194
Elco® - Austria	40	Starflite®	195
Falcon®	41	Sternweber®	196
FangSCO®	42-43	TKD	197
Faro®	44-45	Torx®	198
Futura®	46	W&H®	199-219
Gamborini®	47	Yoshida®	220
Henry Schein®	48-49		
Impact-Air®	50		
Jalenco®	51		
KaVo®	52-114		



To find and order parts online, simply type the part number from this catalogue into the search box on any web page.

Unless specified otherwise, all parts are generic offerings to replace the original parts. Some items are offered with alternative materials. Please check carefully before placing an order.

## 350000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SFR144A-D27-ZSO-KC**



**Replacement Bearing**

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SFR144A-D27-ZSO-TO**



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

# ACETON SATELEC®

## F12200 Airscalers

Position / Description / Application

Manufacturers Part

MBA Part

**Bush / Rotor Repair Kit**

**707091**

**BKU67091**

Brand: MBA



Key



Usually stocked or Short Delivery

\*\* Available while stocks last



Long Delivery or Large Minimum Quantities

*Not in current production*

## 5000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	

## 5100

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	





Key

Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# B.A. International®


## 13Z

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece  Brand: MBA		<b>BKU65587</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece  Brand: MBA		<b>BKU65687</b>	
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece  Brand: MBA		<b>BKU65687X</b>	
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece  Brand: MBA		<b>BKU65587X</b>	

## BA 200 LTS

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Dental replacement drive kit  Brand: MBA	<b>6175587 + 6175603</b>	<b>BKU68703</b>	

## BA 40 LSS Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Dental replacement drive kit  Brand: MBA	<b>6405695 + 6000959</b>	<b>BKU69559</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# B.A. International®

## BA200LTS

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete

**BKU65687**

Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece

Brand: MBA



### Head Gear Complete

**BKU65587X**

Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: MBA



### Head Gear Complete

**BKU65587**

Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece

Brand: MBA



### Head Gear Complete

**BKU65687X**

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: MBA



### Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: myonic

**BKU65687M**



## BA40

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**C875060 + C880022**

**BKU66022**

Dental replacement drive kit

Brand: MBA



## BA45

Position / Description / Application

Manufacturers Part

MBA Part


### Drive Kit

**C875060 + C880022**

**BKU66022**

Dental replacement drive kit

Brand: MBA



## BA45 LS

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**6000983 + 6000959**

**BKU68359**

Dental replacement drive kit

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# B.A. International®

## Diameter 8.7 BA200LTS

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: myonic

**BKU65687M**



## Optima® BA525

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine Excluding Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA

**BKUR310000**




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

**BKURC310000**




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

**BKUR3100004**




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA

**BKUR3100005**




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

**BKUR3100003**



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# B.A. International®

## Optima® BA528

Position / Description / Application

Manufacturers Part


MBA Part

### Replacement Turbine Excluding Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA

BKUR310001



### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

BKUR3100013



### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA

BKUR3100015



### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

BKUR3100014




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic

BKURC310001



## Ultima® 650L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

High Speed Bearing - Phenolic Retainer  
Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer  
2 Required  
Brand: myonic

SR144-D106-ZS-KC



### Bearing for Cap End - Alternative

High Speed Ceramic-Hybrid Bearing - Torlon Retainer  
Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38)  
2 Required  
Brand: MBA

SR144-D106-CSN-ZS-TO



### Bearing for Cap End - Alternative

High Speed Bearing - Torlon Retainer  
Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer  
2 Required  
Brand: myonic

SR144-D106-ZS-TO



### Bearing for Cap End - Alternative

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer  
Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38)  
2 Required  
Brand: MBA

SR144-D106-CSN-ZS-KC



Key

Usually stocked or Short Delivery





\*\* Available while stocks last

Long Delivery or Large Minimum Quantities





Not in current production

# B.A. International®

## Ultimate 635

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200243</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200245</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200244</b>	

## Ultimate 670

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220022</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200225</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200224</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200223</b>	

Key

Usually stocked or Short Delivery





\*\* Available while stocks last

Long Delivery or Large Minimum Quantities





*Not in current production*

# B.A. International®

## Ultimate 675

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220023</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200234</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200235</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200233</b>	

## Ultimate 680

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220023</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200234</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200235</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200233</b>	

Key

Usually stocked or Short Delivery





\*\* Available while stocks last

Long Delivery or Large Minimum Quantities





*Not in current production*

# B.A. International®

## Ultimate 685

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220023</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200233</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200234</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200235</b>	

## Ultimate 755

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR220024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR2200245</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200244</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR2200243</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

# B.A. International®

## Ultimate 758

Position / Description / Application

Manufacturers Part


MBA Part

### Replacement Turbine Excluding Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

**BKUR220021**

Brand: MBA




### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

**BKUR2200213**

Brand: myonic



### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

**BKUR2200215**

Brand: MBA



### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

**BKUR2200214**

Brand: myonic



## Ultimate BA 250LT Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**6414457 + 6175603**

**BKU65703**

Brand: MBA



## Ultimate BA688 New Type

Position / Description / Application

Manufacturers Part


MBA Part

### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

**BKUR2200363**

Brand: myonic



## Ultimate BA755 New Type

Position / Description / Application

Manufacturers Part


MBA Part

### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

**BKUR2200353**

Brand: myonic





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# B.A. International®

## Ultimate Power

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>6175587 + 6175603</b>	<b>BKU68703</b>	
Dental replacement drive kit			
Brand: MBA			
<b>Drive Kit</b>	<b>6405695 + 6000959</b>	<b>BKU69559</b>	
Dental replacement drive kit			
Brand: MBA			

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Banner

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144A-D27-ZS-TO

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic



Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Slow Speed 40000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x7x3)

2 Required

S683A-ZZ-MC3-JP6



Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# BIEN AIR®

## Black Pearl®

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Thrust Washer</b> Stainless Steel  Brand: MBA	009.34.58	BKU63458	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Press Button</b>  Brand: MBA	009.40.25+009.40.30+009.40.24+009.40.23	BKU64025	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b>  Turbine sub assembly excluding bearings  Brand: MBA	1500164	BKUR240036	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1500164	BKUR2400365	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1500164	BKUR2400363	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Turbine Sub Assembly With Bearings</b>  With Full stainless bearings, Deep Groove, Phenolic retainer, myonic  Brand: myonic	1500164	BKUR2400364	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental ceramic-hybrid bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-CSN-ZS-KC	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Black Pearl® Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement O-Ring**

**705.01.81**

**OR-00804-100-N70**

High Temperature

Deleted Item. Available while stocks last.



## Black Pearl® Friction Grip

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**701.01.00**

**SR2-5A-D17-ZS-TO**

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Bora®

Position / Description / Application	Manufacturers Part	MBA Part	
<b>LED Bulb</b>  Brand: MBA	<b>204.33.07</b>	<b>BKU63307LED</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>009.34.46+209.34.11+009.40.24+009.34.58</b>	<b>BKU63446</b>	
<b>Thrust Washer</b> Stainless Steel  Brand: MBA	<b>009.34.58</b>	<b>BKU63458</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1500181</b>	<b>BKUR240037</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1500181</b>	<b>BKUR2400375</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1500181</b>	<b>BKUR2400373</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1500181</b>	<b>BKURMU240037</b>	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Phenolic retainer, myonic  Brand: myonic	<b>1500181</b>	<b>BKUR2400374</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>1500181</b>	<b>BKURM240037</b>	

Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# BIEN AIR®

<p><b>Xenon Lamp</b></p> <p>Brand: MBA</p>	<p>204.33.07 (1500000)</p>	<p>BKU63307Xenon</p>	 <p>LED Bulb, Couplers &amp; Dental Motors</p>
<p><b>Bearing for Rotor Cartridge</b></p> <p>High Speed Bearing - Phenolic Retainer</p> <p>Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)</p>	<p>701.01.00</p>	<p>SR2-5A-D17-ZS-KC</p>	
<p><b>Bearing for Rotor Cartridge</b></p> <p>High Speed Bearing - Torlon Retainer</p> <p>Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)</p> <p>Brand: MBA</p>	<p>701.01.00</p>	<p>SR2-5A-D17-ZS-TO</p>	
<p><b>Bearing for Rotor Cartridge</b></p> <p>High Speed Ceramic-Hybrid Bearing - Phenolic Retainer</p> <p>Single shielded Abec 7 precision special dental ceramic-hybrid bearing with high speed phenolic retainer (3.18x7.94x3.57)</p> <p>Brand: MBA</p>	<p>701.01.00</p>	<p>SR2-5A-D17-CSN-ZS-KC</p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Bearing - Phenolic Retainer</p> <p>Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)</p> <p>2 Required</p>	<p>701.01.00</p>	<p>SR2-5A-D17-ZS-KC</p>	

## Bora® Air Driven Handpiece





Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Replacement O-Ring</b></p> <p>High Temperature</p> <p>Deleted Item. Available while stocks last.</p>	<p>705.01.81</p>	<p>OR-00804-100-N70</p>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
 Not in current production

# BIEN AIR®


## Bora® lina TU

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1500164	BKUR240036	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1500164	BKUR2400365	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1500164	BKUR2400363	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Phenolic retainer, myonic  Brand: myonic	1500164	BKUR2400364	

## CA 1:1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b>  Brand: MBA	1500289 / 1501083	BKU60289	
<b>Head Gear Complete 1:1</b>  Brand: MBA	1500549 / 1501085	BKU61085	
<b>Replacement Press Button</b>  Brand: MBA	1500550	BKU60550	

## CA 1:1 L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	1501085 + 1501083	BKU68583	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## CA 1:1 NOVA

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 1:1**

**1503014**

**BKU63014**

Brand: MBA



## CA 1:5

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**701.01.75**

**SMR4.762-2.35A-ZZS-MC3-GR7-TO**

High Speed Bearing - Torlon Retainer

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)

Brand: MBA



**Head Gear Complete - Premium**

**1501933**

**BKU61933**

For Handpiece with Wedge Type Chuck

Premium head gear assembled with ceramic ball bearings and tested in master electric handpiece




**Head Gear Complete - Standard**

**1501933**

**BKU61933X**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings and tested in master electric handpiece




**Replacement Middle Gear Complete 1:5**

**1500057**

**BKU60057**

Brand: MBA



**Replacement Press Button**

**1500169**

**BKU60169**

Brand: MBA



**Bearing for Headgear**

**701.01.68**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



## CA 1:5 L

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**1501933 + 1500057**

**BKU63357**

Brand: MBA






Key	Usually stocked or Short Delivery
	** Available while stocks last


Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# BIEN AIR®





## CA 20:1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 20.1</b>  Brand: MBA	<b>1501026</b>	<b>BKU61026</b>	
<b>Head Gear Complete 20.1</b>  Brand: MBA	<b>1501898</b>	<b>BKU61898</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>1501027</b>	<b>BKU61027</b>	

## CA 20:1 L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>1501898 + 1501026</b>	<b>BKU69826</b>	

## Croma M

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>209.30.19</b>	<b>BKUR240039</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>209.30.19</b>	<b>BKUR2400393</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>209.30.19</b>	<b>BKUR2400394</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>209.30.19</b>	<b>BKUR2400395</b>	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# BIEN AIR®

## HP11

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

701.01.65

S623-ZZ-MC3-JP6

High Precision Bearing Shielded Option  
Double shielded P6 precision stainless steel bearing (3x10x4)



### Bearing for Shaft

701.01.53

SMR85C-ZZS-MC34-TH5-L23

High Precision Bearing Standard Cage Option  
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



Brand: MBA

## LaBora® tory

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

701.01.00

SR2-5A-D17-ZS-KC

High Speed Bearing - Phenolic Retainer  
Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)  
2 Required



## MC2

Position / Description / Application

Manufacturers Part


MBA Part

### LED Bulb

211.60.16

BKU66016LED

Brand: MBA




### Replacement Motor Brush

211.45.06

BKU64506

Brand: MBA




### Xenon Lamp

211.60.16 (1500007)

BKU66016Xenon

Brand: MBA



## MC2 Intra-Matic

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement O-Ring

705.01.74

OR-00800-070-N70

High Temperature  
Deleted Item. Available while stocks last. Special o-ring for dental applications  
2 Required



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## MC2 Isolite 300

Position / Description / Application

Manufacturers Part

MBA Part


**Dental Motor Set with 15 parts**

**1502219-001**

**BKU62219**

Includes parts - 3x 1300967 - 1x 705.03.15 - 1x 652.02.02 - 1x 701-01-14 - 1x 701.01.23 - 2x 705.02.07 - 2x 211.60.13 - 2x 211.45.06 - 1x 011.35.37 - 1x 1500391

Brand: MBA



## MC2 Isolite100

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**701.01.23**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



## MC2 Isolite300

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**701.01.23**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities


*Not in current production*

# BIEN AIR®


## MC3

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)	701.01.70	SR1-5A-ZZS-MC3-JP6	
<b>Dental Motor Set with 21 parts</b> Includes parts - 3x 01862300 - 1x 04062500 - 1x 04061800 - 1x 04062300 - 1x 04058800 - 2x 04062100 - 2x 04062200 - 2x 04062400 - 2x 04058800 - 3x 04059200 - 2x 04059300 - 1x 04069600 Brand: MBA	1502217-001	BKU62217	
<b>LED Bulb</b> Brand: MBA	211.60.16	BKU66016LED	
<b>Replacement Motor Brush</b> Brand: MBA	211.75.08	BKU67508	
<b>Xenon Lamp</b> Brand: MBA	211.60.16 (1500007)	BKU66016Xenon	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	701.01.71	SMR104A-ZZ-MC3-JP6	

## Micromotor Aquilon 830

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front and Brush End Bearings</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.14	SR2-5A-D17-ZS-KC	

## Micromotor Aquilon 850

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front and Brush End Bearings</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.14	SR2-5A-D17-ZS-KC	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Micromotor Isolite 100 and 300

Position / Description / Application

Manufacturers Part

MBA Part

### Front Bearing

701.01.14

SR2-5A-D17-ZS-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)



## Micromotor MC2 20

Position / Description / Application

Manufacturers Part

MBA Part

### Front Bearing

701.01.14

SR2-5A-D17-ZS-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)



### Replacement O-Ring Small

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

OR-00386-091-N70



### Replacement O-Ring

Medium

Special sized o-ring used in dental Applications

OR-00400-100-N70



### Replacement O-Ring Large

High Temperature

Deleted Item. Available while stocks last.

OR-00804-100-N70



Key





Usually stocked or Short Delivery

\*\* Available while stocks last





Long Delivery or Large Minimum Quantities

Not in current production

## Micromotor MC2 20GT

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
<b>Replacement O-Ring Small</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
<b>Replacement O-Ring</b> Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
<b>Replacement O-Ring Large</b> High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	





## Micromotor MC2 20GTAV

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
<b>Replacement O-Ring Small</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
<b>Replacement O-Ring</b> Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
<b>Replacement O-Ring Large</b> High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	





Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Micromotor MC2 40E

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
<b>Replacement O-Ring Small</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
<b>Replacement O-Ring</b> Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
<b>Replacement O-Ring Large</b> High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	





## Micromotor MC2 40GT

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
<b>Replacement O-Ring Small</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
<b>Replacement O-Ring</b> Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
<b>Replacement O-Ring Large</b> High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	


Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*


## Micromotor MC2 40GTAV

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
<b>Replacement O-Ring Small</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
<b>Replacement O-Ring</b> Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
<b>Replacement O-Ring Large</b> High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	


## Motor Isolite 300

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	

## Motor Isolite100

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	



## Motor MC2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## MX

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3.97x7.94x3.18)	701.01.94	SR155A-ZZS-MC3-JP6	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4.76x7.94x3.18)	701.01.93	SR156A-ZZS-MC3-JP6	



## Ondine

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	209.30.19	BKUR240039	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	209.30.19	BKUR2400395	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	209.30.19	BKUR2400394	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	209.30.19	BKUR2400393	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required Brand: Myonic		SR144A-D27-ZS-TO	


Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*


## PA11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	701.01.65	S623-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	701.01.53	SMR85C-ZZS-MC34-TH5-L23	

## PM AM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)	701.01.77	SR133A-ZZS-MC3-JP6	

## PM RM 1:1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)	701.01.77	SR133A-ZZS-MC3-JP6	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# BIEN AIR®

## Prestige





Position / Description / Application	Manufacturers Part	MBA Part	
<b>LED Bulb</b>  Brand: MBA	204.33.07	BKU63307LED	
<b>Replacement Press Button</b>  Brand: MBA	009.30.82+209.31.34+009.30.63+009.31.45	BKU63082	
<b>Xenon Lamp</b>  Brand: MBA	204.33.07 (1500000)	BKU63307Xenon	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required Brand: Myonic		SR144A-D27-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)	SG1329X-D	SR144A-D31-CSN-WZ-TO	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**




Long Delivery or Large Minimum Quantities  
*Not in current production*

# BIEN AIR®

## Prestige S32L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1500485	BKUR240040	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1500485	BKUR2400404	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1500485	BKUR2400405	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1500485	BKUR2400403	




## TD783

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Rotor / Cartridge</b>  Brand: MBA	208.55.01	BKUR2400411	
<b>Replacement Rotor / Cartridge</b>  Brand: MBA	208.55.01	BKUR240041	
<b>Replacement Rotor / Cartridge</b>  Brand: MBA	208.55.01	BKUR2400412	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TDS890

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Rotor / Cartridge</b>	<b>208.55.01</b>	<b>BKUR240041</b>	
Brand: MBA			
<b>Replacement Rotor / Cartridge</b>	<b>208.55.01</b>	<b>BKUR240042</b>	
Brand: MBA			
<b>Replacement Rotor / Cartridge</b>	<b>208.55.01</b>	<b>BKUR240041</b>	
Brand: MBA			

## Tornado

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	<b>1501968</b>	<b>BKUR240042</b>	
Brand: MBA			
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	<b>1501968</b>	<b>BKUR2400423</b>	
Brand: myonic			
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	<b>1501968</b>	<b>BKUR2400424</b>	
Brand: myonic			
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	<b>1501968</b>	<b>BKUR2400425</b>	
Brand: MBA			
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	<b>SG1302X-D</b>	<b>SR144A-D27-ZS-TO</b>	
Brand: Myonic			

Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Unifix 4V Coupler

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement O-Ring**

**1300380+1300391**

**BKU60391**

Special O-Ring for Dental Applications

Special dental applications O-Ring

Brand: MBA



Key



Usually stocked or Short Delivery

\*\* Available while stocks last



Long Delivery or Large Minimum Quantities

*Not in current production*



## 60

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-KC</b>	

## 70

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-KC</b>	





## Airotor Silent Cartridge

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-KC</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## Silencer

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		SR144-D106-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		SR144-D106-CSN-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		SR144-D106-ZS-TO	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		SR144-D106-CSN-ZS-TO	 <b>Image Coming Soon</b>

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Brio 440 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA		<b>BKUR310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR310003</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA		<b>BKUR310005</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKURC310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic		<b>BKUR310004</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Bravo Bedetrice®

## Drill 8 12000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**S685A-ZZ-MC3-JP6**

Stainless Bearing

Double shielded P6 precision stainless steel bearing (5x11x5)

2 Required



Key

Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## EM3000 Nose Cone

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing Motor End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (4x13x5)

S624-ZZ-MC3-JP6



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Goldspeed EVO M5 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

97230612 + 97230601

BKU61201

Brand: MBA



## Goldspeed EVO R20 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

97230522 + 97230527

BKU62227

Brand: MBA



## Goldspeed EVO S1 L

Position / Description / Application

Manufacturers Part


MBA Part

Drive Kit

97230630 + 97230623

BKU63023

Brand: MBA



## Goldspeed S1 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

FR400N103 + FR400N104

BKU60304





Brand: MBA







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## High Power M

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	L0000874	BKUR270103	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0000874	BKUR2701034	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	L0000874	BKUR2701035	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0000874	BKUR2701033	

## Silent Power

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	L0001857	BKUR270105	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0001857	BKUR2701054	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	L0001857	BKUR2701055	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0001857	BKUR2701053	

Key

Usually stocked or Short Delivery





\*\* Available while stocks last

Long Delivery or Large Minimum Quantities





*Not in current production*

# Castellini®

## Silent Power Evo

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	95529230	BKUR270107	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	95529230	BKUR2701073	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	95529230	BKUR2701074	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	95529230	BKUR2701075	

## Titanium Gold 2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	L0001370	BKUR270104	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0001370	BKUR2701043	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	L0001370	BKUR2701044	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	L0001370	BKUR2701045	

Key

**Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TG 656

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	T031 2560 4/1	BKUR320000	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	T031 2560 4/1	BKURC3200013	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	T031 2560 4/1	BKUR3200004	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	T031 2560 4/1	BKUR3200005	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	T031 2560 4/1	BKUR3200003	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	T031 2560 4/1	BKURC3200015	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	T031 2560 4/1	BKURC320001	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## CX207 Head

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**210-G SP**

**BKUR350001**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



## CX207 Standard Head

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**210-G SP**

**BKUR350001**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



## CX207 Torque Head

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**210-G TP**

**BKUR350002**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

# Dabi Atlante®

## L350

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		SR144-D106-ZS-KC	 Image Coming Soon
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		SR144-D106-CSN-ZS-KC	 Image Coming Soon
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		SR144-D106-CSN-ZS-TO	 Image Coming Soon
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		SR144-D106-ZS-TO	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last



Long Delivery or Large Minimum Quantities  
*Not in current production*

# Dabi Atlante®

## MS350

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required Brand: myonid		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonid		<b>SR144-D106-ZS-KC</b>	
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonid		<b>SR144-D106-ZS-TO</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-KC</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-TO</b>	



## RS350

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonid		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonid		<b>SFR144A-D27-ZSO-TO</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Dublin DSLCH11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete Dental Head 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>DSSP2005</b>	<b>BKU62005</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Intermediate Piece Dental Head 1:1</b>  Brand: MBA	<b>DSSP2024</b>	<b>BKU62024</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>

## Monaco Lux

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKURC3301083</b>	 <p>Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKUR3301083</b>	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKURC3301085</b>	 <p>Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKUR3301085</b>	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>DSTC2034</b>	<b>BKUR330108</b>	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Monaco Lux Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	DSTC2035	BKUR330109	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	DSTC2035	BKURC3301093	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	DSTC2035	BKUR3301095	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	DSTC2035	BKUR3301093	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	DSTC2035	BKURC3301095	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	DSTC2035	BKURC330109	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Monaco Lux Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKUR3301083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKUR3301085</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKURC3301083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKURC3301085</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck   Brand: MBA	<b>DSTC2034</b>	<b>BKURC330108</b>	




Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Montreal Lux

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2041</b>	<b>BKURC3301105</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2041</b>	<b>BKURC3301103</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2041</b>	<b>BKUR3301105</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2041</b>	<b>BKUR3301103</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>DSTC2041</b>	<b>BKUR330110</b>	

## Montreal Lux Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2042</b>	<b>BKURC3301203</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2042</b>	<b>BKURC3301205</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>DSTC2042</b>	<b>BKURC330120</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Montreal Lux Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>DSTC2034</b>	<b>BKUR330108</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>DSTC2041</b>	<b>BKUR330110</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2041</b>	<b>BKURC3301103</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2041</b>	<b>BKUR3301105</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2041</b>	<b>BKUR3301103</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2041</b>	<b>BKURC3301105</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>DSTC2041</b>	<b>BKURC330110</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Vienna



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>DSTC2034</b>	<b>BKUR330108</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKUR3301085</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>DSTC2034</b>	<b>BKURC3301085</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKURC3301083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>DSTC2034</b>	<b>BKUR3301083</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>DSTC2034</b>	<b>BKURC330108</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last



Long Delivery or Large Minimum Quantities  
*Not in current production*

# DME

## Aden

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	



## True Speed

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Hi-Torq Surgical Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Bearing for Cap End</b></p> <p>Stainless Flanged High Speed Bearing</p> <p>Single shielded flanged Abec 7 precision stainless steel bearing with high speed phenolic retainer (4.76x9.52x3.18)</p> <p>Brand: Myonic</p>		SFR166-ZS-MC5-GR7-KB-J260	
<p><b>Bearing for Chuck End</b></p> <p>Stainless Flanged High Speed Bearing</p> <p>Single shielded flanged Abec 7 precision stainless steel bearing with high speed phenolic retainer (4.76x9.52x3.18)</p> <p>Brand: Myonic</p>		SFR166-ZS-MC5-GR7-KB-J260	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Elco® - Austria

## Elco - Austria

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**BB 7-14-5**

**S687A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)

2 Required

Brand: MBA





Key

Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 11 - 4000 RPM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Precision Bearing Deleted Item. Available while stocks last. (6x10x3) 2 Required Brand: Myonic		SMR106A-ZZS-MC34-GR7-R-G48	
<b>Replacement Bearing</b> High Precision Bearing Deleted Item. Available while stocks last. (6x10x3) 2 Required Brand: Myonic		SMR106A-ZZS-MC34-P4P-R-AF12	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## A5 Small Head

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>CCHPZ4</b>	<b>BKUR370001</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPZ4</b>	<b>BKUR3700014</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPZ4</b>	<b>BKUR3700015</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPZ4</b>	<b>BKUR3700013</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPZ4</b>	<b>BKURC3700015</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPZ4</b>	<b>BKURC3700013</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>CCHPZ4</b>	<b>BKURC370001</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## A6 Head

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKURC3700023</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPC4</b>	<b>BKURC3700025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKUR3700024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPC4</b>	<b>BKUR3700025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKUR3700023</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>CCHPC4</b>	<b>BKUR370002</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## A6 Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>CCHPC4</b>	<b>BKUR370002</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPC4</b>	<b>BKURC3700025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKUR3700024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKUR3700023</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>CCHPC4</b>	<b>BKUR3700025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>CCHPC4</b>	<b>BKURC3700023</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>CCHPC4</b>	<b>BKURC370002</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## BTC77

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

Brand: Myonic

SFR144A-D27-ZSO-TO



### Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

Brand: Myonic

SFR144A-D27-ZSO-KC



## Delta

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

SR144A-D27-ZS-TO



## LC81

Position / Description / Application

Manufacturers Part

MBA Part

### Bearings - Small 3.175-9.525-3.969mm

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3.18x9.52x3.97)

2 Required

SR2-ZZ-MC3-JP6



### Medium Bearing Replacement

BB 7.938-12.7-3.969

SR1810-ZZS-MC3

Stainless Bearing

Double shielded standard stainless steel bearing (7.94x12.7x3.97)



### Large Bearing Replacement 8-16-6mm

BB 8-16-6

S688B-ZZ-MC3-GR3-T9H-G36

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



## MC01R

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: myonic

BKU64040M




### Head Gear Complete - Standard

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA

BKU64040X



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

Not in current production

## MC01RL

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA

**BKU64040**



### Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: myonic

**BKU64040M**



### Head Gear Complete - Standard

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA

**BKU64040X**



## S380

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required

**SR144-D28-ZS-KC**



## S400

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SR144A-D27-ZS-TO**



## S405

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SR144A-D27-ZS-TO**



Key

**Usually stocked or Short Delivery**

**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities

*Not in current production*

## Star

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

**SR2-5A-D17-ZS-TO**

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



Key




Usually stocked or Short Delivery

\*\* Available while stocks last




Long Delivery or Large Minimum Quantities

*Not in current production*

## MKQ

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Armature Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		<b>S694-ZZ-MC3-JP6</b>	
<b>Bearing for Chuck End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		<b>S695-ZZ-MC3-JP6</b>	
<b>Bearing for Cap End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		<b>S687A-ZZ-MC3-JP6</b>	




## MKS

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Armature Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		<b>S694-ZZ-MC3-JP6</b>	
<b>Bearing for Chuck End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		<b>S695-ZZ-MC3-JP6</b>	
<b>Bearing for Cap End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		<b>S687A-ZZ-MC3-JP6</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## MKT

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Armature Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		<b>S694-ZZ-MC3-JP6</b>	
<b>Bearing for Chuck End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		<b>S695-ZZ-MC3-JP6</b>	
<b>Bearing for Cap End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		<b>S687A-ZZ-MC3-JP6</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 35 II

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Armature**

**101-26**

**BKU61103**

Brand: MBA



## 35 Z

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**101-26**

**BKU61103**

Brand: MBA



## L35 Z

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**101-26**

**BKU61103**

Brand: MBA



## Maxima® ELITE

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine Excluding Bearings**

**570-0136**

**BKUR310001**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine With Bearings**

**570-0136**

**BKURC310001**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**570-0136**

**BKUR3100014**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**570-0136**

**BKUR3100015**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




**Replacement Turbine With Bearings**

**570-0136**

**BKUR3100013**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# Henry Schein®

## Maxima® ELITE L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	570-0136	BKUR310001	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	570-0136	BKUR3100015	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0136	BKUR3100013	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0136	BKURC310001	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0136	BKUR3100014	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Henry Schein®

## Maxima® ELITE Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	570-0137	BKUR310000	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKUR3100003	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKURC310000	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKUR3100004	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	570-0137	BKUR3100005	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
*Not in current production*

# Henry Schein®

## Maxima® ELITE Mini L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	570-0137	BKUR310000	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKUR3100003	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKUR3100004	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	570-0137	BKUR3100005	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	570-0137	BKURC310000	

## Maxima® PRO 45

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	112-6773	BKURC210022	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	112-6773	BKURC2100223	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	112-6773	BKURC2100225	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Henry Schein®

## Maxima® PRO 45L

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine Excluding Bearings

112-6773

BKURC210022

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

112-6773

BKURC2100225

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




### Replacement Turbine With Bearings

112-6773

BKURC2100223

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



Key

Usually stocked or Short Delivery




\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

# Impact-Air®

## High Speed Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Lab Handpiece Air

Position / Description / Application

Manufacturers Part

MBA Part

### Both Turbine Bearings

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 0768LHC

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0139**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x10x4)



## 10B

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0126**

**S685A-ZZS-MC34-TH7-R-G48**

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)

Brand: MBA



## 10C

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**220-0219 / 0126**

**S685A-ZZ-MC3-JP6**

Stainless Bearing

Double shielded P6 precision stainless steel bearing (5x11x5)

2 Required



## 10CHC

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**0.524.2942**

**BKU62942**

True running accuracy less than 0.03, bore diameter 1.6mm

Brand: MBA



## 10CN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Cap End**

**0.220.1302**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



Key




Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 10E

Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Bearing for Shaft</b></p> <p>High Precision Bearing Shielded Option</p> <p>Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer. Band around OD is a manufacturers precision marker. (4x9x4)</p> <p>Brand: Myonic</p>	<b>0.220.1301</b>	<b>S684A-ZZS-MC34-GR7-R-G48</b>	
<p><b>Bearing for Shaft</b></p> <p>High Precision Bearing Standard Cage Option</p> <p>Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)</p> <p>Brand: MBA</p>	<b>0.220.1302</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	
<p><b>Replacement Chuck</b></p> <p>True running accuracy less than 0.03, bore diameter 1.6mm</p> <p>Brand: MBA</p>	<b>0.524.2942</b>	<b>BKU62942</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 10LH

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer. Band around OD is a manufacturers precision marker. (4x9x4)  Brand: Myonic	<b>0.220.1301</b>	<b>S684A-ZZS-MC34-GR7-R-G48</b>	
<b>Bearing for Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>0.220.1302</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	
<b>Replacement Chuck</b>  True running accuracy less than 0.03, bore diameter 1.6mm  Brand: MBA	<b>0.524.2942</b>	<b>BKU62942</b>	
<b>Replacement Spindle</b>  Brand: MBA	<b>0.524.1421</b>	<b>BKU61421</b>	 <p style="font-size: 8px; color: red; position: absolute; top: 5px; left: 5px;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Pin</b>  Brand: MBA	<b>0.524.3262</b>	<b>BKU63262</b>	 <p style="font-size: 8px; color: red; position: absolute; top: 5px; left: 5px;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Roll Pin</b>  Brand: MBA	<b>0.524.2742</b>	<b>BKU62742</b>	 <p style="font-size: 8px; color: red; position: absolute; top: 5px; left: 5px;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production


## 10LP

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer. Band around OD is a manufacturers precision marker. (4x9x4)  Brand: Myonic	<b>0.220.1301</b>	<b>S684A-ZZS-MC34-GR7-R-G48</b>	
<b>Bearing for Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>0.220.1302</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	
<b>Replacement Chuck</b>  True running accuracy less than 0.03, bore diameter 1.6mm  Brand: MBA	<b>0.524.2942</b>	<b>BKU62942</b>	
<b>Replacement Spindle</b>  Brand: MBA	<b>0.524.1421</b>	<b>BKU61421</b>	
<b>Replacement Pin</b>  Brand: MBA	<b>0.524.3262</b>	<b>BKU63262</b>	
<b>Replacement Roll Pin</b>  Brand: MBA	<b>0.524.2742</b>	<b>BKU62742</b>	

## 11B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)  Brand: MBA	<b>0.220.0126</b>	<b>S685A-ZZS-MC34-TH7-R-G48</b>	

## 121

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)  Brand: MBA	<b>0.220.0021</b>	<b>S683-GR7-Y</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 121A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0021**

**S683-GR7-Y**

High Precision Bearing  
Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)

Brand: MBA



## 13A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0103**

**S684A-ZZ-MC3-JP6**

High Precision Bearing  
Double shielded P6 precision stainless steel bearing (4x9x4)



**Bearing for Centre**

**220-0024 / 0103 / 1301**

**S684A-ZZ-MC3-JP6**

4x9x4 Stainless Steel Bearing  
Double shielded P6 precision stainless steel bearing (4x9x4)



**Bearing for Cap End**

**220-0068**

**S624-MC3-JP6-AF12**

Stainless Bearing  
Deleted Item. Available while stocks last. (4x13x5)



**Bearing for Chuck End**

**220-0106**

**SMR85-MC3-JP6-AF12**

Stainless Bearing  
Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



## 13B

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0139**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option  
Double shielded P6 precision stainless steel bearing (4x10x4)




**Bearing for Shaft**

**0.220.0126**

**S685A-ZZS-MC34-TH7-R-G48**

High Precision Bearing Shielded Option  
Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

## 14E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.2004**

**S683B-ZZS-MC3-GR3**

High Precision Bearing Shielded Option

Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)

Brand: MBA



**Bearing for Shaft**


**0.220.0098**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



**Bearing for Headgear**

**0.220.2005**

**SMR74-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



## 14ES

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.2004**

**S683B-ZZS-MC3-GR3**

High Precision Bearing

Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)

Brand: MBA



## 15CN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.1306**

**SMR9525-050-MC3-L23**

High Precision Bearing

Double shielded standard bearing (5x9.52x3.18)

Brand: Myonic



## 15G

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0126**

**S685A-ZZS-MC34-TH7-R-G48**

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)

Brand: MBA



## 181CB

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0097**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 181DB

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0097**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## 181L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0021**

**S683-GR7-Y**

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)

Brand: MBA



## 186B

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0136**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



**Bearing for Rotor**

**0.220.0089**

**SMR104A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)

Brand: MBA



## 187

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0136**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



**Bearing for Rotor**

**0.220.0089**

**SMR104A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



Brand: MBA





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## 188

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)  Brand: MBA	<b>0.220.0136</b>	<b>S692-ZS-GR3-KC</b>	
<b>Bearing for Rotor</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	<b>0.220.0089</b>	<b>SMR104A-ZZS-MC3-GR3-T9H</b>	



## 190

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)  Brand: MBA	<b>0.220.0136</b>	<b>S692-ZS-GR3-KC</b>	
<b>Bearing for Rotor</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	<b>0.220.0089</b>	<b>SMR104A-ZZS-MC3-GR3-T9H</b>	

## 191

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	<b>0.220.0089</b>	<b>SMR104A-ZZS-MC3-GR3-T9H</b>	

## 192

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)  Brand: MBA	<b>0.220.0136</b>	<b>S692-ZS-GR3-KC</b>	
<b>Bearing for Rotor</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	<b>0.220.0089</b>	<b>SMR104A-ZZS-MC3-GR3-T9H</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 193

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0136**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



**Bearing for Rotor**

**0.220.0089**

**SMR104A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)

Brand: MBA



## 196

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**0.220.0136**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



**Bearing for Rotor**

**0.220.0089**

**SMR104A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)

Brand: MBA



## 19E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0175**

**S683-GR7-Y**

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)

Brand: MBA



**Bearing for Shaft**


**0.220.0174**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## 19ES

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0175**

**S683-GR7-Y**

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 200

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**1.000.7185**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option  
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



## 200XD

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete**

**1.001.3520**

**BKU63520X**

For Handpiece with Wedge Type Chuck  
Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA




**Head Gear Complete**

**1.001.3520**

**BKU63520**

For Handpiece with Wedge Type Chuck  
Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



## 201

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**1.000.7185**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option  
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



## 2068D

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**0.540.2491 + 0.540.1151**

**BKU69151**

Brand: MBA



## 20E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.2004**

**S683B-ZZS-MC3-GR3**

High Precision Bearing  
Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)

Brand: MBA



**Bearing for Shaft**

**0.220.0098**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing  
Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 20LH

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0098**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



**Replacement Spindle**

**0.534.0681**

**BKU60681**

Brand: MBA



## 20LP

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0098**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



**Replacement Spindle**

**0.534.0681**

**BKU60681**

Brand: MBA



## 2304LD - 553 5121

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**220 0735**

**SMR74-ZS-T4HB-MC3-GR5-AF12**

High Speed Bearing - Torlon Retainer

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)





2 Required




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



## 2307

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)  Brand: Myonic	<b>220-0157 / 0779</b>	<b>S692-ZZS-MC34-GR5-R-G48</b>	
<b>Bearing for Chuck End</b> Stainless Bearing Deleted Item. Available while stocks last. (2x7x3)		<b>S602B-ZZS-MC3-JP6</b>	
<b>Replacement Drive Ball</b> High Precision Stainless Ball Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball 6 Required		<b>** BL-00210-S440-G3</b>	
<b>Bearing for Centre</b> Stainless Bearing Deleted Item. Available while stocks last. (4x11x4)	<b>220-0158</b>	<b>S694-MC3-JP6-AF12</b>	

## 2307L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)  Brand: MBA	<b>0.220.0144</b>	<b>S692-ZS-GR3-KC</b>	

## 2307LN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	<b>0.220.0161</b>	<b>S683-MC3-JP6-AF12</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)  Brand: MBA	<b>0.220.0174</b>	<b>SMR104A-ZZS-MC3-GR7-T9H</b>	

## 2313

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x10x4)	<b>0.220.0139</b>	<b>SMR104A-ZZ-MC3-JP6</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 2313L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0144**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



## 2313LN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0161**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## 2320

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Cap End**

**220-0157 / 0779**

**S692-ZZS-MC34-GR5-R-G48**

High Precision Bearing

Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)

Brand: Myonic



**Bearing for Chuck End**

Stainless Bearing

Deleted Item. Available while stocks last. (2x7x3)

**S602B-ZZS-MC3-JP6**




**Replacement Drive Ball**

High Precision Stainless Ball

Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball

6 Required

**\*\* BL-00210-S440-G3**



**Bearing for Centre**

**220-0158**

**S694-MC3-JP6-AF12**

Stainless Bearing

Deleted Item. Available while stocks last. (4x11x4)



## 2320L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0144**

**S692-ZS-GR3-KC**

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)

Brand: MBA



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



\*\* Available while stocks last

*Not in current production*


## 2320LN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)  Brand: Myonic	220-0157 / 0779	S692-ZZS-MC34-GR5-R-G48	
<b>Bearing for Chuck End</b> Stainless Bearing Deleted Item. Available while stocks last. (2x7x3)		S602B-ZZS-MC3-JP6	
<b>Replacement Drive Ball</b> High Precision Stainless Ball Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball 6 Required		** BL-00210-S440-G3	
<b>Bearing for Centre</b> Stainless Bearing Deleted Item. Available while stocks last. (4x11x4)	220-0158	S694-MC3-JP6-AF12	


## 2325L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)  Brand: MBA	0.220.0144	S692-ZS-GR3-KC	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	0.220.0161	S683-MC3-JP6-AF12	

## 2325LN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)  Brand: MBA	0.220.0175	S683-GR7-Y	

## 2329

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)  Brand: MBA	0.220.0175	S683-GR7-Y	

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 2329LN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0174**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## 23ES

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0098**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## 24LH

Position / Description / Application

Manufacturers Part

MBA Part

**Front Gear 1:5 Middle Gear 0.553.7201**

**0.553.8282**

**BKU68282**

Brand: MBA



## 24LN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0734**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



## 25LCA

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete**

**0.553.7341**

**BKU67341**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



**Head Gear Complete - Standard**

**0.553.7341**

**BKU67341X**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 25LH

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)	<b>0.220.0730</b>	<b>SMR5-0235-T4HB-MC3-AF12</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	<b>1.001.2836</b>	<b>SMR74-CSN-TH7-TO</b>	
<b>Drive Kit</b>  Brand: MBA	<b>553.7341 + 0.553.7201</b>	<b>BKU64101</b>	
<b>Drive Kit</b>  Brand: MBA	<b>553.7341 + 0.553.7201 + 0.553.7251</b>	<b>BKU64151</b>	
<b>Front Gear 1:5 Middle Gear 0.553.7201</b>  Brand: MBA	<b>0.553.8282</b>	<b>BKU68282</b>	
<b>Intermediate Gear 1:5</b>  Brand: MBA	<b>0.553.7251</b>	<b>BKU67251</b>	
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>0.553.7201</b>	<b>BKU67201</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.7031 + 0.553.9932</b>	<b>BKU67031</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 25LHA

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)	<b>0.220.0730</b>	<b>SMR5-0235-T4HB-MC3-AF12</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	<b>0.220.0734</b>	<b>SMR74-CSN-TH7-TO</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)	<b>0.220.0171</b>	<b>S684-MC3-JP6-AF12</b>	
<b>Drive Kit</b>  Brand: MBA	<b>553.7341 + 0.553.7201</b>	<b>BKU64101</b>	
<b>Drive Kit</b>  Brand: MBA	<b>553.7341 + 0.553.7201 + 0.553.7251</b>	<b>BKU64151</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>0.553.7341</b>	<b>BKU67341</b>	
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>0.553.7341</b>	<b>BKU67341X</b>	
<b>Intermediate Gear 1:5</b>  Brand: MBA	<b>0.553.7251</b>	<b>BKU67251</b>	
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>0.553.7201</b>	<b>BKU67201</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

<b>Replacement Press Button</b>	0.553.7031 + 0.553.9932	BKU67031	
Brand: MBA			
<b>25LHC</b>			
Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b>	0.220.0171	S684-MC3-JP6-AF12	
High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)			
<b>25LHP</b>			
Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Gear 1:5 Middle Gear 1.001.3515</b>	1.001.3516	BKU63516	
Brand: MBA			
<b>Intermediate Gear 1:5</b>	1.001.3519	BKU63519	
Brand: MBA			
<b>Replacement Middle Gear Complete 1:5</b>	1.001.3515	BKU63515	
Brand: MBA			

**Key** Usually stocked or Short Delivery \*\* Available while stocks last

Long Delivery or Large Minimum Quantities *Not in current production*

## 25LP

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Idler Gear</b> High Speed Bearing - Torlon Retainer Duplex tube bearing Abec 7 precision with high speed Torlon retainer (2x4.7x6)  Brand: MBA	1.001.2835	SMR4.762-200-ZZS-MC3-GR7-TO	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open Abec 7 precision bearing with Torlon retainer and light oil lubricant (2x5x1.5)  Brand: MBA	1.001.2834	S682-MC34-GR7-TO-L23	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open, stainless steel bearing, with extended inner ring and very narrow width, for dental applications. (3.18x6.35x1.9)  Brand: MBA	1.001.3167	SWR144C-D38W-CSN-TO	 <b>Image Coming Soon</b>
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	1.001.2836	SMR74-CSN-TH7-TO	
<b>Replacement Dog Gear 1:5</b> Gear  Brand: MBA	1.001.3071	BKU63071	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Drive Kit</b>  Brand: MBA	1.001.3520 + 1.001.3515	BKU62015	 <small>Drive Kits for Electric Driven Handpieces and Heads</small>
<b>Drive Kit</b>  Brand: MBA	1.001.3520 + 1.001.3515 + 1.001.3519	BKU62019	 <small>Drive Kits for Electric Driven Handpieces and Heads</small>
<b>Front Gear 1:5 Middle Gear 1.001.3515</b>  Brand: MBA	1.001.3516	BKU63516	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: myonic	1.001.3520	BKU63520M	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces with Carbide Insert</small>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Bearing</b></p> <p>Special Integrated Bearing Special Size - Checking Details</p> <p>Brand: MBA</p>	<p><b>1.001.2833</b></p>	<p><b>BKU22833</b></p>	
<p><b>Intermediate Gear 1:5</b></p> <p>Brand: MBA</p>	<p><b>1.001.3519</b></p>	<p><b>BKU63519</b></p>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Middle Gear Complete 1:5</b></p> <p>Brand: MBA</p>	<p><b>1.001.3515</b></p>	<p><b>BKU63515</b></p>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Press Button</b></p> <p>Brand: MBA</p>	<p><b>1.001.3284 + 0.201.8157 + 1.001.3283</b></p>	<p><b>BKU63284</b></p>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<p><b>Replacement Spring Washer</b></p> <p>Brand: MBA</p>	<p><b>1.001.2987</b></p>	<p><b>BKU62987</b></p>	 <p style="font-size: 8px; color: red;">Spring Washers for Air Driven &amp; Electric Driven Handpieces</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 25LPA


Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Idler Gear</b> High Speed Bearing - Torlon Retainer Duplex tube bearing Abec 7 precision with high speed Torlon retainer (2x4.7x6)  Brand: MBA	1.001.2835	SMR4.762-200-ZZS-MC3-GR7-TO	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open Abec 7 precision bearing with Torlon retainer and light oil lubricant (2x5x1.5)  Brand: MBA	1.001.2834	S682-MC34-GR7-TO-L23	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open, stainless steel bearing, with extended inner ring and very narrow width, for dental applications. (3.18x6.35x1.9)  Brand: MBA	1.001.3167	SWR144C-D38W-CSN-TO	 <b>Image Coming Soon</b>
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	1.001.2836	SMR74-CSN-TH7-TO	
<b>Drive Kit</b>  Brand: MBA	1.001.3520 + 1.001.3515	BKU62015	 <small>Drive Kits for Electric Driven Handpieces and Heads</small>
<b>Drive Kit</b>  Brand: MBA	1.001.3520 + 1.001.3515 + 1.001.3519	BKU62019	 <small>Drive Kits for Electric Driven Handpieces and Heads</small>
<b>Head Gear Complete</b> Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	1.001.3520	BKU63520	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces</small>
<b>Head Gear Complete</b> Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	1.001.3520	BKU63520X	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces</small>
<b>Head Gear Complete</b> Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: myonic	1.001.3520	BKU63520M	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces with Carbide Insert</small>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last


Long Delivery or Large Minimum Quantities  
Not in current production

<b>Replacement Bearing</b> Special Integrated Bearing Special Size - Checking Details  Brand: MBA	<b>1.001.2833</b>	<b>BKU22833</b>	
<b>Intermediate Gear 1:5</b>  Brand: MBA	<b>1.001.3519</b>	<b>BKU63519</b>	
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>1.001.3515</b>	<b>BKU63515</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.001.3284 + 0.201.8157 + 1.001.3283</b>	<b>BKU63284</b>	
<b>Replacement Spring Washer</b>  Brand: MBA	<b>1.001.2987</b>	<b>BKU62987</b>	


## 282A Slow Speed

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x16x5)	<b>247-6053 / 0063</b>	<b>S625-ZZ-MC3-JP6</b>	

## 29LH

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	<b>0.220.0138</b>	<b>SMR104A-ZZS-MC3-GR7-T9H</b>	



## 29LP

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	<b>0.220.0138</b>	<b>SMR104A-ZZS-MC3-GR7-T9H</b>	



<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 31E

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 3 precision bearing (2x5x1.5)	<b>1.003.5383</b>	<b>S682-MC3-GR3-GPR</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	<b>0.220.0104</b>	<b>SMR63-MC3-JP6-AF12</b>	

## 31ES

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 3 precision bearing (2x5x1.5)	<b>1.003.5383</b>	<b>S682-MC3-GR3-GPR</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)  Brand: MBA	<b>0.220.0098</b>	<b>S693-ZZS-MC34-TH5-R-L23</b>	

## 3548

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>0.220.0103</b>	<b>S684A-ZZ-MC3-JP6</b>	

## 3549

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>0.220.0103</b>	<b>S684A-ZZ-MC3-JP6</b>	

## 3550

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>0.220.0103</b>	<b>S684A-ZZ-MC3-JP6</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 3554

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0103**

**S684A-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



## 3554 Surgical Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

**Nosecone Bearing**

**SR133-MC34-GR5-R-L23**

Extra High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (2.38x4.76x1.59)

Brand: Myonic



## 3555

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0103**

**S684A-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



## 3610

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**0.220.1302**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

2 Required

Brand: MBA



## 3610N1 Staright Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**0-220-1306**

**SMR9525-050-MC3-L23**

Stainless Bearing

Double shielded standard bearing (5x9.52x3.18)

2 Required

Brand: Myonic



## 3640L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0020**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## 38A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0106**

**SMR85-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 4005

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft 7-14-5mm**

**0.220.0112**

**S687A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)

Brand: MBA



## 4500

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Spring Washer**

**1.001.6391**

**BKU66391**

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## 4500B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.006.1079</b>	<b>BKU61079</b>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.006.4564</b>	<b>BKUR210010</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.006.4564</b>	<b>BKURC210010</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.006.4564</b>	<b>BKUR2100104</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.006.4564</b>	<b>BKUR2100103</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.006.4564</b>	<b>BKUR2100105</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.006.4564</b>	<b>BKURC2100105</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.006.4564</b>	<b>BKURC2100103</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

## 465LRN

Position / Description / Application

Manufacturers Part

MBA Part

**Multiflex Coupler O-ring**

**0.200.6037+0.200.6113**

**BKU66113**

Special dental applications O-Ring

Brand: MBA



**Multiflex Coupler O-ring**

**0.553.5262**

**BKU65262**

Special dental applications O-Ring

Brand: MBA



## 465LRN Coupling

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Ball**

**0.257.4013**

**\*\* 0.257.4013**

Deleted Item. Available while stocks last. High precision ball for dental applications

Brand: KaVo



**Replacement O-Ring**

**0.200.6038**

**\*\* 0.200.6038**

Deleted Item. Available while stocks last. special o-ring for dental applications

Brand: KaVo




**Replacement Valve**

**0.553.9452**

**0.553.9452**

high precision valve for dental applications

Brand: KaVo



## 465R

Position / Description / Application

Manufacturers Part

MBA Part


**Multiflex Coupler O-ring**

**0.200.6037+0.200.6113**

**BKU66113**

Special dental applications O-Ring

Brand: MBA



## 465RN

Position / Description / Application

Manufacturers Part

MBA Part

**Multiflex Coupler O-ring**

**0.553.1872**

**BKU61872**

Special dental applications O-Ring

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 5000

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Spring Washer**

**1.001.6391**

**BKU66391**

Brand: MBA



## 5000B

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**1.006.1079**

**BKU61079**

Brand: MBA




**Replacement Turbine Excluding Bearings**

**1.006.4564**

**BKUR210010**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine Excluding Bearings**

**1.006.4564**

**BKURC210010**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine With Bearings**

**1.006.4564**

**BKUR2100105**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




**Replacement Turbine With Bearings**

**1.006.4564**

**BKUR2100104**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**1.006.4564**

**BKUR2100103**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**1.006.4564**

**BKURC2100105**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




**Replacement Turbine With Bearings**

**1.006.4564**

**BKURC2100103**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 6000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-CSN-Z-KC	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	SG1304X	SWR144-D28W-Z-KC	 <b>Image Coming Soon</b>
<b>Replacement Spring Washer</b>   Brand: MBA	1.001.5032	BKU65032	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last







Long Delivery or Large Minimum Quantities  
Not in current production

## 6000B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.000.4770</b>	<b>BKU64770</b>	 <p style="font-size: 8px;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.001.3888</b>	<b>BKUR210011</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.001.3888</b>	<b>BKURC210011</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.001.3888</b>	<b>BKUR2100115</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.3888</b>	<b>BKUR2100114</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.3888</b>	<b>BKUR2100113</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.3888</b>	<b>BKURC2100113</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.3888</b>	<b>BKURMU210011</b>	 <p style="font-size: 8px;">Ultra precision sub-assemblies, cartridges</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.001.3888</b>	<b>BKURC2100115</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

<p><b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck</p> <p>Brand: MBA</p>	<p><b>1.001.3888</b></p>	<p><b>BKURM210011</b></p>	
<p><b>614A</b></p>			
<p>Position / Description / Application</p>	<p>Manufacturers Part</p>	<p>MBA Part</p>	
<p><b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB</p>	<p><b>0.220.2007, 220-0793K</b></p>	<p><b>SR144A-D107-CSN-ZS-TO</b></p>	
<p><b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-TO</b></p>	
<p><b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-TO</b></p>	
<p><b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-KC</b></p>	
<p><b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-KC</b></p>	
<p><b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required</p>	<p><b>OD of Bearing</b></p>	<p><b>OR-00600-100-N70</b></p>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 624

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 625

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
*Not in current production*

## 625B SuperTorque

Position / Description / Application

Manufacturers Part

MBA Part

<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer</p> <p>Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required</p> <p>Brand: MPB</p>	<p><b>0.220.2007, 220-0793K</b></p>	<p><b>SR144A-D107-CSN-ZS-TO</b></p>	 <b>Image Coming Soon</b>
<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Phenolic Retainer</p> <p>Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces</p> <p>2 Required</p> <p>Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-KC</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Bearing - Torlon Retainer</p> <p>Stepped OD bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required</p> <p>Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-TO</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Bearing - Phenolic Retainer</p> <p>Stepped OD bearing with high speed phenolic retainer for precision dental handpieces</p> <p>2 Required</p> <p>Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-KC</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer</p> <p>Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required</p> <p>Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-TO</b></p>	
<p><b>Replacement O-Ring</b></p> <p>Special Dental Application O-Ring</p> <p>Special sized o-ring used in dental Applications</p> <p>2 Required</p>	<p><b>OD of Bearing</b></p>	<p><b>OR-00600-100-N70</b></p>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last




Long Delivery or Large Minimum Quantities

*Not in current production*

## 625C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

## 625D

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production





## 627

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	





Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 628 - 629 FG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-TO	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-TO	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-KC	 <b>Image Coming Soon</b>

## 628 BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-TO	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-TO	 <b>Image Coming Soon</b>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 630

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last








Long Delivery or Large Minimum Quantities  
*Not in current production*

## 630B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.5881</b>	<b>BKUR210001</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.5881</b>	<b>BKURC210001</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.5881</b>	<b>BKUR2100015</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>0.553.5881</b></p>	<p><b>BKUR2100014</b></p>	
<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>0.553.5881</b></p>	<p><b>BKUR2100013</b></p>	
<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>0.553.5881</b></p>	<p><b>BKUR210001T</b></p>	<p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>  <p>Triple pack - Unit with 3 pieces</p>
<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer</p> <p>Brand: MBA</p>	<p><b>0.553.5881</b></p>	<p><b>BKURC2100015</b></p>	<p>Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p> 
<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>0.553.5881</b></p>	<p><b>BKURC2100013</b></p>	<p>Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p> 
<p><b>Replacement Turbine With Bearings</b>                  For Handpiece with Wedge Type Chuck                  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>0.553.5881</b></p>	<p><b>BKURMU210001</b></p>	<p>Ultra precision sub-assemblies, cartridges</p> 
<p><b>Turbine Sub Assembly Without Bearings</b>                  For Handpiece with Wedge Type Chuck</p> <p>Brand: MBA</p>	<p><b>0.553.5881</b></p>	<p><b>BKURM210001</b></p>	<p>Ultra precision sub-assemblies, cartridges</p> 

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 630C PB SuperTorque

Position / Description / Application

Manufacturers Part

MBA Part

<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required Brand: MPB</p>	<p><b>0.220.2007, 220-0793K</b></p>	<p><b>SR144A-D107-CSN-ZS-TO</b></p>	 <b>Image Coming Soon</b>
<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-TO</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces</p> <p>2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-KC</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces</p> <p>2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-ZS-TO</b></p>	
<p><b>Replacement Bearing</b></p> <p>High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces</p> <p>2 Required Brand: Myonic</p>	<p><b>220-0125</b></p>	<p><b>SR144A-D105-CSN-ZS-KC</b></p>	
<p><b>Replacement O-Ring</b></p> <p>Special Dental Application O-Ring Special sized o-ring used in dental Applications</p> <p>2 Required</p>	<p><b>OD of Bearing</b></p>	<p><b>OR-00600-100-N70</b></p>	


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 632D

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

## 632D - 632E

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 632E

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

## 634

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	<b>0.220.0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	<b>0.220.0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	<b>0.220.0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	<b>0.220.0125</b>	<b>SR144A-D105-ZS-KC</b>	

Key

**Usually stocked or Short Delivery**

**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities

*Not in current production*

## 635

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**0.220.2007, 220-0793K**

**SR144A-D107-CSN-ZS-TO**

High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer  
Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces

Brand: MPB



## 635 Back

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-ZS-TO**

High Speed Bearing - Torlon Retainer  
3.175x6.35x2.381/2.778, stainless steel, Torlon retainer

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-CSN-ZS-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, ceramic balls, and phenolic cage, for dental applications

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-CSN-ZS-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, ceramic balls, and Torlon cage, for dental applications

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-ZS-KC**

High Speed Bearing - Phenolic Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, and phenolic cage, for dental applications

Brand: MBA



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>





## 635B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.8491</b>	<b>BKU68491</b>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8581</b>	<b>BKUR210003</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8581</b>	<b>BKUR2100034</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.8581</b>	<b>BKUR2100035</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8581</b>	<b>BKUR2100033</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKURC210003</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKURC2100033</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKURC2100035</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 635BR

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	0.553.8661	BKUR210005	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	0.553.8661	BKUR2100055	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	0.553.8661	BKUR2100053	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	0.553.8661	BKUR2100054	





## 636

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 636B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.8061</b>	<b>BKU68061</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8911</b>	<b>BKURC210004</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	<b>0.553.8911</b>	<b>BKURC2100045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8911</b>	<b>BKURC2100043</b>	

## 636CP

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8911</b>	<b>BKUR210004</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8911</b>	<b>BKUR2100043</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8911</b>	<b>BKUR210004T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8911</b>	<b>BKUR2100044</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.8911</b>	<b>BKURE2100045</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 637

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**0.220.2007, 220-0793K**

**SR144A-D107-CSN-ZS-TO**

High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer  
Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces

Brand: MPB



## 637 Back

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-ZS-TO**

High Speed Bearing - Torlon Retainer  
3.175x6.35x2.381/2.778, stainless steel, Torlon retainer

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-CSN-ZS-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, ceramic balls, and Torlon cage, for dental applications

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-ZS-KC**

High Speed Bearing - Phenolic Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, and phenolic cage, for dental applications

Brand: MBA



**Bearing for Rotor Cartridge**

**0.220.2015**

**SWR144-D106W-CSN-ZS-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer  
Open, stainless steel bearing, with extended inner ring and stepped OD, ceramic balls, and phenolic cage, for dental applications

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



## 637B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.8491</b>	<b>BKU68491</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8661</b>	<b>BKUR210005</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8661</b>	<b>BKUR2100053</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.8661</b>	<b>BKUR2100055</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8661</b>	<b>BKUR2100054</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKURC210005</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKURC2100055</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKURC2100053</b>	







**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 637B Mini Lux 3 BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.		SR144A-D650-CSN-ZS-TO	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		OR-00600-100-N70	

## 639A Friction Grip

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-ZS-TO	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-ZS-KC	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-CSN-ZS-KC	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-CSN-ZS-TO	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 639B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.6911</b>	<b>BKUR210006</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.6911</b>	<b>BKUR2100063</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.6911</b>	<b>BKUR2100064</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last


Long Delivery or Large Minimum Quantities  
Not in current production

**Replacement Turbine With Bearings**
**0.553.6911**
**BKUR2100065**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



## 639C Push Button

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**
**0.220.2007, 220-0793K**
**SR144A-D107-CSN-ZS-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces

2 Required

Brand: MPB


**Replacement Bearing**
**220-0125**
**SR144A-D105-ZS-TO**

High Speed Bearing - Torlon Retainer

Stepped OD bearing with Torlon retainer for precision dental handpieces

2 Required

Brand: Myonic


**Replacement Bearing**
**220-0125**
**SR144A-D105-CSN-ZS-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces

2 Required

Brand: Myonic


**Replacement Bearing**
**220-0125**
**SR144A-D105-CSN-ZS-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces

2 Required

Brand: Myonic


**Replacement Bearing**
**220-0125**
**SR144A-D105-ZS-KC**

High Speed Bearing - Phenolic Retainer

Stepped OD bearing with high speed phenolic retainer for precision dental handpieces

2 Required

Brand: Myonic


**Replacement O-Ring**
**OD of Bearing**
**OR-00600-100-N70**

Special Dental Application O-Ring

Special sized o-ring used in dental Applications

2 Required


**Key**
**Usually stocked or Short Delivery**

Long Delivery or Large Minimum Quantities

**\*\* Available while stocks last**
*Not in current production*

## 640

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>

## 640A Friction Grip

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 640B



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.6511</b>	<b>BKU66511</b>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Spray Cover</b>  Brand: MBA	<b>0.553.8752</b>	<b>BKU68752</b>	 <p style="font-size: 8px; color: red;">Spray Covers</p>
<b>Replacement Threaded Ring</b>  Brand: MBA	<b>0.553.3322</b>	<b>BKU63322</b>	 <p style="font-size: 8px; color: red;">Spray Covers</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.5881</b>	<b>BKUR210001</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.5881</b>	<b>BKURC210001</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.5881</b>	<b>BKUR2100014</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.5881</b>	<b>BKUR2100013</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.5881</b>	<b>BKUR2100015</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Turbine With Bearings</b>      <b>0.553.5881</b></p> <p>For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>BKUR210001T</b></p>	
<p><b>Replacement Turbine With Bearings</b>      <b>0.553.5881</b></p> <p>For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer</p> <p>Brand: MBA</p>	<p><b>BKURC2100015</b></p>	
<p><b>Replacement Turbine With Bearings</b>      <b>0.553.5881</b></p> <p>For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>BKURMU210001</b></p>	
<p><b>Replacement Turbine With Bearings</b>      <b>0.553.5881</b></p> <p>For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>BKURC2100013</b></p>	
<p><b>Turbine Sub Assembly Without Bearings</b>      <b>0.553.5881</b></p> <p>For Handpiece with Wedge Type Chuck</p> <p>Brand: MBA</p>	<p><b>BKURM210001</b></p>	

## 640C

Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Replacement Spray Cover</b>      <b>0.553.8752</b></p> <p>Brand: MBA</p>	<p><b>BKU68752</b></p>		
<p><b>Replacement Threaded Ring</b>      <b>0.553.3322</b></p> <p>Brand: MBA</p>	<p><b>BKU63322</b></p>		

## 640C PB SuperTorque

Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Replacement O-Ring</b>      <b>OD of Bearing</b></p> <p>Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required</p>	<p><b>OR-00600-100-N70</b></p>		

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## 640FG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 642A FG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 642B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.7391</b>	<b>BKU67391</b>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.6911</b>	<b>BKUR210006</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.6911</b>	<b>BKURC210006</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.6911</b>	<b>BKUR2100065</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.6911</b>	<b>BKUR2100064</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.6911</b>	<b>BKUR2100063</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.6911</b>	<b>BKURC2100065</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.6911</b>	<b>BKURC2100063</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 642B - 642C PB BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	<b>OD of Bearing</b>	<b>OR-00600-100-N70</b>	





Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 645

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>

## 645B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	0.553.7391	BKU67391	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	0.553.6911	BKURC210006	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	0.553.6911	BKURC2100063	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	0.553.6911	BKURC2100065	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 646

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 646B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b> Brand: MBA	<b>0.553.8061</b>	<b>BKU68061</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8911</b>	<b>BKUR210004</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8911</b>	<b>BKURC210004</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8911</b>	<b>BKUR2100044</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8911</b>	<b>BKUR2100043</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8911</b>	<b>BKUR210004T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	<b>0.553.8911</b>	<b>BKURC2100045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>0.553.8911</b>	<b>BKURE2100045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8911</b>	<b>BKURC2100043</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 647

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-ZS-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces  Brand: Myonic	0.220.0125	SR144A-D105-CSN-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces  Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>


Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 647B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b> Brand: MBA	<b>0.553.8061</b>	<b>BKU68061</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8541</b>	<b>BKUR210007</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8541</b>	<b>BKURC210007</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKUR2100073</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKUR2100074</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>0.553.8541</b>	<b>BKUR2100075</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>0.553.8541</b>	<b>BKURC2100075</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKURC2100073</b>	

## 649

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 649B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b> Brand: MBA	<b>0.553.8061</b>	<b>BKU68061</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8541</b>	<b>BKUR210007</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>0.553.8541</b>	<b>BKURC210007</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKUR2100074</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>0.553.8541</b>	<b>BKUR2100075</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKUR2100073</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>0.553.8541</b>	<b>BKURC2100073</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>0.553.8541</b>	<b>BKURC2100075</b>	



## 650

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.	<b>0.220.2009</b>	<b>SR144A-D650-CSN-ZS-TO</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 6500

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-CSN-Z-KC	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	SG1304X	SWR144-D28W-Z-KC	 <b>Image Coming Soon</b>
<b>Replacement Spring Washer</b>  Brand: MBA	1.003.5900	BKU65900	
<b>Replacement Spring Washer</b>  Brand: MBA	1.003.7576	BKU67576	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 6500B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1.003.9664	BKUR210012	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1.003.9664	BKURC210015	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.003.9664	BKUR2100124	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1.003.9664	BKUR2100125	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.003.9664	BKUR2100123	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.003.9664	BKURC2100153	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.003.9664	BKURMU210012	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1.003.9664	BKURC2100155	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	1.003.9664	BKURM210012	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
Not in current production

## 6500BR



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)	<b>1.000.7980</b>	<b>SR144A-D31-CSN-WZ-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Press Button</b>  Brand: MBA	<b>1.003.5525</b>	<b>BKU65525</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.003.6884</b>	<b>BKUR210014</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.003.6884</b>	<b>BKURC210014</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.003.6884</b>	<b>BKUR2100145</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR2100143</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR210014T</b>	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR2100144</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.003.6884</b>	<b>BKURC2100145</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKURC2100143</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKURMU210014</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>1.003.6884</b>	<b>BKURM210014</b>	

## 650A FG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Chuck End</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.	<b>220-2009</b>	<b>SR144A-D650-CSN-ZS-TO</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		<b>OR-00600-100-N70</b>	

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>



Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 650B



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>0.553.8061</b>	<b>BKU68061</b>	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Spray Cover</b>  Brand: MBA	<b>0.553.0573</b>	<b>BKU60573</b>	 <p style="font-size: 8px; color: red;">Spray Covers</p>
<b>Replacement Threaded Ring</b>  Brand: MBA	<b>0.553.0543</b>	<b>BKU60543</b>	 <p style="font-size: 8px; color: red;">Spray Covers</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8051</b>	<b>BKUR210008</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>0.553.8051</b>	<b>BKURC210008</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8051</b>	<b>BKUR2100083</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8051</b>	<b>BKUR2100084</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>0.553.8051</b>	<b>BKUR2100085</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>0.553.8051</b>	<b>BKURC2100083</b>	 <p style="font-size: 8px; color: red;">Precision sub- assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last



Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Turbine With Bearings</b>                      <b>0.553.8051</b></p> <p>For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>BKURMU210008</b></p>	 <p style="font-size: 8px; color: red;">Ultra precision sub-assemblies, cartridges</p>
<p><b>Turbine Sub Assembly Without Bearings</b>                      <b>0.553.8051</b></p> <p>For Handpiece with Wedge Type Chuck</p> <p>Brand: MBA</p>	<p><b>BKURM210008</b></p>	 <p style="font-size: 8px; color: red;">Ultra precision sub-assemblies, cartridges</p>



## 650B PB - 650C PB

Position / Description / Application	Manufacturers Part	MBA Part
<p><b>Bearing for Chuck End</b>                      <b>220-2009</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.</p>	<p><b>SR144A-D650-CSN-ZS-TO</b></p>	
<p><b>Replacement O-Ring</b></p> <p>Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required</p>	<p><b>OR-00600-100-N70</b></p>	

## 650C

Position / Description / Application	Manufacturers Part	MBA Part
<p><b>Replacement Spray Cover</b>                      <b>0.553.0573</b></p> <p>Brand: MBA</p>	<p><b>BKU60573</b></p>	 <p style="font-size: 8px; color: red;">Spray Covers</p>
<p><b>Replacement Threaded Ring</b>                      <b>0.553.0543</b></p> <p>Brand: MBA</p>	<p><b>BKU60543</b></p>	 <p style="font-size: 8px; color: red;">Spray Covers</p>

## 655

Position / Description / Application	Manufacturers Part	MBA Part
<p><b>Bearing for Rotor Cartridge</b>                      <b>1.000.7980</b></p> <p>High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)</p>	<p><b>SR144A-D31-CSN-WZ-TO</b></p>	
<p><b>Replacement Spring Washer</b>                      <b>1.001.5032</b></p> <p>Brand: MBA</p>	<p><b>BKU65032</b></p>	 <p style="font-size: 8px; color: red;">Spring Washers for Air Driven &amp; Electric Driven Handpieces</p>

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>





Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 655B



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.000.2076</b>	<b>BKU62076</b>	
<b>Replacement Spray Cover</b>  Brand: MBA	<b>0.553.0573</b>	<b>BKU60573</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.000.2078</b>	<b>BKUR210009</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.000.2078</b>	<b>BKURC210009</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.000.2078</b>	<b>BKUR2100095</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.000.2078</b>	<b>BKUR2100093</b>	
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.000.2078</b>	<b>BKUR210009T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.000.2078</b>	<b>BKUR2100094</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.000.2078</b>	<b>BKURMU210009</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer</p> <p>Brand: MBA</p>	<p><b>1.000.2078</b></p>	<p><b>BKURC2100095</b></p>	
<p><b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p> <p>Brand: myonic</p>	<p><b>1.000.2078</b></p>	<p><b>BKURC2100093</b></p>	
<p><b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck</p> <p>Brand: MBA</p>	<p><b>1.000.2078</b></p>	<p><b>BKURM210009</b></p>	
<p><b>Replacement Threaded Ring</b></p> <p>Brand: MBA</p>	<p><b>1.000.3391</b></p>	<p><b>BKU63391</b></p>	

## 655C

Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Replacement Spray Cover</b></p> <p>Brand: MBA</p>	<p><b>0.553.0573</b></p>	<p><b>BKU60573</b></p>	
<p><b>Replacement Threaded Ring</b></p> <p>Brand: MBA</p>	<p><b>1.000.3391</b></p>	<p><b>BKU63391</b></p>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 659BR

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Spray Cover</b> Brand: MBA	<b>0.553.0573</b>	<b>BKU60573</b>	 <p style="font-size: 8px; text-align: center;">Spray Covers</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>1.000.2078</b>	<b>BKUR210009</b>	 <p style="font-size: 8px; text-align: center;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>1.000.2078</b>	<b>BKUR210009T</b>	 <p style="font-size: 8px; text-align: center;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>1.000.2078</b>	<b>BKUR2100094</b>	 <p style="font-size: 8px; text-align: center;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>1.000.2078</b>	<b>BKUR2100093</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>1.000.2078</b>	<b>BKUR2100095</b>	 <p style="font-size: 8px; text-align: center;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>1.000.2078</b>	<b>BKURMU210009</b>	 <p style="font-size: 8px; text-align: center;">Ultra precision sub-assemblies, cartridges</p>
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck Brand: MBA	<b>1.000.2078</b>	<b>BKURM210009</b>	 <p style="font-size: 8px; text-align: center;">Ultra precision sub-assemblies, cartridges</p>

## 660

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Spring Washer</b> Brand: MBA	<b>1.001.5032</b>	<b>BKU65032</b>	 <p style="font-size: 8px; text-align: center;">Spring Washers for Air Driven &amp; Electric Driven Handpieces</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 660B

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**1.000.2076**

**BKU62076**

Brand: MBA




**Replacement Turbine Excluding Bearings**

**1.000.2078**

**BKURC210009**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine With Bearings**

**1.000.2078**

**BKURC2100095**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




**Replacement Turbine With Bearings**

**1.000.2078**

**BKURC2100093**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## 66GD

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0105**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing  
Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## 66LU

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.2005**

**SMR74-MC3-JP6-AF12**

High Precision Bearing  
Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



## 67G

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0104**

**SMR63-MC3-JP6-AF12**

High Precision Bearing  
Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 67GD

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0105**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## 67L

Position / Description / Application

Manufacturers Part

MBA Part


**Bearing for Headgear**

**0.220.0051**

**SWA683B-ZZ-C69-GR5-TO**

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



## 67LDN

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0713**

**S683A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (3x7x3)



**Bearing for Headgear**

**0.220.0105**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



\*\* Available while stocks last

*Not in current production*

## 67LH

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	0.220.2005	SMR74-MC3-JP6-AF12	
<b>Drive Kit</b>  Brand: MBA	0.540.3671 + 0.540.1161	BKU67151	
<b>Replacement Head Housing</b>  Brand: MBA	0.540.3641	BKU63641	
<b>Head Gear Complete 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	0.540.3671	BKU63671	
<b>Replacement Head Gear 2:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	0.540.4022	BKU64022	
<b>Intermediate Piece Dental Head 2:1</b>  Brand: MBA	0.540.1161	BKU61161	

## 67RIC

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	0.220.0713	S683A-ZZ-MC3-JP6	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0023	SMR85-MC3-JP6-AF12	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 68B - 68C - 68CN - 68D

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)

2 Required

**SMR63-MC3-JP6-AF12**



## 68CN

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**0.540.1781 + 0.540.1151**

**BKU68151**

Brand: MBA



### Head Gear Complete 1:1

For Handpiece with Latch Type Chuck

**0.540.1781**

**BKU61781**

Brand: MBA



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 68G

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	<b>0.220.0104</b>	<b>SMR63-MC3-JP6-AF12</b>	
<b>Bush</b> Intermediate Piece 0.540.1151  Brand: MBA	<b>0.540.2392</b>	<b>BKU62392</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Drive Kit</b>  Brand: MBA	<b>0.540.1781 + 0.540.1151</b>	<b>BKU68151</b>	 <p style="font-size: 8px; color: red;">Drive Kits for Electric Driven Handpieces and Heads</p>
<b>Fitting set</b> Includes 0.540.4072, 0.540.4102, 0.540.4522  Brand: MBA	<b>0.540.4072, 0.540.4102, 0.540.4522</b>	<b>BKU64102</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Head Gear Complete 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.1781</b>	<b>BKU61781</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Intermediate Piece Dental Head 1:1</b>  Brand: MBA	<b>0.540.1151</b>	<b>BKU61151</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Screw</b>  Brand: MBA	<b>0.201.0005</b>	<b>BKU60005</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Spring</b>  Brand: MBA	<b>0.540.3162</b>	<b>BKU63162</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 68GDN

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 1:1**

**0.540.2491**

**BKU62491**

For Handpiece with Latch Type Chuck

Brand: MBA




**Replacement Head Gear 1:1**

**0.540.4092**

**BKU64092**

For Handpiece with Latch Type Chuck

Brand: MBA



## 68L

Position / Description / Application

Manufacturers Part

MBA Part


**Bearing for Headgear**

**0.220.0051**

**SWA683B-ZZ-C69-GR5-TO**

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 68LDN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	<b>0.220.0713</b>	<b>S683A-ZZ-MC3-JP6</b>	
<b>Bearing for Cap End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (3x7x3)	<b>220-0713 / 0102</b>	<b>S683A-ZZ-MC3-JP6</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)  Brand: Myonic	<b>0.220.0105</b>	<b>SMR74-MC34-GR7-R-L23</b>	
<b>Bearing for Chuck End</b> Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>220-0105</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Bush</b> Intermediate Piece 0.540.1151  Brand: MBA	<b>0.540.2392</b>	<b>BKU62392</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Drive Kit</b>  Brand: MBA	<b>0.540.2491 + 0.540.1151</b>	<b>BKU69151</b>	 <p style="font-size: 8px; color: red;">Drive Kits for Electric Driven Handpieces and Heads</p>
<b>Fitting set</b> Includes 0.540.4072, 0.540.4102, 0.540.4522  Brand: MBA	<b>0.540.4072, 0.540.4102, 0.540.4522</b>	<b>BKU64102</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Head Housing</b>  Brand: MBA	<b>0.540.2231</b>	<b>BKU62231</b>	 <p style="font-size: 8px; color: red;">Head Housings</p>
<b>Head Gear Complete 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.2491</b>	<b>BKU62491</b>	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Head Gear 1:1</b> For Handpiece with Latch Type Chuck</p> <p>Brand: MBA</p>	<p><b>0.540.4092</b></p>	<p><b>BKU64092</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Intermediate Piece Dental Head 1:1</b></p> <p>Brand: MBA</p>	<p><b>0.540.1151</b></p>	<p><b>BKU61151</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Press Button</b></p> <p>Brand: MBA</p>	<p><b>0.540.2501</b></p>	<p><b>BKU62501</b></p>	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<p><b>Replacement Screws</b></p> <p>Brand: MBA</p>	<p><b>0.201.0005</b></p>	<p><b>BKU60005</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Spring</b></p> <p>Brand: MBA</p>	<p><b>0.540.3162</b></p>	<p><b>BKU63162</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
*Not in current production*

## 68LH

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Front Bearing</b> High Speed Bearing - Torlon Retainer Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)		<b>SWA683B-ZZ-C69-GR5-TO</b>	
<b>Replacement Rear Bearing</b> Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)		<b>S683-MC3-JP6-AF12</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>0.220.2005</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Bush</b> Intermediate Piece 0.540.1151  Brand: MBA	<b>0.540.2392</b>	<b>BKU62392</b>	
<b>Drive Kit</b>  Brand: MBA	<b>0.540.3661 + 0.540.1151</b>	<b>BKU66151</b>	
<b>Fitting set</b> Includes 0.540.4072, 0.540.4102, 0.540.4522  Brand: MBA	<b>0.540.4072, 0.540.4102, 0.540.4522</b>	<b>BKU64102</b>	
<b>Replacement Head Housing</b>  Brand: MBA	<b>0.540.3631</b>	<b>BKU63631</b>	
<b>Head Gear Complete 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.3661</b>	<b>BKU63661</b>	
<b>Replacement Head Gear 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.4012</b>	<b>BKU64012</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

<p><b>Intermediate Piece Dental Head 1:1</b></p> <p>Brand: MBA</p>	<p><b>0.540.1151</b></p>	<p><b>BKU61151</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Press Button</b></p> <p>Brand: MBA</p>	<p><b>0.540.3741</b></p>	<p><b>BKU63741</b></p>	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<p><b>Replacement Screws</b></p> <p>Brand: MBA</p>	<p><b>0.201.0005</b></p>	<p><b>BKU60005</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<p><b>Replacement Spring</b></p> <p>Brand: MBA</p>	<p><b>0.540.3162</b></p>	<p><b>BKU63162</b></p>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 68LU

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>0.220.2005</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Bush</b> Intermediate Piece 0.540.1151  Brand: MBA	<b>0.540.2392</b>	<b>BKU62392</b>	
<b>Fitting set</b> Includes 0.540.4072, 0.540.4102, 0.540.4522  Brand: MBA	<b>0.540.4072, 0.540.4102, 0.540.4522</b>	<b>BKU64102</b>	
<b>Head Gear Complete 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.3661</b>	<b>BKU63661</b>	
<b>Replacement Head Gear 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>0.540.4012</b>	<b>BKU64012</b>	
<b>Intermediate Piece Dental Head 1:1</b>  Brand: MBA	<b>0.540.1151</b>	<b>BKU61151</b>	
<b>Replacement Screws</b>  Brand: MBA	<b>0.201.0005</b>	<b>BKU60005</b>	
<b>Replacement Spring</b>  Brand: MBA	<b>0.540.3162</b>	<b>BKU63162</b>	

## 6N

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>0.220.0103</b>	<b>S684A-ZZ-MC3-JP6</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 7000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	<b>1.000.7979</b>	<b>SWR144-D28W-CSN-Z-KC</b>	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	<b>1.000.7979</b>	<b>SWR144-D28W-Z-TO</b>	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	<b>1.000.7979</b>	<b>SWR144-D28W-CSN-Z-TO</b>	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	<b>SG1304X</b>	<b>SWR144-D28W-Z-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Spring Washer</b>   Brand: MBA	<b>1.001.5032</b>	<b>BKU65032</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last





Long Delivery or Large Minimum Quantities  
Not in current production

## 7000B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	1.000.4770	BKU64770	 <p style="font-size: 8px;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1.000.4772	BKUR210013	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1.000.4772	BKURC210013	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1.000.4772	BKUR2100135	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.000.4772	BKUR2100133	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.000.4772	BKUR2100134	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1.000.4772	BKURC2100135	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.000.4772	BKURMU210013	 <p style="font-size: 8px;">Ultra precision sub-assemblies, cartridges</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.000.4772	BKURC2100133	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
Not in current production

<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck	<b>1.000.4772</b>	<b>BKURM210013</b>	
Brand: MBA			
<b>700KL</b>			
Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)	<b>0.220.0218</b>	<b>S693-ZZS-MC34-TH5-R-L23</b>	
Brand: MBA			
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)	<b>0.220.0219</b>	<b>S685A-ZZS-MC34-TH7-R-G48</b>	
Brand: MBA			
<b>79G</b>			
Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	<b>0.220.0105</b>	<b>SMR74-MC34-GR7-R-L23</b>	
Brand: Myonic			
<b>7LH</b>			
Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	<b>0.220.0139</b>	<b>SMR104A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	<b>0.220.0174</b>	<b>SMR104A-ZZS-MC3-GR7-T9H</b>	
Brand: MBA			
<b>Replacement Gear Unit</b>	<b>0.534.1501</b>	<b>BKU61501</b>	
Brand: MBA			

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production


## 7LP

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	<b>0.220.0139</b>	<b>SMR104A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)  Brand: MBA	<b>0.220.0174</b>	<b>SMR104A-ZZS-MC3-GR7-T9H</b>	
<b>Replacement Gear Unit</b>  Brand: MBA	<b>0.534.1501</b>	<b>BKU61501</b>	 <p style="font-size: small; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>

## 8000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Spring Washer</b>  Brand: MBA	<b>1.003.5900</b>	<b>BKU65900</b>	 <p style="font-size: small; color: red;">Spring Washers for Air Driven &amp; Electric Driven Handpieces</p>

## 8000 Front

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)	<b>1.000.7980</b>	<b>SR144A-D31-CSN-WZ-TO</b>	 <p style="font-size: small; color: red;">Image Coming Soon</p>

Key Usually stocked or Short Delivery \*\* Available while stocks last



Long Delivery or Large Minimum Quantities *Not in current production*

## 8000B



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.003.5525</b>	<b>BKU65525</b>	 <p style="font-size: 8px;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.003.6884</b>	<b>BKUR210014</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.003.6884</b>	<b>BKURC210014</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR2100143</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR2100144</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKUR210014T</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p> <p style="font-size: 8px;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.003.6884</b>	<b>BKUR2100145</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKURMU210014</b>	 <p style="font-size: 8px;">Ultra precision sub-assemblies, cartridges</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.003.6884</b>	<b>BKURC2100143</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last


Long Delivery or Large Minimum Quantities  
Not in current production

<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.003.6884</b>	<b>BKURC2100145</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>1.003.6884</b>	<b>BKURM210014</b>	


## 80L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Single shielded, Abec 7 precision stainless steel bearing with special clearance and bronze cage (4x7x2)  Brand: MBA	<b>0.220.0735</b>	<b>SMR74-ZS-MC11-15-GR5-Y-L23</b>	
<b>Drive Kit</b>  Brand: MBA	<b>1.009.1782 + 0.540.1151</b>	<b>BKU68251</b>	


## 80LD

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Single shielded, Abec 7 precision stainless steel bearing with special clearance and bronze cage (4x7x2)  Brand: MBA	<b>0.220.0735</b>	<b>SMR74-ZS-MC11-15-GR5-Y-L23</b>	

## 81

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	<b>0.220.0106</b>	<b>SMR85-MC3-JP6-AF12</b>	



## C09

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	<b>0.220.0020</b>	<b>S683-MC3-JP6-AF12</b>	

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## CL09

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	<b>0.220.0023</b>	<b>SMR85-MC3-JP6-AF12</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x11x3)	<b>1.001.2581</b>	<b>S685-MC3-JP6-AF12</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## CL3

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	0.220.0051	SWA683B-ZZ-C69-GR5-TO	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0023	SMR85-MC3-JP6-AF12	
<b>Drive Kit</b>  Brand: MBA	1.011.5595 + 1.003.4207	BKU69507	 <p style="font-size: 8px; color: red;">Drive Kits for Electric Driven Handpieces and Heads</p>
<b>Head Gear Complete 3:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	1.001.5005	BKU65005	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Head Housing</b>  Brand: MBA	1.001.9257	BKU69257	 <p style="font-size: 8px; color: red;">Head Housings</p>
<b>Head Gear Complete 3:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	1.011.5595	BKU65595	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Intermediate Piece Dental Head 3:1</b>  Brand: MBA	1.003.4207	BKU64207	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Press Button</b>  Brand: MBA	1.001.4486+0.540.3752+0.201.8958	BKU64486	 <p style="font-size: 8px; color: red;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>

## Coupling - Low Speed

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> High Temperature Special sized o-ring used in dental Applications 2 Required	200-6159	OR-00830-075-N70	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## E 25C

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete

**1.008.4040**

**BKU64040**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



## E120

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Headgear

**0.220.2005**

**SMR74-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



## E15C

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**0.540.3671 + 0.540.1161**

**BKU67151**

Brand: MBA



## E15L

Position / Description / Application

Manufacturers Part

MBA Part


### Bearing for Headgear

**0.220.0051**

**SWA683B-ZZ-C69-GR5-TO**

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



### Replacement Gear Unit

**0.534.1501**

**BKU61501**

Brand: MBA



### Head Gear Complete 1:1

**0.540.3671**

**BKU63671**

For Handpiece with Latch Type Chuck

Brand: MBA



### Replacement Head Gear 2:1

**0.540.4022**

**BKU64022**

For Handpiece with Latch Type Chuck

Brand: MBA



### Intermediate Piece Dental Head 2:1

**0.540.1161**

**BKU61161**

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## E20C

Position / Description / Application

Manufacturers Part


MBA Part

**Drive Kit**

**0.540.3661 + 0.540.1151**

**BKU66151**

Brand: MBA



## E20L

Position / Description / Application

Manufacturers Part

MBA Part


**Bearing for Headgear**

**0.220.0051**

**SWA683B-ZZ-C69-GR5-TO**

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



## E25

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**1.008.4040 + 1.007.9916**

**BKU64016**

Brand: MBA



**Drive Kit**

**1.014.5555 + 1.014.5557**

**BKU65557**

Brand: MBA




**Front Gear 1:5 Middle Gear**

**1.007.9916**

**BKU69999**

Brand: MBA



## E25C

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete - Standard**

**1.008.4040**

**BKU64040X**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



**Replacement Press Button**

**1.008.1488**

**BKU61488**

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## E25L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open, stainless steel bearing, with extended inner ring and very narrow width, for dental applications. (3.18x6.35x1.9)  Brand: MBA	1.008.4038	SWR144C-D38W-CSN-TO	 <b>Image Coming Soon</b>
<b>Bearing for Headgear</b> Japanese High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 3 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	1.008.4039	SMR74-CSN-J3-TO	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	1.008.4039	SMR74-CSN-TH7-TO	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	1.008.4040	BKU64040	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces</small>
<b>Head Gear Complete - Premium</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: myonic	1.008.4040	BKU64040M	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces with Carbide Insert</small>
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	1.008.4040	BKU64040X	 <small>Head Gears Complete for 1:5 Electric Driven Handpieces</small>
<b>Replacement Dog Gear 1:5</b> Gear  Brand: MBA	1.007.5555	BKU65555	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	1.007.9916	BKU69916	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Press Button</b>  Brand: MBA	1.008.1488	BKU61488	 <small>Press Button cpl. for Heads &amp; Air Driven Handpieces</small>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

**Replacement Spring Washer**

**1.001.2987**

**BKU62987**

Brand: MBA



## E31

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**0.220.0104**

**SMR63-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



## E40

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.220.0020**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



**Bearing for Headgear**

**0.220.2005**

**SMR74-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



## E675

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**1.000.7979**

**SWR144-D28W-CSN-Z-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA



**Bearing for Rotor Cartridge**

**1.000.7979**

**SWR144-D28W-Z-TO**

High Speed Bearing - Torlon Retainer

Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)

Brand: MBA



**Bearing for Rotor Cartridge**

**1.000.7979**

**SWR144-D28W-CSN-Z-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)



**Bearing for Rotor Cartridge**

**SG1304X**

**SWR144-D28W-Z-KC**

High Speed Bearing - Phenolic Retainer

Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## E675L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.7051</b>	<b>BKU67051</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9457</b>	<b>BKUR210017</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9457</b>	<b>BKURC210017</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKUR2100174</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKUR2100173</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	<b>1.007.9457</b>	<b>BKUR2100175</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9457</b>	<b>BKURC2100175</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKURC2100173</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## E677

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-CSN-Z-KC	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	1.000.7979	SWR144-D28W-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	SG1304X	SWR144-D28W-Z-KC	 <b>Image Coming Soon</b>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## E677L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.7051</b>	<b>BKU67051</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9457</b>	<b>BKUR210017</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9457</b>	<b>BKURC210017</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKUR2100174</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKUR2100173</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	<b>1.007.9457</b>	<b>BKUR2100175</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9457</b>	<b>BKURC2100173</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9457</b>	<b>BKURC2100175</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## E679

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**1.000.7980**

**SR144A-D31-CSN-WZ-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)



**Replacement Spring Washer**

**1.003.7576**

**BKU67576**

Brand: MBA



Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## E679L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.8903</b>	<b>BKU68903</b>	 <p style="font-size: 8px;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKUR210015</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKURC210018</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR210015T</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p> <p style="font-size: 8px;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR2100154</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKUR2100155</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR2100153</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKURC2100185</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKURC2100183</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## E680

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**1.000.7980**

**SR144A-D31-CSN-WZ-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision dental handpiece ceramic-hybrid angular contact bearing with Torlon retainer (3.18x6.35x2.78)



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*



## E680L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.8903</b>	<b>BKU68903</b>	 <p style="font-size: 8px;">Press Button cpl. for Heads &amp; Air Driven Handpieces</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKUR210015</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.007.9313</b>	<b>BKURC210018</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR2100153</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKUR2100155</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR2100154</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKUR210015T</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p> <p style="font-size: 8px;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.007.9313</b>	<b>BKURC2100183</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.007.9313</b>	<b>BKURC2100185</b>	 <p style="font-size: 8px;">Precision sub-assemblies, cartridges with carbide inserts &amp; WCC coated actuators</p>


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production


## ERGOgrip

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 7-14-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)  Brand: MBA	<b>1.004.0207</b>	<b>S687A-ZZS-MC3-GR3-T9H</b>	
<b>Replacement Bearing Set</b> Bearing Set with bearing, washer, and two o-rings Includes Bearing, Washer, and two O-Rings  Brand: MBA	<b>1.004.0203</b>	<b>BKU20203</b>	

## Intra 24CN - 24LN - 553 5761

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2) 2 Required	<b>220 0735</b>	<b>SMR74-ZS-T4HB-MC3-GR5-AF12</b>	


## Intra 24LS - 553 5641

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Open precision dental handpiece bearing with Torlon retainer (2.35x5.5x2) 2 Required	<b>2200185</b>	<b>SMR55-0235-MC3-TO-AF12</b>	

## Intra 24LZ - 24LN - 553 5641

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Precision Bearing Open Abec 3 precision bearing (2.35x5.5x2) 2 Required	<b>2200010</b>	<b>SMR55-0235-MC3-GR3-AF12</b>	

## Intra 25CH Contra Angle

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2) 2 Required	<b>2201305 / 1307</b>	<b>SMR74-ZS-T4HB-MC3-GR5-AF12</b>	

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Intra 25LN

Position / Description / Application

Manufacturers Part

MBA Part

**Stem Bearing**

**2200730 in 553 7201 Stem**

**SMR5-0235-T4HB-MC3-AF12**

High Speed Bearing - Torlon Retainer

Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)



**Replacement Bearing**

**220 1307 in 553 7341 Cartridge**

**SMR74-ZS-T4HB-MC3-GR5-AF12**

High Speed Bearing - Torlon Retainer

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)

2 Required



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## K10

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 3-8-4mm</b> High Speed Bearing - Polyamide Retainer Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	0.220.0778	S693A-ZZS-MC3-T9H	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	0.220.0779	SMR104A-ZZS-MC3-GR3-T9H	
<b>Bearing for Shaft 8-16-6mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)	0.220.0113	S688B-ZZ-MC3-GR3-T9H-G36	
<b>Replacement Chuck</b>  True running accuracy less than 0.03, bore diameter 2.35mm  Brand: MBA	0.674.7752	BKU67752	
<b>Replacement Spring Washer</b> Laboratory Drive  Brand: MBA	0.222.4018	BKU64018	
<b>Replacement Spring Washer</b> Laboratory Drive  Brand: MBA	0.222.4045	BKU64045	
<b>Replacement Lock Ring</b> Laboratory Drive  Brand: MBA	0.222.4126	BKU64126	
<b>Replacement Quad Ring</b> Laboratory Drive  Brand: MBA	0.224.2408	BKU62408	
<b>Replacement Rod</b> Laboratory Drive  Brand: MBA	0.674.1191	BKU61191	

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

<p><b>Replacement Cap</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.2861</b></p>	<p><b>BKU62861</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Screw nut</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.2871</b></p>	<p><b>BKU62871</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Roll Pin</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4152</b></p>	<p><b>BKU64152</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Distance Bush</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4162</b></p>	<p><b>BKU64162</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Screw</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4172</b></p>	<p><b>BKU64172</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Driver</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4192</b></p>	<p><b>BKU64192</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Bush</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.6302</b></p>	<p><b>BKU66302</b></p>	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
<p><b>Replacement Motor Brush</b></p> <p>Brand: MBA</p>	<p><b>0.674.2481</b></p>	<p><b>BKU62481</b></p>	 <p style="font-size: 8px; color: red;">Brushes, Electric Dental Motors &amp; Laboratory Drives</p>


<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## K11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 3-8-4mm</b> High Speed Bearing - Polyamide Retainer Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	<b>0.220.0778</b>	<b>S693A-ZZS-MC3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)  Brand: MBA	<b>0.220.0779</b>	<b>SMR104A-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x6)	<b>0.220.0044</b>	<b>SR11010A-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed polaymide retainer (7.94x15.88x7.5)	<b>0.220.0045</b>	<b>SR11010B-ZZS-MC3-GR7-T9H</b>	
<b>Replacement Roll Pin</b> Laboratory Drive  Brand: MBA	<b>0.674.4152</b>	<b>BKU64152</b>	
<b>Replacement Distance Bush</b> Laboratory Drive  Brand: MBA	<b>0.674.4162</b>	<b>BKU64162</b>	
<b>Replacement Screw</b> Laboratory Drive  Brand: MBA	<b>0.674.4172</b>	<b>BKU64172</b>	
<b>Replacement Driver</b> Laboratory Drive  Brand: MBA	<b>0.674.4192</b>	<b>BKU64192</b>	

## K12 PowerGrip®

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rear Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x6)	<b>0.220.0044</b>	<b>SR11010A-ZZS-MC3-GR3-T9H</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## K186

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Motor Brush**

**0.674.0241**

**BKU60241**

Brand: MBA



## K190

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Dental Motor**

**0.676.0871**

**BKU60871**

Brand: MBA



**Replacement Motor O-Ring**

**0.200.6120**

**BKU66120**

Special dental applications O-Ring

Brand: MBA




**Replacement Motor Brush**

**0.676.0881**

**BKU60881**

Brand: MBA



## K192

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Dental Motor**

**0.676.0971**

**BKU60971**

Brand: MBA



## K196

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Dental Motor**

**0.676.0871**

**BKU60871**



Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>




## K200

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Dental Motor</b>	<b>1.000.6859</b>	<b>BKU66859</b>	
Brand: MBA			
<b>Replacement Motor O-Ring</b>	<b>0.200.6120</b>	<b>BKU66120</b>	
Special dental applications O-Ring			
Brand: MBA			

## K201

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Motor Brush</b>	<b>0.676.0881</b>	<b>BKU60881</b>	
Brand: MBA			




## K4

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Rear Shaft Bearing 10x19x7</b>	<b>0.220.0056</b>	<b>S6800B-ZZS-MC3-GR3-KC</b>	
High Speed Bearing - Phenolic Retainer Double shielded Abec 3 stainless steel bearing with a high speed phenolic cage (10x19x7)			
Brand: MBA			
<b>Replacement Chuck</b>	<b>0.674.4871</b>	<b>BKU64871</b>	
True running accuracy less than 0.03, bore diameter 2.35mm			
Brand: MBA			
<b>Front Shaft Bearing Set</b>	<b>0.220.0058</b>	<b>BKU20058</b>	
Includes Bearing, Spacer, Felt Washer and 2 Orings Special double grooved OD bearing with o-rings, spacer, and washer			
Brand: MBA			

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## K5

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Rear Shaft Bearing 10x19x7</b> High Speed Bearing - Phenolic Retainer Double shielded Abec 3 stainless steel bearing with a high speed phenolic cage (10x19x7)  Brand: MBA	<b>0.220.0056</b>	<b>S6800B-ZZS-MC3-GR3-KC</b>	
<b>Replacement Chuck</b>  True running accuracy less than 0.03, bore diameter 2.35mm  Brand: MBA	<b>0.674.4871</b>	<b>BKU64871</b>	
<b>Front Shaft Bearing Set</b> Includes Bearing, Spacer, Felt Washer and 2 Orings Special double grooved OD bearing with o-rings, spacer, and washer  Brand: MBA	<b>0.220.0058</b>	<b>BKU20058</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## K9



Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)  Brand: Myonic	<b>220-0158</b>	<b>S692-ZZS-MC34-GR5-R-G48</b>	
<b>Bearing for Centre</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)	<b>220-0157 / 0779</b>	<b>S694-ZZ-MC3-JP6</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6) 2 Required	<b>220-0113</b>	<b>S688B-ZZ-MC3-GR3-T9H-G36</b>	
<b>Replacement Chuck</b>  True running accuracy less than 0.03, bore diameter 2.35mm  Brand: MBA	<b>0.674.7752</b>	<b>BKU67752</b>	
<b>Replacement motor</b> Laboratory Drive  Brand: MBA	<b>0.674.1261</b>	<b>BKU61261</b>	 Armature, Rotors, Electric Dental Motors & Laboratory Drives
<b>Replacement Spring Washer</b> Laboratory Drive  Brand: MBA	<b>0.222.4018</b>	<b>BKU64018</b>	 Spare Parts for Laboratory Drives
<b>Replacement Spring Washer</b> Laboratory Drive  Brand: MBA	<b>0.222.4045</b>	<b>BKU64045</b>	 Spare Parts for Laboratory Drives
<b>Replacement Lock Ring</b> Laboratory Drive  Brand: MBA	<b>0.222.4126</b>	<b>BKU64126</b>	 Spare Parts for Laboratory Drives
<b>Replacement Quad Ring</b> Laboratory Drive  Brand: MBA	<b>0.224.2408</b>	<b>BKU62408</b>	 Spare Parts for Laboratory Drives

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

<p><b>Replacement Roll Pin</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4152</b></p>	<p><b>BKU64152</b></p>	
<p><b>Replacement Distance Bush</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4162</b></p>	<p><b>BKU64162</b></p>	
<p><b>Replacement Screw</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4172</b></p>	<p><b>BKU64172</b></p>	
<p><b>Replacement Driver</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.4192</b></p>	<p><b>BKU64192</b></p>	
<p><b>Replacement Bush</b> Laboratory Drive</p> <p>Brand: MBA</p>	<p><b>0.674.6302</b></p>	<p><b>BKU66302</b></p>	
<p><b>Replacement Motor Brush</b></p> <p>Brand: MBA</p>	<p><b>0.674.0241</b></p>	<p><b>BKU60241</b></p>	

## K9 Blau




Position / Description / Application	Manufacturers Part	MBA Part	
<p><b>Bearing for Rotor</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (5x11x4)</p>	<p><b>0.220.0223</b></p>	<p><b>SMR115-ZZ-MC3-JP6</b></p>	
<p><b>Replacement Bearing</b> Includes Bearing, Washer, and two O-Rings Includes Bearing, Washer, and two O-Rings</p> <p>Brand: MBA</p>	<p><b>1.000.3586</b></p>	<p><b>BKU23586</b></p>	

Key




Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## K9 Blue

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rear Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x6)	<b>0.220.0044</b>	<b>SR11010A-ZZS-MC3-GR3-T9H</b>	
<b>Replacement Chuck</b> True running accuracy less than 0.03, bore diameter 2.35mm Brand: MBA	<b>0.674.4871</b>	<b>BKU64871</b>	
<b>Replacement Dental Motor</b> Brand: MBA	<b>0.674.4811</b>	<b>BKU61261</b>	Armature, Rotors, Electric Dental Motors & Laboratory Drives 




## K9 EWL 950

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Rod</b> Laboratory Drive Brand: MBA	<b>0.674.1191</b>	<b>BKU61191</b>	Spare Parts for Laboratory Drives 
<b>Replacement Cap</b> Laboratory Drive Brand: MBA	<b>0.674.2861</b>	<b>BKU62861</b>	Spare Parts for Laboratory Drives 
<b>Replacement Screw nut</b> Laboratory Drive Brand: MBA	<b>0.674.2871</b>	<b>BKU62871</b>	Spare Parts for Laboratory Drives 


Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## K9 Umbra

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)  Brand: Myonic	<b>0.220.0158</b>	<b>S692-ZZS-MC34-GR5-R-G48</b>	
<b>Bearing for Rotor</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	<b>0.220.0157</b>	<b>SMR104A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft 8-16-6mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)	<b>0.220.0113</b>	<b>S688B-ZZ-MC3-GR3-T9H-G36</b>	



## K9 Umbra EWL 955

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)  Brand: ECO	<b>0.220.0148</b>	<b>SMR148A-ZZS-MC3-GR3-T9H</b>	

## KL702

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)  Brand: MBA	<b>1.003.5818</b>	<b>S693-ZZS-MC34-TH5-R-L23</b>	
<b>Bearing for Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>1.003.6023</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	



## KL703 LED

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	<b>1.007.0966</b>	<b>S683A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>1.003.6023</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	





<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## L80

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>1.009.1782</b>	<b>BKU61782</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>=L890:L891+U890:U891=C874</b>	<b>BKU63606</b>	



## M05L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Precision ceramic-hybrid dental handpiece bearing with Torlon retainer (3.18x5.79x2270)  Brand: MBA	<b>1.010.9766</b>	<b>SMR5.79-3.175-CSN-ZZS-MC3-TO</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.5x6.5x2)  Brand: MBA	<b>1.010.9767</b>	<b>SMR6.50-3.50A-CSN-ZZS-MC3-TO</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>1.011.1374</b>	<b>BKU61374</b>	
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>1.011.3026</b>	<b>BKU63206</b>	

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## M25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	1.011.0129 + 1.001.3515	BKU62915	
<b>Drive Kit</b>  Brand: MBA	1.011.0129 + 1.001.3515 + 1.001.3519	BKU62919	
<b>Replacement Head Housing</b>  Brand: MBA	1.009.3635	BKU63635	
<b>Replacement Shaft</b>  Brand: MBA	1.001.3518	BKU63518	
<b>Replacement Sprayinset</b>  Brand: MBA	1.009.3629	BKU63629	




Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## M25L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Precision ceramic-hybrid dental handpiece bearing with Torlon retainer (3.18x5.79x2270)  Brand: MBA	<b>1.010.9766</b>	<b>SMR5.79-3.175-CSN-ZZS-MC3-TO</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.5x6.5x2)  Brand: MBA	<b>1.010.9767</b>	<b>SMR6.50-3.50A-CSN-ZZS-MC3-TO</b>	
<b>Replacement Dog Gear 1:5</b> Gear  Brand: MBA	<b>1.001.3071</b>	<b>BKU63071</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>1.011.0129</b>	<b>BKU60129</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.009.3634 + 0.201.8157 + 1.009.3638</b>	<b>BKU63634</b>	

## M8700L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>3.001.9915</b>	<b>BKURC210020</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>3.001.9915</b>	<b>BKURC2100203</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>3.001.9915</b>	<b>BKURC2100205</b>	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

## M8700L Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Application O-Ring  Brand: MBA	<b>3.002.0256</b>	<b>BKU60256</b>	
<b>M8900L</b> Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>3.000.1874</b>	<b>BKU61874</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2.000.2266</b>	<b>BKUR210019</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2.000.2266</b>	<b>BKURC210019</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2.000.2266</b>	<b>BKUR2100195</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKUR2100193</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKUR2100194</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKURC2100193</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2.000.2266</b>	<b>BKURC2100195</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production








## M9000L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>3.000.1874</b>	<b>BKU61874</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2.000.2266</b>	<b>BKUR210019</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2.000.2266</b>	<b>BKURC210019</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKUR2100193</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2.000.2266</b>	<b>BKUR2100195</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKUR2100194</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2.000.2266</b>	<b>BKURC2100193</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2.000.2266</b>	<b>BKURC2100195</b>	
<b>M9000L Air Driven Handpiece</b>			
<b>Replacement O-Ring</b>  Brand: MBA	<b>3.000.2581</b>	<b>BKU62581</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## Magno 634

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic	220-0128 / 0145	SR144-D106-ZS-TO	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-KC	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA	220-0128 / 0145	SR144-D106-CSN-ZS-TO	 <b>Image Coming Soon</b>
<b>Miscellaneous Refer OEM Part Number</b>			
<b>Bearing for Rotor Cartridge</b> High Precision Bearing Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer Brand: myonic	0.220.0128	SR144-D106-ZS-KC	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) Brand: MBA	0.220.0128	SR144-D106-CSN-ZS-TO	 <b>Image Coming Soon</b>
<b>Bearing for Rotor Cartridge - Alternative</b> High Precision Bearing Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer Brand: myonic	0.220.0128	SR144-D106-ZS-TO	
<b>Bearing for Rotor Cartridge - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) Brand: MBA	0.220.0128	SR144-D106-CSN-ZS-KC	 <b>Image Coming Soon</b>



Key Usually stocked or Short Delivery \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## Multiflex® Coupling

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Small - Silicone Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required	<b>200-6113 (Part of Set)</b>	<b>OR-00370-099-S60-R</b>	
<b>Replacement O-Ring</b> Large Special sized o-ring used in dental Applications 3 Required	<b>200-6113 (Part of Set)</b>	<b>OR-00640-080-N70</b>	


## Multiflex® Coupling Motor K190

Position / Description / Application	Manufacturers Part	MBA Part	
<b>LED Bulb</b>  Brand: MBA	<b>0.553.3881</b>	<b>BKU63881LED</b>	
<b>Xenon Lamp</b>  Brand: MBA	<b>0.553.3881</b>	<b>BKU63881Xenon</b>	

## Multiflex® Coupling Motor K201

Position / Description / Application	Manufacturers Part	MBA Part	
<b>LED Bulb</b>  Brand: MBA	<b>0.553.3881</b>	<b>BKU63881LED</b>	
<b>Xenon Lamp</b>  Brand: MBA	<b>0.553.3881</b>	<b>BKU63881Xenon</b>	

## S201 C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>1.001.4490 + 1.001.4501</b>	<b>BKU69001</b>	




Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*




## S201XL

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 20:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>3.000.5177</b>	<b>BKU65177</b>	
<b>Replacement Head Housing</b>  Brand: MBA	<b>3.000.5315</b>	<b>BKU65315</b>	

## S459C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.001.2518</b>	<b>BKURC210022</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.001.2518</b>	<b>BKURC2100225</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.2518</b>	<b>BKURC2100223</b>	

## S459L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.001.2518</b>	<b>BKURC210022</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.001.2518</b>	<b>BKURC2100225</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.001.2518</b>	<b>BKURC2100223</b>	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## S550

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**1.000.8101**

**S693-CSN-ZBS-MC34-GR7-L23**

High Precision Bearing

Double shielded with special bronze shield one side. Abec 7 precision stainless ceramic-hybrid bearing (3x8x3)

Brand: MBA



## S605

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**1.000.7979**

**SWR144-D28W-CSN-Z-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA



**Bearing for Rotor Cartridge**

**1.000.7979**

**SWR144-D28W-CSN-Z-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)

Brand: MBA



**Bearing for Rotor Cartridge**

**SG1304X**

**SWR144-D28W-Z-KC**

High Speed Bearing - Phenolic Retainer

Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA



**Replacement Press Button**

**1.007.7277**

**BKU67277**

Brand: MBA



**Replacement Spring Washer**

**1.001.5032**

**BKU65032**




Brand: MBA





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## S605C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.008.4164</b>	<b>BKURC210021</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.008.4164</b>	<b>BKURC2100215</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.008.4164</b>	<b>BKURC2100213</b>	

## S609

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.0854</b>	<b>BKU60854</b>	
<b>Replacement Spring Washer</b>  Brand: MBA	<b>1.001.5032</b>	<b>BKU65032</b>	

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>





## S609C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.008.4163</b>	<b>BKUR210016</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.008.4163</b>	<b>BKUR210016S</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.008.4163</b>	<b>BKUR2100163</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.008.4163</b>	<b>BKUR210016T</b>	 <p style="font-size: small; text-align: right;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.008.4163</b>	<b>BKUR2100164</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly without ball bearings  Brand: myonic	<b>1.008.4163</b>	<b>BKURC210016</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## S615L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>1.007.7277</b>	<b>BKU67277</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>1.008.4164</b>	<b>BKURC210021</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>1.008.4164</b>	<b>BKURC2100215</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>1.008.4164</b>	<b>BKURC2100213</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## S619C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	1.008.4163	BKUR210016	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.008.4163	BKUR2100163	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.008.4163	BKUR2100164	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	1.008.4163	BKUR2100165	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	1.008.4163	BKUR210016T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly without ball bearings  Brand: myonic	1.008.4163	BKURC210016	

## S619L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	1.007.0854	BKU60854	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## SF4004

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft 7-14-5mm**

**0.220.0112**

**S687A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)

Brand: MBA



**Bearing for Shaft**

**0.247.6056**

**S608-ZZS-MC3-GR3-KC**

High Speed Bearing - Phenolic Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (8x22x7)

Brand: MBA



**Replacement Chuck**

**0.671.1562**

**BKU61562**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## SF4005

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**0.247.6057**

**S608-ZZS-MC3-GR3-KC**

High Speed Bearing - Phenolic Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (8x22x7)

Brand: MBA



**Replacement Chuck**

**0.671.1562**

**BKU61562**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## SONICflex 2000 Airscaler

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Bush**

**0.571.5032**

**BKU65032S**

Refer OEM Part Number

Brand: MBA



## SONICflex 2003 Airscaler

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Bush**

**0.571.5032**

**BKU65032S**

Refer OEM Part Number



Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Surgical Instrument 3554

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Lower Bearing</b> 4x9x4 Stainless Steel Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>220 0103</b>	<b>S684A-ZZ-MC3-JP6</b>	
<b>Bearing for Centre</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>0.220.1302</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Colibri

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

# Kool Torque®

## USA

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**BB 6-13-5**

**S686B-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)

2 Required

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 330E

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## 557 CX ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


10632

BKUR360001

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## 557 Push Button

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

SR144C-D26-KC

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x1.91)



### Bearing for Chuck End

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)



### Replacement O-Ring

OR-00600-075-N70

Special Dental Application O-Ring

Special sized o-ring used in dental Applications

2 Required



## 557 Turbo ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


10633

BKUR360003

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## 557SR

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 757

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x1.91)

SR144C-D26-KC



### Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-ZS-KC



### Replacement O-Ring

Special Dental Application O-Ring

Special sized o-ring used in dental Applications

2 Required

OR-00600-075-N70



## 757 CX ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


10629

BKUR360002

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## 757 Multi Port

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


10630

BKUR360004

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## 757 Single Port

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


10631

BKUR360005

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



Key




Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Workhorse

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x1.91)		<b>SR144C-D26-KC</b>	 <b>Image Coming Soon</b>
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)		<b>SR144-D28-ZS-KC</b>	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		<b>OR-00600-075-N70</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## BM50S1

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**SUS420J2**

**BKU64202**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## H35LSP

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Motor Brush**

**BKU60108**

Brand: MBA



## H35SP1

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**BKU61107**

Brand: MBA



## H37L

Position / Description / Application


Manufacturers Part

MBA Part

**Replacement Motor Brush**

**BKU60108**

Brand: MBA



## Handy 700

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**SUS420J2**

**BKU64202**





True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA





# Marathon®


## Motor DC Type

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	<b>830ZZ</b>	<b>S693A-ZZS-MC3-T9H</b>	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded bearing with a high speed polyamide retainer and light oil lubricant (3x10x4)  Brand: MBA	<b>1030ZZ</b>	<b>S623-ZZS-MC3-T9H-AF12</b>	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>1360ZZ</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>188ZZ</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	

## Motor E Type

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	<b>830ZZ</b>	<b>S693A-ZZS-MC3-T9H</b>	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded bearing with a high speed polyamide retainer and light oil lubricant (3x10x4)  Brand: MBA	<b>1030ZZ</b>	<b>S623-ZZS-MC3-T9H-AF12</b>	

## Motor ETC Type

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>1360ZZ</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	

## Multi 600

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Chuck</b> True running accuracy less than 0.03, bore diameter 2.35mm  Brand: MBA	<b>SUS420J2</b>	<b>BKU64202</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Marathon®

## SDEH37L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61107

Brand: MBA



## SDESH40C

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61108

Brand: MBA



## SH40C

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60108

Brand: MBA



## SH45C

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61108


Brand: MBA



Replacement Motor Brush

BKU60108

Brand: MBA



## Spindle DC Type

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Shaft

1260ZZ

SMR126A-ZZS-MC3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded high speed stainless steel bearing with Polyamide retainer (6x12x4)



Bearing for Shaft

1480ZZ

SMR148A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)



Brand: ECO



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Spindle ETC Type

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (6x12x4)	<b>1260ZZ</b>	<b>SMR126A-ZZS-MC3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)  Brand: ECO	<b>1480ZZ</b>	<b>SMR148A-ZZS-MC3-GR3-T9H</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Micro Mega®

**350**

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**SFR144A-D27-ZSO-KC**

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

Brand: Myonic



Key



Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 1000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D28-CSN-KC</b>	

## 400

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-KC</b>	
<b>Bearing for Cap End</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 5 Hole Hose

Position / Description / Application

Manufacturers Part

MBA Part

**Xenon Lamp**

**790257R**

**BKU60257Xenon**

Brand: MBA



## 5 Hole Hoses

Position / Description / Application

Manufacturers Part


MBA Part

**LED Bulb**

**790257R**

**BKU63307LED**

Brand: MBA



## 6 Pin Coupling

Position / Description / Application

Manufacturers Part

MBA Part

**LED Bulb**

**790257R**

**BKU63307LED**

Brand: MBA




**Xenon Lamp**

**790257R**

**BKU60257Xenon**

Brand: MBA



## 8000 - 8000i

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA

**SR144-D28-KC**



**Replacement Bearing**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA

**SR144-D28-CSN-KC**



## Bright ISO

Position / Description / Application

Manufacturers Part

MBA Part

**Xenon Lamp**

**262522**

**BKU62522Xenon**

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## C6

Position / Description / Application

Manufacturers Part

MBA Part

**Xenon Lamp**

**262522**

**BKU62522Xenon**

Brand: MBA



## Easycare Mini

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**SG1223X.2mq**

**SR145-D75-ZZE-KC**

High Speed Bearing - Phenolic Retainer Shielded Option

Double shielded Abec 7 precision bearing with stepped OD and high speed Phenolic retainer

Brand: MBA



**Bearing for Rotor Cartridge**

**SG1223X-D**

**SR145-D75-CSN-ZZE-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Shielded Option

Double shielded precision special ceramic-hybrid dental bearing with stepped OD and high speed Phenolic retainer

Brand: MBA



## eStylus Motor

Position / Description / Application

Manufacturers Part


MBA Part

**LED Bulb**

**774150**

**BKU66016LED**

Brand: MBA



## Quiet Air

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Cap End**

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA

**SR144-D28-KC**



**Bearing for Cap End**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA

**SR144-D28-CSN-KC**



**Bearing for Chuck End**

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

Brand: Myonic

**SFR144A-D27-ZSO-KC**




**Bearing for Chuck End**

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

Brand: Myonic





**SFR144A-D27-ZSO-TO**





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Quiet Air Compact

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-KC</b>	
<b>Bearing for Cap End</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	





## Stylus

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Shielded Option Double shielded Abec 7 precision bearing with stepped OD and high speed Phenolic retainer  Brand: MBA	<b>SG1223X.2mq</b>	<b>SR145-D75-ZZE-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Shielded Option Double shielded precision special ceramic-hybrid dental bearing with stepped OD and high speed Phenolic retainer  Brand: MBA	<b>SG1223X-D</b>	<b>SR145-D75-CSN-ZZE-KC</b>	




Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*


## Tradition

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	

## Tradition Lever

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece ceramic-hybrid bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	<b>SG1222X-D</b>	<b>SR144-D28-CSN-ZS-KC</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>780118</b>	<b>BKURMU260101</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>780118</b>	<b>BKURM260101</b>	

## Tradition Lever Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> Daily Lube - High Speed Bearing - Phenolic Retainer Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>SG1301X-D</b>	<b>SFR144A-D1301-CSN-ZSO-KC</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## Tradition XGT

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


790118

BKURMU260100

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



### Turbine Sub Assembly Without Bearings

790118

BKURM260100

For Handpiece with Wedge Type Chuck

Brand: MBA



## Traditon Lever

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge

SG1222X

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)



## Traditon Lever Push Button

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge

SG1301X

SFR144A-D1301-ZSO-KC

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)

Brand: MBA



## Tru-Line

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA

SR144-D28-KC



### Bearing for Cap End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

Brand: MBA

SR144-D28-CSN-KC



### Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

Brand: Myonic

SFR144A-D27-ZSO-KC




### Bearing for Chuck End

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

Brand: Myonic

SFR144A-D27-ZSO-TO



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## XGT Lever Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement O-Ring**

**BKU69998**

Special Dental Application O-Ring

Brand: MBA





Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# Miscellaneous



## Dental Washers

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13)		W0030-F-007-001-PT	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## EM1190K

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>BE5015</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	
<b>Replacement Motor Brush</b>  Brand: MBA	<b>SP1219</b>	<b>BKU60881</b>	

## HB21

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC1034</b>	<b>BKUR280108</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC1034</b>	<b>BKUR2801085</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC1034</b>	<b>BKUR2801083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2034</b>	<b>BKURC2800033</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2034</b>	<b>BKURC2800035</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2034</b>	<b>BKURC280003</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## HC20

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC2038	BKUR280120	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2038	BKURC2800075	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2038	BKURC2800073	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2038	BKUR2801205	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2038	BKUR2801203	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2038	BKURC280007	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HC21

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC2038</b>	<b>BKUR280120</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2038</b>	<b>BKURC2800073</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2038</b>	<b>BKURC2800075</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2038</b>	<b>BKUR2801205</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2038</b>	<b>BKUR2801203</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2038</b>	<b>BKURC280007</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HC4021

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1031 / TC2031	BKUR280106	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801063	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801064	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1031 / TC2031	BKUR2801065	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2031	BKURC2800015	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2031	BKURC2800013	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2031	BKURC280001	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HC4022

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1032 / TC2032	BKUR280107	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801074	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1032 / TC2032	BKUR2801075	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801073	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2032	BKURC2800025	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2032	BKURC2800023	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2032	BKURC280002	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

# MK-Dent®

## HC5021

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1031 / TC2031	BKUR280106	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801064	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR280106T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1031 / TC2031	BKUR2801065	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801063	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2031	BKURC2800013	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2031	BKURC2800015	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2031	BKURC280001	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HC5022

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1032 / TC2032	BKUR280107	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801074	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1032 / TC2032	BKUR2801075	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR280107T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801073	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2032	BKURC2800025	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2032	BKURC2800023	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2032	BKURC280002	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## HC60012

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1032 / TC2032	BKUR280107	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801073	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR280107T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1032 / TC2032	BKUR2801075	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1032 / TC2032	BKUR2801074	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2032	BKURC2800023	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2032	BKURC2800025	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2032	BKURC280002	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HC6021

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	TC1031 / TC2031	BKUR280106	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801063	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR2801064	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC1031 / TC2031	BKUR2801065	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC1031 / TC2031	BKUR280106T	 <p>Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	TC2031	BKURC2800013	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	TC2031	BKURC2800015	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	TC2031	BKURC280001	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## HC8022

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**TC2042**

**BKURC2800063**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



### Replacement Turbine With Bearings

**TC2042**

**BKURC2800065**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




### Turbine Sub Assembly Without Bearings

**TC2042**

**BKURC280006**

For Handpiece with Wedge Type Chuck

Brand: MBA



## HC9021

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine Excluding Bearings


**TC2041**

**BKUR280110**

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA



### Replacement Turbine With Bearings


**TC2041**

**BKUR2801103**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



### Replacement Turbine With Bearings


**TC2041**

**BKURC2800055**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



### Replacement Turbine With Bearings


**TC2041**

**BKURC2800053**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



### Replacement Turbine With Bearings


**TC2041**

**BKUR2801105**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




### Turbine Sub Assembly Without Bearings

**TC2041**

**BKURC280005**

For Handpiece with Wedge Type Chuck

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## HE17

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC1034</b>	<b>BKUR280108</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC1034</b>	<b>BKUR2801085</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC1034</b>	<b>BKUR2801083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2034</b>	<b>BKURC2800033</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2034</b>	<b>BKURC2800035</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2034</b>	<b>BKURC280003</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HE18

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC1035</b>	<b>BKUR280109</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC1035</b>	<b>BKUR2801093</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC1035</b>	<b>BKUR2801095</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2035</b>	<b>BKURC2800045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2035</b>	<b>BKURC2800043</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2035</b>	<b>BKURC280004</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## HE22

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC1035</b>	<b>BKUR280109</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC1035</b>	<b>BKUR2801095</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC1035</b>	<b>BKUR2801093</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2035</b>	<b>BKURC2800043</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2035</b>	<b>BKURC2800045</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2035</b>	<b>BKURC280004</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last




Long Delivery or Large Minimum Quantities  
*Not in current production*

# MK-Dent®

## HP21

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TC2041</b>	<b>BKUR280110</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2041</b>	<b>BKUR2801103</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2041</b>	<b>BKUR2801105</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2041</b>	<b>BKURC2800053</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2041</b>	<b>BKURC2800055</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2041</b>	<b>BKURC280005</b>	

## HP22

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TC2042</b>	<b>BKURC2800065</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TC2042</b>	<b>BKURC2800063</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TC2042</b>	<b>BKURC280006</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## LC01L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**BE5014**

**S684A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x9x4)



**Bearing for Shaft**

**BE5015**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



## LC11L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**


**BE5013**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## LC27L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**BE5018**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x10x4)



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities




\*\* Available while stocks last

*Not in current production*

## LCH11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	<b>BE5011</b>	<b>SWA683B-ZZ-C69-GR5-TO</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>BE5012</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Head Gear Complete Dental Head 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>SP2005</b>	<b>BKU62005</b>	
<b>Intermediate Piece Dental Head 1:1</b>  Brand: MBA	<b>SP2024</b>	<b>BKU62024</b>	
<b>Drive Kit</b>  Brand: MBA	<b>SP2005 + SP2024</b>	<b>BKU60524</b>	

## LCH21

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>SP2105 + SP2124</b>	<b>BKU60254</b>	
<b>Head Gear Complete Dental Head 2:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>SP2105</b>	<b>BKU62105</b>	
<b>Intermediate Piece Dental Head 2:1</b>  Brand: MBA	<b>SP2124</b>	<b>BKU62124</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## LCH21P

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Headgear

**BE6055**

**SMR63-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



## LE11

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**BE5013**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## LE148L

Position / Description / Application

Manufacturers Part

MBA Part


### Bearing for Shaft

**BE5019**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



### Bearing for Shaft

**BE5017**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## LE15L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Middle Gear

**BE5027**

**SMR4.762-2.35A-ZZS-MC3-GR7-TO**

High Speed Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)

Brand: MBA



### Bearing for Middle Gear

**BE5027**

**SR133A-ZZS-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)



### Bearing for Headgear

**BE5029**

**SWR144C-D38W-CSN-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open, stainless steel bearing, with extended inner ring and very narrow width, for dental applications. (3.18x6.35x1.9)

Brand: MBA



### Bearing for Headgear

**BE5021B**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## LE21

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**


**BE5013**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## LE27L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**BE5018**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



**Bearing for Shaft**

**BE5016**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## LE74L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**BE5019**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



**Bearing for Shaft**

**BE5017**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## LP01L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**BE5014**

**S684A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x9x4)



**Bearing for Shaft**


**BE5015**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## LP11L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft


**BE5013**

**S693-ZZS-MC34-TH5-R-L23**

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## LP15L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Middle Gear

**BE5027**

**SMR4.762-2.35A-ZZS-MC3-GR7-TO**

High Speed Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)

Brand: MBA




### Bearing for Middle Gear

**BE5027**

**SR133A-ZZS-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)



### Bearing for Headgear

**BE5029**

**SWR144C-D38W-CSN-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open, stainless steel bearing, with extended inner ring and very narrow width, for dental applications. (3.18x6.35x1.9)

Brand: MBA



### Bearing for Headgear


**BE5021B**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



## LP54L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**BE5018**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



### Bearing for Shaft

**BE5016**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



## LP74L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**BE5016**

**SMR104A-ZZS-MC3-GR7-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## LPH11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	<b>BE5011</b>	<b>SWA683B-ZZ-C69-GR5-TO</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>BE5012</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Head Gear Complete Dental Head 1:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>SP2005</b>	<b>BKU62005</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Intermediate Piece Dental Head 1:1</b>  Brand: MBA	<b>SP2024</b>	<b>BKU62024</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Drive Kit</b>  Brand: MBA	<b>SP2005 + SP2024</b>	<b>BKU60524</b>	 <p>Drive Kits for Electric Driven Handpieces and Heads</p>



Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*


## LPH21

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	<b>BE5011</b>	<b>SWA683B-ZZ-C69-GR5-TO</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>BE5012</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Drive Kit</b>  Brand: MBA	<b>SP2105 + SP2124</b>	<b>BKU60254</b>	
<b>Head Gear Complete Dental Head 2:1</b> For Handpiece with Latch Type Chuck  Brand: MBA	<b>SP2105</b>	<b>BKU62105</b>	
<b>Intermediate Piece Dental Head 2:1</b>  Brand: MBA	<b>SP2124</b>	<b>BKU62124</b>	

## LS2115L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (2.38x4.76x1588)	<b>BE5022</b>	<b>SR133-MC3-TH7-TO</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	<b>BE5021A</b>	<b>SMR74-CSN-TH7-TO</b>	

## Miscellaneous Refer OEM Part Number

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)	<b>BE5010</b>	<b>SR1-5A-ZZS-MC3-JP6</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## QC6016K

Position / Description / Application

Manufacturers Part


MBA Part

LED Bulb and Cap

BU8012

BKU68012LED

Brand: MBA



## QC6016KW

Position / Description / Application

Manufacturers Part

MBA Part

LED Bulb and Cap

BU8012

BKU68012LED

Brand: MBA



Key



Usually stocked or Short Delivery

\*\* Available while stocks last



Long Delivery or Large Minimum Quantities

*Not in current production*


## Falcon

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Chuck End</b> High Precision Bearing Shielded Option Deleted Item. Available while stocks last. (6x10x3) 3 Required Brand: Myonic		<b>SMR106A-ZZS-MC34-P4P-R-AF12</b>	
<b>Bearing for Chuck End</b> High Precision Bearing Shielded Option Deleted Item. Available while stocks last. (6x10x3) 3 Required Brand: Myonic		<b>SMR106A-ZZS-MC34-GR7-R-G48</b>	

## Micromotor

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (2.5x8x4)	<b>BB 2.5-8-4</b>	<b>S602XA-ZZ-MC3-JP6</b>	
<b>Replacement Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (3x9x4)	<b>BB 3-9-4</b>	<b>SMR93A-ZZ-MC3-JP6</b>	

## PAR4HEXO

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)  Brand: MBA	<b>418NA-CER</b>	<b>SR144-D28-CSN-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	<b>SG1357X-D</b>	<b>SR144-D28-CSN-ZE-TO</b>	

## PAR4HUEXO

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)  Brand: MBA	<b>418NA-CER</b>	<b>SR144-D28-CSN-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)  Brand: MBA	<b>SG1357X-D</b>	<b>SR144-D28-CSN-ZE-TO</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Torx 1000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

Stainless Bearing Shielded Option

Special anti-corrosion stainless TBN-6H bearing designed primarily for fishing reels and salt-water applications (5x8x2.5)

2 Required

**SACB-MR85A-ZZ-MC3**



### Replacement Bearing

Stainless Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (5x8x2.5)

2 Required

**SMR85A-ZZS-MC3-JP6**



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

# MTC

**105**

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**SR144-D28-CSN-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA



**Replacement Bearing**

**SR144-D28-KC**

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA



Key






Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Lynx

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic		SR144A-D105-ZS-TO	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic		SR144A-D105-ZS-KC	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	0.220.2007, 220-0793K	SR144A-D107-CSN-ZS-TO	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-CSN-ZS-TO	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	220-0125	SR144A-D105-CSN-ZS-KC	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Micromat-Athlet Type 200 Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Chuck End**

**0229**

**\*\* S687A-ZZ-MC3-2OR**

Stainless Bearing

Deleted Item. Available while stocks last. Special bearing with 2 o-rings in OD grooves. (7x14x5)



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 3 405

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft 6-13-5mm

**B11223-1306**

**S686B-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)

Brand: MBA



### Bearing for Shaft

**B20215-2763**

**SR188A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)

Brand: MBA



## 450L

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NCH-SU03**

**BKUR250059**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



## 455L Phatelus 45S

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NCH-TU03**

**BKUR250060**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



## A200 N

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-MU03**

**BKUR250052**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## A300 N

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-SU03**

**BKUR250053**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## A400 N

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-TU03**

**BKUR250054**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## A500

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**Ti-MU03**

**BKUR250045**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## A600

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**Ti-SU03**

**BKUR250046**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## A700

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**Ti-TU03**

**BKUR250047**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



Key







Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Alpha Air

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Alpha Light

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic	<b>418WIF</b>	<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic	<b>418WIF</b>	<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIF-CER</b>	<b>SFR144A-D27-CSN-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Polyamide Retainer Precision flanged dental handpiece ceramic-hybrid bearing with polyamide retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIF-CER</b>	<b>SFR144A-D27-CSN-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Replacement O-Ring

OR-00607-079-N70

High Temperature  
Deleted Item. Available while stocks last. Special o-ring for dental applications  
2 Required



## BB Y

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

C037015 + C032012

BKU61512

Brand: MBA



### Head Gear Complete 1:1

C037015

BKU67015

Brand: MBA



### Replacement Middle Gear Complete 1:1

C032012

BKU62012

Brand: MBA



## CHQD

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-ZS-KC



### Bearing for Chuck End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-CSN-KC

Brand: MPB



### Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-KC

Brand: MPB



### Bearing for Chuck End

High Speed Bearings - Polyamide Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-ZSO

Brand: MPD



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

Not in current production

## Coupling

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement O-Ring

Small  
Special sized o-ring used in dental Applications  
2 Required

*OR-00460-070-N70*



### Replacement O-Ring

Large  
Special sized o-ring used in dental Applications  
4 Required

*OR-00635-089-N70*



## Coupling KCL

Position / Description / Application

Manufacturers Part


MBA Part

### LED Bulb

Brand: MBA

10000198

BKU60198LED



## Coupling KCL SCL

Position / Description / Application

Manufacturers Part


MBA Part

### LED Bulb

Brand: MBA

10000198

BKU60198LED



## DYNAled

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic

SX-SU03

BKUR250073



## E4R Slow Speed - Contra Angle Green Band

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

Stainless Bearing  
Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)

SMR63-MC3-JP6-AF12



### Bearing for Chuck End

Stainless Bearing  
Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)

SMR105-MC3-JP6-AF12



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## EX100

Position / Description / Application

Manufacturers Part

MBA Part


**Bearing for Shaft**

**B22227-7923**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



## EX103

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B22227-7923**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



## EX12

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B10321-1005**

**SMR105-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



## EX200

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B22227-7923**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



## EX5

Position / Description / Application

Manufacturers Part

MBA Part


**Bearing for Shaft**

**B13213-0804**

**SMR84A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x8x3)



**Bearing for Shaft**

**B10321-1005**

**SMR105-MC3-JP6-AF12**

High Precision Bearing



Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



## EX6

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	<b>B13213-0804</b>	<b>SMR84A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	<b>B10321-1005</b>	<b>SMR105-MC3-JP6-AF12</b>	



## EX6 Blue Band

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (4x8x3)	<b>BB 4-8-3</b>	<b>SMR84A-ZZ-MC3-JP6</b>	
<b>Replacement Bearing</b> Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	<b>BB 5-10-3</b>	<b>SMR105-MC3-JP6-AF12</b>	

## EX7

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	<b>B13213-0804</b>	<b>SMR84A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	<b>B10321-1005</b>	<b>SMR105-MC3-JP6-AF12</b>	

## EX8

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	<b>B13213-0804</b>	<b>SMR84A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	<b>B10321-1005</b>	<b>SMR105-MC3-JP6-AF12</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## EX9

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**B13213-0804**

**SMR84A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x8x3)



### Bearing for Shaft

**B10321-1005**

**SMR105-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



## FPB

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

**C875050**

**BKU65050**

Brand: MBA



## FPB Y

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**C875060 + C875012**

**BKU66012**

Brand: MBA



### Replacement Middle Gear Complete 1:1

**C875012**

**BKU65012**

Brand: MBA



## FPB03

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete 1:1

**C875060**

**BKU65060**

Brand: MBA



## FX15

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**C1052700 + C10000946**

**BKU60046**




Brand: MBA







Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## FX23

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b>	<b>C1052700</b>	<b>BKU62700</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 1:1</b>	<b>10000946</b>	<b>BKU60946</b>	
Brand: MBA			
<b>Replacement Press Button</b>	<b>C1052500 - FX25-04</b>	<b>BKU62500X</b>	
Brand: MBA			




## FX25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>C1052700 + C10000946</b>	<b>BKU60046</b>	
Brand: MBA			
<b>Head Gear Complete 1:1</b>	<b>C1056700</b>	<b>BKU62700</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 1:1</b>	<b>10000946</b>	<b>BKU60946</b>	
Brand: MBA			
<b>Replacement Press Button</b>	<b>C1052500 - FX25-04</b>	<b>BKU62500X</b>	
Brand: MBA			


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## FX57

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement End Cap 4:1</b>	<b>20002065</b>	<b>BKU62065</b>	
Brand: MBA			
<b>Replacement End Cap 4:1</b>	<b>20002066</b>	<b>BKU62066</b>	
Brand: MBA			
<b>Head Gear Complete 4:1</b>	<b>10001571</b>	<b>BKU61571</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 4:1</b>	<b>10000996</b>	<b>BKU60996</b>	
Brand: MBA			


## IS205

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b>	<b>B22227-7923</b>	<b>SR1-5A-ZZS-MC3-JP6</b>	
High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)			

## IS65

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b>	<b>B13213-0804</b>	<b>SMR84A-ZZ-MC3-JP6</b>	
High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)			
<b>Bearing for Shaft</b>	<b>B10321-1005</b>	<b>SMR105-MC3-JP6-AF12</b>	
High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)			

## M25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>	<b>C601050</b>	<b>BKU61050</b>	
Brand: MBA			

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## M25L New

Position / Description / Application

Manufacturers Part


MBA Part

**Head Gear Complete 1:1**

**C1038015**

**BKU68015**

Brand: MBA



**Replacement Middle Gear Complete 1:1**

**C1142022**

**BKU62022**

Brand: MBA



## M25L Old

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Middle Gear Complete 1:1**

**10000574**

**BKU60574**

Brand: MBA



## M500

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**P1001500 - SX-MU04**

**BKU61500**

Brand: MBA



## M500 Pana Max M Plus

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**SX-MU03**

**BKUR250048**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## M600

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**P1002500 - SX-SU04**

**BKU62500M**

Brand: MBA



## M65

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**H190180**

**BKU60180**

True running accuracy less than 0.03

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities




\*\* Available while stocks last

*Not in current production*

## M95

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	<b>740-CER</b>	<b>SMR74-CSN-TH7-TO</b>	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427062</b>	<b>BKU67062</b>	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427168</b>	<b>BKU67168</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>C600015</b>	<b>BKU60015</b>	
<b>Replacement Press Button</b>   Brand: MBA	<b>10000664</b>	<b>BKU60664</b>	

## M95L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Dog Gear 1:5 Dog ClutchC600C013</b> Clutch  Brand: MBA		<b>BKU69993</b>	
<b>Drive Gear Assembly - Dog Clutch Complete 1:5</b> Complete Assembly  Brand: MBA	<b>10000560</b>	<b>BKU60560</b>	
<b>Replacement Press Button</b>   Brand: MBA	<b>10000664</b>	<b>BKU60664</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## M95L New

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Middle Gear Complete 1:5**

**C1141022**

**BKU61022N**

Brand: MBA



## M95L Old

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Middle Gear Complete 1:5**

**10000559**

**BKU60559**

Brand: MBA



## Mach - Mach-LITE P045SU03 - P046TU03 Green-Black

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA

**SR144-D28-KC**



**Replacement Bearing**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA

**SR144-D28-CSN-KC**



**Replacement Wave Washer**

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)

**W0048-T-007-008-W3-S2**



**Replacement O-Ring**

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

Brand: NSK

**OR-00620-100-N70**



## Mach Lite 65 T

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**NCH-TU03**

**BKUR250060**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Mach Lite XT M

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-MU03**

**BKUR250052**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## Mach Lite XT S

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-SU03**

**BKUR250053**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## Mach Lite XT T

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**NMC-TU03**

**BKUR250054**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## Mach Mini- Mach-LITE Mini P047MU03 Red Band

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

Brand: MPB

**SFR144-D28-KC**



### Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

Brand: MPB

**SFR144-D28-CSN-KC**



### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

Brand: MPD

**SFR144-D28-ZSO**



### Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

Brand: NSK





**OR-00620-100-N70**




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## Mach N75

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-KC	


## Mach QD 45

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer  Brand: MBA	NCH-SU03	BKUR250059	 <p>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

## Motor KCL Coupler

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Xenon Lamp</b>  Brand: MBA	Y900132	BKU60132Xenon	 <p>Xenon Lamp for Dental Couplers &amp; Dental Motors</p>

## Motor Ultimate XL 450 UHR

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (4x9x4)  Brand: MBA	B11253-0904	S684A-ZZS-MC3-GR3-T9H	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*



## MP F

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>C521015 + C547022</b>	<b>BKU61522</b>	
Brand: MBA			
<b>Head Gear Complete 1:1</b>	<b>C521015</b>	<b>BKU61015M</b>	
Brand: myonic			

## MP Y

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b>	<b>C521015</b>	<b>BKU62101</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 1:1</b>	<b>C547022</b>	<b>BKU67022</b>	
Brand: MBA			

## MPA Y

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>C521015 + C547022</b>	<b>BKU61522</b>	
Brand: MBA			
<b>Head Gear Complete 1:1</b>	<b>C521015</b>	<b>BKU61015M</b>	
Brand: myonic			

Key




Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


\*\* Available while stocks last

*Not in current production*



## NAC Y

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C030011 + C032012</b>	<b>BKU61112</b>	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>C030011</b>	<b>BKU60011</b>	
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>C032012</b>	<b>BKU62012</b>	

## Nano 25LS

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C1038015 + C1038022</b>	<b>BKU61580</b>	

## Nano 95LS

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>C1034015</b>	<b>BKU64015</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>C1034017</b>	<b>BKU64017</b>	

Key





Usually stocked or Short Delivery

\*\* Available while stocks last





Long Delivery or Large Minimum Quantities

*Not in current production*

## NL45

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-KC	





## NL65 CH QD

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA	SR144-D28-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)  Brand: MBA	418NA-CER	SR144-D28-CSN-TO	





Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## NL9000L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TiX-TU03</b>	<b>BKUR250044</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TiX-TU03</b>	<b>BKUR2500445</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR2500444</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR2500443</b>	

## NL9000M

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TiX-MU03</b>	<b>BKUR250042</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-MU03</b>	<b>BKUR2500424</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TiX-MU03</b>	<b>BKUR2500425</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-MU03</b>	<b>BKUR2500423</b>	

Key


Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## NL9000S

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TiX-SU03</b>	<b>BKUR250043</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR2500433</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR2500434</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR250043T</b>	 Triple pack - Unit with 3 pieces
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TiX-SU03</b>	<b>BKUR2500435</b>	

## NL9000S Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b>  Brand: MBA	<b>0312-063100</b>	<b>BKU63100</b>	

## NL9000T

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR250044T</b>	 Triple pack - Unit with 3 pieces

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## NL95S

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**Ti-MU03**

**BKUR250045**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**Ti-SU03**

**BKUR250046**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## NL95T

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**Ti-TU03**

**BKUR250047**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## NPB Slow Speed

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**BB 4-7-2**

**SMR74-MC3-JP6-AF12**

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



**Replacement Bearing**

**BB 4-8-2**

**SMR84-MC3-JP6-AF12**

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)



## NSK Presto Rotor

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**418WIF**

**SFR144A-D27-ZSO-KC**

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

Brand: Myonic



**Bearing for Rotor Cartridge**

**418WIF-CER**

**SFR144A-D27-CSN-ZSO-KC**

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)







Brand: MBA



**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production






## NTF TU03

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production







## Optica

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic		<b>SR144A-D27-ZS-TO</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	


Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Pana Air Mini MU03 Red Band

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic		<b>SR144A-D27-ZS-TO</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	





## Pana Air Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer  Brand: MBA	<b>NPA-SU03</b>	<b>BKUR250056</b>	 <p>Precision sub- assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>



Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Pana Air SU03 - TU03 Green-Black

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic		<b>SR144A-D27-ZS-TO</b>	
<b>Bearing for Chuck End</b> High Speed Bearings - Polyamide Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)  Brand: MPD		<b>SFR144-D28-ZSO</b>	
<b>Bearing for Chuck End</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)  Brand: MPB		<b>SFR144-D28-CSN-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)  Brand: MPB		<b>SFR144-D28-KC</b>	




## Pana Max 2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>P1190050 - PAX2-SU03</b>	<b>BKU60050</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer  Brand: MBA	<b>PAX2-SU03</b>	<b>BKUR250068</b>	


Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production




## Pana Max S

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>PAX-SU</b>	<b>BKUR250051</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-SU</b>	<b>BKUR2500514</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-SU</b>	<b>BKUR2500513</b>	

## Pana Max S Plus

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>SX-SU03</b>	<b>BKUR250073</b>	





## Pana Max T

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>PAX-TU</b>	<b>BKUR250050</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-TU</b>	<b>BKUR2500503</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-TU</b>	<b>BKUR2500504</b>	

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>






## Phatelus

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic	<b>418WIF</b>	<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic	<b>418WIF</b>	<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Polyamide Retainer Precision flanged dental handpiece ceramic-hybrid bearing with polyamide retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIF-CER</b>	<b>SFR144A-D27-CSN-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIF-CER</b>	<b>SFR144A-D27-CSN-ZSO-KC</b>	







Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Phatelus II

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic		<b>SR144A-D27-ZS-TO</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	







## Phatelus III

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required Brand: NSK		<b>W0030-F-007-001-PT</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) Brand: Myonic		<b>SR144A-D27-ZS-TO</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Phatelus MU03

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	







## Phatelus Set Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b>      Brand: MBA	<b>Y900580</b>	<b>BKU60580</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production







## Phatelus SU03

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	


Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Phatelus TU03

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	

## Presto AQUA Cartridge

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Rotor / Cartridge</b>   Brand: MBA	<b>T816-016</b>	<b>BKUR250067</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

## Presto Aqua II

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Chuck</b>  True running accuracy less than 0.03  Brand: MBA	<b>T815019</b>	<b>BKU65019</b>	

**Key**

<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
Not in current production

## Presto Aqua Lux

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**T815019**

**BKU65019**

True running accuracy less than 0.03

Brand: MBA



## Presto Cartridge

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Rotor / Cartridge**

**T815-016**

**BKUR250066**

Brand: MBA



## PTL

Position / Description / Application

Manufacturers Part


MBA Part

**LED Bulb**

**Y900529**

**BKU60529LED**

Brand: MBA




**Xenon Lamp**

**Y900529**

**BKU60529Xenon**

Brand: MBA



## PTL 950TI

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B22227-7923**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



## QD J

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement O-Ring**

**0312-021083/063080**




**BKU61083**

Special dental applications O-Ring




Brand: MBA




## QD S

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>PAX-SU</b>	<b>BKUR250051</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-SU</b>	<b>BKUR2500514</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-SU</b>	<b>BKUR2500513</b>	


## QD T

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>PAX-TU</b>	<b>BKUR250050</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-TU</b>	<b>BKUR2500504</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>PAX-TU</b>	<b>BKUR2500503</b>	

## S Max M25 NEW

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C1038015 + C1142022</b>	<b>BKU61502</b>	

## S Max M25 OLD

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C601015 + 10000574</b>	<b>BKU61574</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## S Max M95L NEW

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**C600015 + C1141022**

**BKU61522**

Brand: MBA



## S Max pico

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**SX-PU03**

**BKUR250065**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## S Max SG20

Position / Description / Application

Manufacturers Part

MBA Part


**Drive Gear Assembly - Dog Clutch Complete 20:1**

**6 Parts**

**BKU63001**

Complete Assembly

Brand: MBA



**Drive Kit**

**C1010150 + C1010022**

**BKU65021**

Brand: MBA



**Replacement Middle Gear Complete 20:1**

**C1010022**

**BKU60022SG**

Brand: MBA



## S Max-Reg tm- M500

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**SX-MU03**

**BKUR250074**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## S Max-Reg tm- M600

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**SX-SU03**

**BKUR250073**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## S Max-Reg tm- M800

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**SX-MU03**

**BKUR250074**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## S Max-Reg tm- M900

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


**SX-SU03**

**BKUR250073**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



## S900 Airscaler

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Bush

**T695104B**

**BKU65104**

Refer OEM Part Number

Brand: MBA



## S950 Airscaler

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Bush

**T695104B**

**BKU65104**

Refer OEM Part Number

Brand: MBA



## SG20

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

**10000228**

**BKU60228**

Brand: MBA



## SG20L

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**C1010150 + C1003022**

**BKU65022**

Brand: MBA



### Head Gear Complete 20:1

**C1010150002**

**BKU60002**

Brand: MBA



Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## SG20X

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**10000228**

**BKU60228**

Brand: MBA



## SGAE

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B12213-0704**

**SMR74A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



## SGAE2G

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B12213-0704**

**SMR74A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



## SGM

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Middle Gear Complete 1:1**

**C293B031**

**BKU63031**

Brand: MBA



## SGMS

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Middle Gear Complete 1:1**

**C293B031**

**BKU63031**

Brand: MBA



## SGR2E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft**

**B10313-0703**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



**Bearing for Shaft**

**B12213-0704**

**SMR74A-ZZS-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x7x2.5)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## SGSE

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**B12213-0704**

**SMR74A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



## SGT2E0

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**B12213-0704**

**SMR74A-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



## SGT2E1

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

**B10313-0703**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## TEP Y

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete 1:1

**C837015**

**BKU67015TEP**

Brand: MBA




### Replacement Middle Gear Complete 1:1

**C837012**

**BKU67012**

Brand: MBA



## Ti Max® M95L

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

**C600015 + 10000664**

**BKU61564**

Brand: MBA



## Ti Max® Ti25L

Position / Description / Application

Manufacturers Part

MBA Part


### Drive Kit

**C875060 + C880022**

**BKU66022**

Dental replacement drive kit

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## Ti Max® X

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**C1010150 + C1003022**

**BKU65022**

Brand: MBA



## Ti Max® X25L

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**C601015 + C601022**

**BKU61522X**

Brand: MBA



## Ti Max® X500

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**

**TiX-MU03**

**BKUR250042**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA



**Replacement Turbine With Bearings**

**TiX-MU03**

**BKUR2500425**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



**Replacement Turbine With Bearings**

**TiX-MU03**

**BKUR2500424**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**

**TiX-MU03**

**BKUR2500423**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



\*\* Available while stocks last

*Not in current production*

## Ti Max® X600

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TiX-SU03</b>	<b>BKUR250043</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR250043T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TiX-SU03</b>	<b>BKUR2500435</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR2500433</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKUR2500434</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKURMU250043</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TiX-SU03</b>	<b>BKURM250043</b>	

## Ti Max® X650

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-SU03</b>	<b>BKURMU250043</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>TiX-SU03</b>	<b>BKURM250043</b>	


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production


## Ti Max® X700

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>TiX-TU03</b>	<b>BKUR250044</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>TiX-TU03</b>	<b>BKUR2500445</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR2500444</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR2500443</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>TiX-TU03</b>	<b>BKUR250044T</b>	

## Ti Max® X95L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C600015 + C600022</b>	<b>BKU61522</b>	

## Ti Max® Z25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C1038015 + C1038022</b>	<b>BKU61580</b>	

## Ti Max® Z45L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	<b>C1064015 + C1064022</b>	<b>BKU61522Z</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Ti Max® Z800L

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine Excluding Bearings

TiZ-SU03

BKUR250072

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

TiZ-SU03

BKUR2500723

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## Ti Max® Z900L

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine Excluding Bearings

TiZ-SU03

BKUR250069

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

TiZ-SU03

BKUR2500693

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## Ti Max® Z95

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

C1034D017 + C1034022

BKU61722

Brand: MBA



## Ti15

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Middle Gear Complete 1:1

C880022

BKU68022

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Ti25

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Middle Gear Complete 1:1**

**C880022**

**BKU68022**

Brand: MBA



**Replacement Press Button**

**C875050**

**BKU65050**

Brand: MBA



## Ti25L

Position / Description / Application

Manufacturers Part


MBA Part

**Head Gear Complete 1:1**

**C875060**

**BKU65060**

Brand: MBA



## Ti85

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**B19157-4723**

**SR133-MC3-TH7-TO**

High Speed Bearing - Torton Retainer

Open Abec 7 precision stainless steel bearing with high speed Torton retainer (2.38x4.76x1588)

Brand: myonic



## Ti95

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**B19157-4723**

**SR133-MC3-TH7-TO**

High Speed Bearing - Torton Retainer

Open Abec 7 precision stainless steel bearing with high speed Torton retainer (2.38x4.76x1588)

Brand: myonic



**Bearing for Dog Clutch**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)

**SMR84-MC3-JP6-AF12**



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



\*\* Available while stocks last

*Not in current production*



## UHR 50T C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11223-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20215-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	



## UHR45C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11243-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20245-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	

## UHR50C

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Dust Proof Cover</b> Laboratory Drive  Brand: MBA	<b>H063B106A</b>	<b>BKU63106</b>	
<b>Replacement Spindle Washer</b> Laboratory Drive  Brand: MBA	<b>H063B109A</b>	<b>BKU63109</b>	



## UHR50T

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Dust Proof Cover</b> Laboratory Drive  Brand: MBA	<b>H063B106A</b>	<b>BKU63106</b>	
<b>Replacement Spindle Washer</b> Laboratory Drive  Brand: MBA	<b>H063B109A</b>	<b>BKU63109</b>	



Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*



## UHRXL T

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11243-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20245-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	



## Ultimate 450

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11223-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20215-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	

## Ultimate 450-VH 2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11223-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20215-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	

## Ultimate 450-VH 3

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11223-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20215-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Ultimate XL

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft 6-13-5mm

**B11243-1306**

**S686B-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)

Brand: MBA



### Bearing for Shaft

**B20245-2763**

**SR188A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)

Brand: MBA



## Ultimate XL 450

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Chuck

**H203180**

**BKU63180**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Ultimate XL 500

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Chuck

**H203180**

**BKU63180**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## VH 2

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft 6-13-5mm

**B11223-1306**

**S686B-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)

Brand: MBA



### Bearing for Shaft

**B20215-2763**

**SR188A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)

Brand: MBA



Key









Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## Vip-Air

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required Brand: NSK	<b>Quad Ring</b>	<b>** QR-00600-080-N70-NSK</b>	
<b>Replacement O-Ring</b> High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00607-079-N70</b>	
<b>VOLVERE 8 30000 RPM</b>			
<b>Bearing for Chuck End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		<b>S695-ZZ-MC3-JP6</b>	
<b>Bearing for Cap End</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (6.35x12.7x4.76)		<b>SR188A-ZZ-MC3-JP6</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## VOLVERE GX 40M

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Cap End

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)

**SMR74-MC3-JP6-AF12**



### Bearing for Chuck End

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)

**SMR84-MC3-JP6-AF12**



## VOLVERE GX 6000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Chuck End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x6x2.5)

**SMR63A-ZZ-MC3-JP6**



### Bearing for Cap End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x7x3)

**S683A-ZZ-MC3-JP6**



## VOLVERE GX35M VH

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Chuck End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (6x13x5)

**S686B-ZZ-MC3-JP6**



### Bearing for Cap End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (6.35x12.7x4.76)

**SR188A-ZZ-MC3-JP6**



## Volvere i7

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Chuck

**H203180**

**BKU63180**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*



## VR E

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 6-13-5mm</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)  Brand: MBA	<b>B11243-1306</b>	<b>S686B-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)  Brand: MBA	<b>B20245-2763</b>	<b>SR188A-ZZS-MC3-GR3-T9H</b>	



## X

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>10000228</b>	<b>BKU60228</b>	 Press Button cpl. for Heads & Air Driven Handpieces

## X SG20L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Head Housing</b>  Brand: MBA	<b>10000254001</b>	<b>BKU64001</b>	 Head Housings
<b>Replacement Middle Gear Complete 20:1</b>  Brand: MBA	<b>C1003022</b>	<b>BKU63022</b>	 Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads

## X15L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>C601015</b>	<b>BKU61015</b>	 Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>C601022</b>	<b>BKU61022</b>	 Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads




Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## X25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427062</b>	<b>BKU67062</b>	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427064</b>	<b>BKU67064</b>	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427168</b>	<b>BKU67168</b>	
<b>Replacement Head Housing</b>  Brand: MBA	<b>C601011</b>	<b>BKU61011</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>C601050</b>	<b>BKU61050</b>	

## X25L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Head Housing</b>  Brand: MBA	<b>C601011</b>	<b>BKU61011</b>	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>C601015</b>	<b>BKU61015</b>	
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>C601022</b>	<b>BKU61022</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## X450

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Press Button**

**P1078500 - TiX45-SU04**

**BKU68500**

Brand: MBA



## X500

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**P067B050 - TiX-MU04**

**BKU67050**

Brand: MBA



## X600

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**SG1441X-D**

**SR144-D29-CSN-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision angular contact ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)

Brand: MBA



**Replacement Press Button**

**P068B050 - TiX-SU04**

**BKU68050**

Brand: MBA



## X600 Air Driven Handpiece

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement O-Ring**

**0312-063100**

**BKU63100**

Brand: MBA



## X65

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Chuck**

**H190180**

**BKU60180**

True running accuracy less than 0.03

Brand: MBA



## X650

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**SG1441X-D**

**SR144-D29-CSN-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision angular contact ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## X700

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**P069B050 - TiX-TU04**

**BKU69050**

Brand: MBA



## X95

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**740-CER**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



**Bearing for Dog Clutch**

**SMR84-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)



**Replacement O-Ring**

**C427062**

**BKU67062**

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring

Brand: MBA



**Replacement O-Ring**


**C427064**

**BKU67064**

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring

Brand: MBA



**Replacement O-Ring**

**C427168**

**BKU67168**

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring

Brand: MBA



**Head Gear Complete**

**C600015**

**BKU60015**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



Key




Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


\*\* Available while stocks last

*Not in current production*





## X95L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Dog Gear 1:5 Dog Clutch C600C013</b> Clutch  Brand: MBA		<b>BKU69993</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Head Housing</b>  Brand: MBA	<b>C600011</b>	<b>BKU60011A</b>	 <small>Head Housings</small>
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>C600022</b>	<b>BKU60022</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>

## XSG20L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 20:1</b>  Brand: MBA	<b>C1010150002</b>	<b>BKU60002</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>

## Z25

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427062</b>	<b>BKU67062</b>	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring  Brand: MBA	<b>C427168</b>	<b>BKU67168</b>	
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>C1038022</b>	<b>BKU68022N</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Press Button</b>  Brand: MBA	<b>C1038050001</b>	<b>BKU60001</b>	 <small>Press Button cpl. for Heads &amp; Air Driven Handpieces</small>

Key **Usually stocked or Short Delivery**

**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities

*Not in current production*

## Z25L

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 1:1**

**C1038015**

**BKU68015**

Brand: MBA



## Z45L

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 4:1.2**

**C1064015**

**BKU64015Z**

Brand: MBA



**Replacement Middle Gear Complete 4:1.2**

**C1064022**

**BKU64022Z**

Brand: MBA



## Z800

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**P1051050 - TiZ-MU04**

**BKU61050Z**

Brand: MBA



## Z900

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Press Button**

**P1052050 - TiZ-SU04**

**BKU62050**





Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Z95L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Gear Assembly - Dog Clutch Complete 1:5</b> Complete Assembly  Brand: MBA	<b>C1034013</b>	<b>BKU64013</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>C1034015</b>	<b>BKU64015</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>C1034017</b>	<b>BKU64017</b>	
<b>Replacement Middle Gear Complete 1:5</b>  Brand: MBA	<b>C1034022</b>	<b>BKU64022N</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Accutorq

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

SR144A-D27-ZS-TO



## Borden

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA

SR144-D28-KC



### Replacement Bearing


High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

Brand: MBA

SR144-D28-CSN-KC



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## ACCESS PRO NL

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2009914</b>	<b>BKUR310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKURC310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKUR3100004</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKUR3100003</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2009914</b>	<b>BKUR3100005</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## KL Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2009914</b>	<b>BKUR310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKUR3100004</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKUR3100003</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009914</b>	<b>BKURC310000</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2009914</b>	<b>BKUR3100005</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## ROTAMAX PRO KL Torque

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2009915</b>	<b>BKUR310001</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2009915</b>	<b>BKUR3100015</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKUR3100013</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKUR3100014</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKURC310001</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## ROTAMAX PRO NL

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>2009915</b>	<b>BKUR310001</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKUR3100013</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKUR3100014</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>2009915</b>	<b>BKURC310001</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>2009915</b>	<b>BKUR3100015</b>	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## 102L Strong

Position / Description / Application


Manufacturers Part

MBA Part

**Replacement Armature**

**BKU61102**

Brand: MBA



## 103

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**BKU61103**

Brand: MBA



## 106

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**BKU61106**

Brand: MBA



## 204 Strong 90

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Armature**

**BKU61102**

Brand: MBA



## Forte 100

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**BKU68882**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 100 II

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU61073**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 100a III

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**BKU68882**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Forte 102

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU68883**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 103L

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU61073**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 106

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU68882**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 106ZO

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU68882**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 300 II

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU68883**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA



## Forte 300 III

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**107-3**

**BKU68883**

True running accuracy less than 0.03, bore diameter 2.35mm

Brand: MBA





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Lab Drive

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	<b>830ZZ</b>	<b>S693A-ZZS-MC3-T9H</b>	
<b>Bearing for Motor Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (4x9x4)  Brand: MBA	<b>940ZZ</b>	<b>S684A-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	<b>1040ZZ</b>	<b>SMR104A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (5x10x4)  Brand: MBA	<b>1050ZZ</b>	<b>SMR105A-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide Cage (7x13x4)  Brand: MBA	<b>1370ZZ</b>	<b>SMR137A-ZZS-MC3-GR3-T9H</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)  Brand: ECO	<b>1480ZZ</b>	<b>SMR148A-ZZS-MC3-GR3-T9H</b>	

## Micro-Motor

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	<b>SAESHIN003</b>	<b>S693A-ZZS-MC3-T9H</b>	
<b>Replacement Bearing</b> 4x9x4 Stainless Steel Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x9x4)	<b>SAESHIN004</b>	<b>S684A-ZZ-MC3-JP6</b>	

Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Strong 102

Position / Description / Application


Manufacturers Part

MBA Part

**Replacement Motor Brush**

**BKU60107**

Brand: MBA



## Strong 105

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Motor Brush**

**BKU60107**

Brand: MBA



## Strong 105L

Position / Description / Application


Manufacturers Part

MBA Part

**Replacement Motor Brush**

**BKU60107**

Brand: MBA



## ZO

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Chuck**

**BKU68882**

True running accuracy less than 0.03, bore diameter 2.35mm





Brand: MBA







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## C1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)  Brand: MBA	4120	S688B-ZZS-MC3-GR3	
<b>Replacement Felt Washer</b> Laboratory Drive  Brand: MBA	4129	BKU64129	




## C2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)  Brand: MBA	4120	S688B-ZZS-MC3-GR3	
<b>Replacement Felt Washer</b> Laboratory Drive  Brand: MBA	4129	BKU64129	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*



## C3

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Precision Bearing Special design for Schick dental handpiece replacement bearings part number 7120 (5x9.52x3.75)	7212	SWMR952505A-ZZS-MC3-GR7	
<b>Bearing for Shaft 8x15.875x6</b> High Speed Bearing - Phenolic Retainer Double shielded precision stainless steel bearing with high speed polyamide retainer (8x15.88x6)	7111	SMR15875-08-ZZ-MC3-KC	
<b>Bearing for Shaft</b> High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with high speed phenolic cage (10.15x19x7)	7120	SMR19-10.15-ZZS-MC3-GR3-KC	

## C3 7050CN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded Abec 3 precision stainless steel bearing (8x16x6)  Brand: MBA	4120	S688B-ZZS-MC3-GR3	




## CN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
<b>Bearing for Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## Q Profi

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Motor Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	<b>5212</b>	<b>S684A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	<b>611/1</b>	<b>S687A-ZZ-MC3-JP6</b>	
<b>Bearing for Shaft</b> High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)  Brand: MBA	<b>4120</b>	<b>S688B-ZZS-MC3-GR3</b>	




## Qube

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Shaft 5x11x5</b> High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (5x11x5)  Brand: MBA	<b>9160</b>	<b>S685A-ZZS-MC3-GR3-KC</b>	
<b>Bearing for Shaft</b> High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (6x13x5)  Brand: MBA	<b>9162</b>	<b>S686B-ZZS-MC3-GR3-KC</b>	


Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*





## 001

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6710425</b>	<b>BKUR2200313</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6710433</b>	<b>BKUR2200283</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6710441</b>	<b>BKUR2200343</b>	

## 3000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required  Brand: myonic		<b>SR144-D28-ZS-KC</b>	

## 4000 MS-ML

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-KC</b>	
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-TO</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-KC</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-TO</b>	





Key Usually stocked or Short Delivery \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 4000 S-SL Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-TO</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>





## P Control

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly excluding bearings Brand: MBA	<b>6084979</b>	<b>BKUR220022</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>6084979</b>	<b>BKUR2200224</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>6084979</b>	<b>BKUR2200225</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>6084979</b>	<b>BKUR2200223</b>	





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## P Racer

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6084987</b>	<b>BKUR2200233</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6084987</b>	<b>BKUR2200234</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>6084987</b>	<b>BKUR2200235</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly excluding bearings  Brand: MBA	<b>6084987</b>	<b>BKUR220023</b>	

## Siroboost

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>6334598</b>	<b>BKUR220024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>6334598</b>	<b>BKUR2200245</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6334598</b>	<b>BKUR2200243</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>6334598</b>	<b>BKUR2200244</b>	

Key





Usually stocked or Short Delivery

\*\* Available while stocks last





Long Delivery or Large Minimum Quantities

Not in current production

## T1 Control

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>4175217</b>	<b>BKUR220020</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>4175217</b>	<b>BKUR2200205</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4175217</b>	<b>BKUR2200204</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4175217</b>	<b>BKUR2200203</b>	



## T1 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b>  Turbine sub assembly excluding bearings  Brand: MBA	<b>4178591</b>	<b>BKUR220018</b>	
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4178591</b>	<b>BKUR2200184</b>	
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4178591</b>	<b>BKUR2200183</b>	
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>4178591</b>	<b>BKUR2200185</b>	





<b>Usually stocked or Short Delivery</b>
<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T1 Mini TMC3LK

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-TO</b>	
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>

## T2 Control

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>5868299</b>	<b>BKUR220022</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>5868299</b>	<b>BKUR220025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5868299</b>	<b>BKUR220024</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5868299</b>	<b>BKUR220023</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## T2 Control Serial Numbers 602 001 to 802 000 with

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings

**6432590**

**BKUR2200303**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: myonic



## T2 Control Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings

**6710425**

**BKUR2200313**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Mini

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine Excluding Bearings

**5957142**

**BKUR220021**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

**5957142**

**BKUR2200215**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




### Replacement Turbine With Bearings

**5957142**

**BKUR2200214**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




### Replacement Turbine With Bearings

**5957142**

**BKUR2200213**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Mini Serial Numbers 602 001 to 702 000 with Wed

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine With Bearings

**6432632**

**BKUR2200323**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T2 Mini Serial Numbers 702 001 to 802 000 with Wed

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6548130**

**BKUR2200333**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Mini Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**

**6710441**

**BKUR2200343**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Mini Serial Numbers 802 001

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6710441**

**BKUR2200343**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Racer

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**

**5868307**

**BKUR220023**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine With Bearings**

**5868307**

**BKUR2200233**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**

**5868307**

**BKUR2200234**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**5868307**

**BKUR2200235**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T2 Serial Numbers 602 001 to 702 000

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine With Bearings**


**6432616**

**BKUR2200263**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**

**6548122**

**BKUR2200273**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6710433**

**BKUR2200283**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Serial Numbers 802 001

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6710433**

**BKUR2200283**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Boost Serial Numbers 602 001 to 702 000 with We

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6432616**

**BKUR2200263**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**6548122**

**BKUR2200273**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Boost Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6710433**

**BKUR2200283**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T3 Boost Serial Numbers 802 001

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6710433**

**BKUR2200283**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Mini

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**

**5957142**

**BKUR220021**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




**Replacement Turbine With Bearings**

**5957142**

**BKUR2200215**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




**Replacement Turbine With Bearings**

**5957142**

**BKUR2200213**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**5957142**

**BKUR2200214**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Mini Serial Numbers 602 001 to 702 000 with Wed

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6432632**

**BKUR2200323**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Mini Serial Numbers 702 001 to 802 000 with Wed

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**6548130**

**BKUR2200333**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T3 Mini Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine With Bearings

**6710441**

**BKUR2200343**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Mini Serial Numbers 802 001

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine With Bearings

**6710441**

**BKUR2200343**

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T3 Racer

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine Excluding Bearings

**5868307**

**BKUR220023**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

**5868307**

**BKUR2200234**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




### Replacement Turbine With Bearings

**5868307**

**BKUR2200235**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA




### Replacement Turbine With Bearings

**5868307**

**BKUR2200233**

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic





Brand: myonic







<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## T3 Racer Serial Numbers 600 000

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>6486240</b>	<b>BKUR220025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>6486240</b>	<b>BKUR2200255</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer  Brand: myonic	<b>6486240</b>	<b>BKUR2200254</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: myonic	<b>6486240</b>	<b>BKUR2200253</b>	





## TC3

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>4175217</b>	<b>BKUR220020</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4175217</b>	<b>BKUR2200204</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>4175217</b>	<b>BKUR2200203</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>4175217</b>	<b>BKUR2200205</b>	





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TM1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded, stepped OD, high precision stainless steel bearing with ceramic balls and Phenolic retainer (3.18x6.35x2.38) 2 Required Brand: MBA		<b>SR144-D106-CSN-ZS-KC</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing - Alternative</b> High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required Brand: myonic		<b>SR144-D106-ZS-TO</b>	

## TMC

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings Brand: MBA	<b>4178591</b>	<b>BKUR220018</b>	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer Brand: MBA	<b>4178591</b>	<b>BKUR2200185</b>	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>4178591</b>	<b>BKUR2200184</b>	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic Brand: myonic	<b>4178591</b>	<b>BKUR2200183</b>	 <small>Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</small>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TS1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic		<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic		<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>5653054</b>	<b>BKUR2200175</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5653054</b>	<b>BKUR2200174</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>5953054</b>	<b>BKUR220017</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b>  Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5953054</b>	<b>BKUR2200173</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TS2

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic		<b>SR144A-D105-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic		<b>SR144A-D105-ZS-TO</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: MPB	<b>0.220.2007, 220-0793K</b>	<b>SR144A-D107-CSN-ZS-TO</b>	 <b>Image Coming Soon</b>
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-KC</b>	
<b>Replacement Bearing</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required Brand: Myonic	<b>220-0125</b>	<b>SR144A-D105-CSN-ZS-TO</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>5653054</b>	<b>BKUR2200175</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5653054</b>	<b>BKUR2200174</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>5953054</b>	<b>BKUR220017</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>5953054</b>	<b>BKUR2200173</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 000

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Press Button**

**5828889**

**BKU68889**

Brand: MBA



## 001

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**


**6554815**

**SWR144B-CSN-ZE-MC3-GR7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)

Brand: MBA



## 001X

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Spring Washer**

**6367648**

**BKU67648**

Brand: MBA



## 13Z

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete**

**6175587 / 6405687**

**BKU65587**

For Handpiece with Wedge Type Chuck

Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece

Brand: MBA



**Head Gear Complete - Standard**

**6175587 / 6405687**

**BKU65587X**

For Handpiece with Wedge Type Chuck

Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: MBA



## 13Z Mini

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete**

**6414457 / 6405679**

**BKU65687**

For Handpiece with Wedge Type Chuck

Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece

Brand: MBA



## 13Z Mini Standard

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete - Standard**

**6414457 / 6405679**

**BKU65687X**

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## 40

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)

Brand: MBA



## 6

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)

Brand: MBA



## 6:1 VDW

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Headgear

5887174

SMR74D-MC3-GR5-L23

High Precision Bearing

Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)

Brand: Myonic



## A170

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete

6414457 / 6405679

BKU65687M

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: myonic



## A200

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Middle Gear Complete 1:5

6175603

BKU65603

Brand: MBA




### Replacement Middle Gear Complete 1:5

6175603

BKU65603WB

Brand: MBA



## A200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

6049337

BKU69337

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## A40L

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Head Gear 1:1**

**6000983**

**BKU60983**

Brand: MBA




**Replacement Middle Gear Complete 1:1**

**6000959**

**BKU60959**

Brand: MBA



## C 1.6

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Drive**

**4181157**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## C 1.6 Serial Number Over 65 000

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Drive**

**6360577**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## C 200 40 6 1.6

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Drive**

**6360601**

**SMR1035-06-MC3-GR7-TO**

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)

Brand: MBA



## C+S 200

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Drive**

**6360601**

**SMR1035-06-MC3-GR7-TO**

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)

Brand: MBA



## C1.6

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 6:1**

**6405703**

**BKU65703**

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## C200

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete</b> Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece Brand: MBA	6175587 / 6405687	BKU65587X	
<b>Head Gear Complete</b> Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece Brand: MBA	6175587 / 6405687	BKU65587	
<b>Head Gear Complete</b> Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece Brand: MBA	6414457 / 6405679	BKU65687X	
<b>Head Gear Complete</b> Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece Brand: MBA	6414457 / 6405679	BKU65687	
<b>Replacement O-Ring</b> Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring Brand: MBA	6328236	BKU68236	
<b>Replacement Middle Gear Complete 1:5</b> Brand: MBA	6175603	BKU65603WB	
<b>Replacement Middle Gear Complete 1:5</b> Brand: MBA	6175603	BKU65603	
<b>C200 T1</b>			
<b>Head Beaing</b> High Speed Bearing - Torlon Retainer Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2) Brand: MBA		SMR7-032-MC3-GR7-TO-AF12	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## C200L 8.7 Diameter

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete - Premium

6414457 / 6405679

BKU65687M

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

Brand: myonic



## C200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

6049337

BKU69337

Brand: MBA



## C40

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Rear Bearing

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)

S683-MC3-JP6-AF12



### Replacement O-Ring

Special Dental Application O-Ring

Special dental applications O-Ring

Brand: MBA

6328236

BKU68236



### Replacement Press Button

4 Parts

BKU65083

Brand: MBA



## C40 Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

6049329

BKU69329

Brand: MBA



Key




Usually stocked or Short Delivery

\*\* Available while stocks last



Long Delivery or Large Minimum Quantities

*Not in current production*


## C40L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Head Gear 1:1</b>	<b>6000983</b>	<b>BKU60983</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
Brand: MBA			
<b>Replacement Head Gear 1:1 Mini</b>	<b>6405695</b>	<b>BKU65695</b>	
Brand: MBA			 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Middle Gear Complete 1:1</b>	<b>6000959</b>	<b>BKU60959</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
Brand: MBA			


## C40L mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b>	<b>6360593</b>	<b>SMR63-MC3-JP6-AF12</b>	
High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)			
Brand: Myonic			
<b>Bearing for Headgear</b>	<b>5887174</b>	<b>SMR74D-MC3-GR5-L23</b>	
High Precision Bearing Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)			

## C6

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 6:1</b>	<b>6405703</b>	<b>BKU65703</b>	 <p>Drive Kits for Electric Driven Handpieces and Heads</p>
Brand: MBA			



## C6L mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b>	<b>6360593</b>	<b>SMR63-MC3-JP6-AF12</b>	
High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)			




Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*


## EL1

Position / Description / Application	Manufacturers Part	MBA Part	
<b>LED Bulb</b>	<b>5945407</b>	<b>BKU60291LED</b>	
Brand: MBA			
<b>Xenon Lamp</b>	<b>5945407</b>	<b>BKU65407Xenon</b>	
Brand: MBA			


## Endo

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 6:1</b>	<b>6000777</b>	<b>BKU60777</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 6:1</b>	<b>5994038</b>	<b>BKU64038</b>	
Brand: MBA			
<b>Replacement Press Button</b>	<b>5997528</b>	<b>BKU67528</b>	
Brand: MBA			

## ENDO 6:1 VDW

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>6000777 + 5994038</b>	<b>BKU67738</b>	
Brand: MBA			

## Endo X

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 6:1</b>	<b>H1033C1051015</b>	<b>BKU61015X</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 6:1</b>	<b>H1033C1051022</b>	<b>BKU61022X</b>	
Brand: MBA			

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## H40 6

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

4181132

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)

Brand: MBA



## H40L

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

4181199

S623-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (3x10x4)



### Bearing for Shaft

4181108

SMR85-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



## Quick Coupling

Position / Description / Application

Manufacturers Part


MBA Part

### LED Bulb and Cap

6314558

BKU64558LED

Brand: MBA



## Quick Coupling Sirosonic L

Position / Description / Application

Manufacturers Part


MBA Part

### LED Bulb

5940291

BKU60291LED

Brand: MBA




### Xenon Lamp

5940291

BKU60291Xenon

Brand: MBA



## S 12000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Chuck End

SF692-ZZS-MC34-GR5-R-L23

Stainless Flanged High Precision Bearing

Double shielded flanged Abec 5 precision stainless steel bearing with NON-ISO Width and low-torque ribbon retainer (2x6x2.3)

Brand: Myonic



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## S 200 40 6 1.6

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)

Brand: MBA



## Siroboost

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


6334598

BKUR220024T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## SL

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor

7036130

S692-ZZS-MC34-GR5-R-G48

High Precision Bearing Shielded Option

Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)

Brand: Myonic



### Bearing for Rotor

7040801

S684A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x9x4)




### LED Bulb

5945407

BKU60291LED

Brand: MBA




### Xenon Lamp

5945407

BKU65407Xenon

Brand: MBA



## Smart Plus

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete 6:1

H1033C1051015

BKU61015X

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## Smart Plus 6:1

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Middle Gear Complete 6:1**

**H1033C1051022**

**BKU61022X**

Brand: MBA



## T1 Boost Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**6554815**

**SWR144B-CSN-ZE-MC3-GR7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)

Brand: MBA



## T1 Classic S200

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**5890855**

**SMR7-032-MC3-GR7-TO-AF12**

High Speed Bearing - Torlon Retainer

Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2)

Brand: MBA



**Bearing for Headgear**

**5890855**

**SMR7-032-CSN-MC3-GR7-TO-AF12**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (3.18x7x2)

Brand: MBA



## T1 CLASSIC S200L

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**6175587 + 6175603**

**BKU68703**

Dental replacement drive kit

Brand: MBA



## T1 CLASSIC S200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

**Drive Kit**

**6405679 + 6175603**

**BKU67903**

Brand: MBA



**Drive Kit**

**6414457 + 6175603**

**BKU65703S**

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## T1 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement O-Ring

7028988

BKU68988

Special Dental Application O-Ring

Brand: MBA



## T1 Control Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Spring Washer

6367648

BKU67648

Brand: MBA



## T1 Control Serial Number Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge


6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)

Brand: MBA



## T1 Control till Serial Number 65

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

5828889

BKU68889

Brand: MBA



## T1 Control till Serial Number 65 000

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

5828889

BKU68889

Brand: MBA



## T1 Line C1.6L

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

6405703 + 6000942

BKU60342




Brand: MBA




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## T1 Line C200

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (3.18x7x2)  Brand: MBA	5890855	SMR7-032-CSN-MC3-GR7-TO-AF12	
<b>Bearing for Headgear</b> High Speed Bearing - Torlon Retainer Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2)  Brand: MBA	5890855	SMR7-032-MC3-GR7-TO-AF12	
<b>Replacement Press Button</b>  Brand: MBA	4 Parts	BKU66967	


## T1 Line C200L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Dental replacement drive kit  Brand: MBA	6175587 + 6175603	BKU68703	


## T1 Line C200L Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Brand: MBA	6405679 + 6175603	BKU67903	

## T1 Line C40L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Dental replacement drive kit  Brand: MBA	6000983 + 6000959	BKU68359	

## T1 Line C40L Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>  Dental replacement drive kit  Brand: MBA	6405695 + 6000959	BKU69559	

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## T1 Line C6L

Position / Description / Application

Manufacturers Part


MBA Part

**Drive Kit**

**6405703 + 6000942**

**BKU60342**

Brand: MBA



## T2

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**5868307**

**BKUR220023T**

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Boost Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**


**6554815**

**SWR144B-CSN-ZE-MC3-GR7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)

Brand: MBA



## T2 Boost Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Spring Washer**

**6367648**

**BKU67648**

Brand: MBA



## T2 Control

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**5994814**

**SR144A-D27-CSN-ZE-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer and encapsulated shield one side (3.18x6.35x2.78)

Brand: MBA



**Replacement Turbine With Bearings**


**5868299**

**BKUR220022T**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## T2 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement O-Ring**

**5829416**

**BKU69416**

Special Dental Application O-Ring

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## T2 Control Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge


6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)

Brand: MBA



## T2 Control Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Spring Washer

6367648

BKU67648

Brand: MBA



## T2 Line A200L

Position / Description / Application

Manufacturers Part


MBA Part

### Drive Kit

6562586 + 6175603

BKU68603

Brand: MBA



## T2 Line A40 Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

6049352 + 6000959

BKU65259

Brand: MBA



## T2 Line A40L

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

6000983 + 6000959

BKU68359

Dental replacement drive kit

Brand: MBA



## T2 Racer

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge

5994814

SR144A-D27-CSN-ZE-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer and encapsulated shield one side (3.18x6.35x2.78)

Brand: MBA




### Replacement Press Button

5828897

BKU68897

Brand: MBA



Key **Usually stocked or Short Delivery**  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## T2 Racer Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement O-Ring

5829416

BKU69416

Special Dental Application O-Ring

Brand: MBA



## T2 Revo

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

4 Parts

BKU65083

Brand: MBA



## T2 Revo Mini

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

6049329

BKU69329

Brand: MBA



## T2 Revo R170

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Middle Gear

7044654

SMR55-0235-T4HB-MC3-AF12

High Speed Bearing - Torlon Retainer

Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)

Brand: MBA



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


\*\* Available while stocks last

*Not in current production*


## T2 Revo R40

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>5943092 + 5943167</b>	<b>BKU69267</b>	
Brand: MBA			
<b>Replacement Head Gear 1:1</b>	<b>5943092</b>	<b>BKU63092</b>	
Brand: MBA			
<b>Replacement Head Gear 1:1 Mini</b>	<b>6049352</b>	<b>BKU69352</b>	
Brand: MBA			
<b>Replacement Middle Gear Complete 1:1</b>	<b>5943167</b>	<b>BKU63167</b>	
Brand: MBA			

## T2 Revo R40 mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>6405695 + 5943167</b>	<b>BKU69567</b>	
Brand: MBA			

## T2 S

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>6562586 + 6175603</b>	<b>BKU68603</b>	
Brand: MBA			

## T3 Line E40 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Drive Kit</b>	<b>6049352 + 6000959</b>	<b>BKU65259</b>	
Brand: MBA			

Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## T3 Racer

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Turbine With Bearings


5868307

BKUR220023T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## TEH40

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft

7025273

SMR105-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



## TS1 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement O-Ring

7028988

BKU68988

Special Dental Application O-Ring

Brand: MBA



## VDW

Position / Description / Application

Manufacturers Part

MBA Part

### Head Gear Complete 6:1

6000777

BKU60777

Brand: MBA



### Replacement Middle Gear Complete 6:1

5994038

BKU64038

Brand: MBA



## VDW 6:1

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Press Button

5997528

BKU67528

Brand: MBA



## X SMART PLUS 6:1

Position / Description / Application

Manufacturers Part

MBA Part

### Drive Kit

H1033C1051015 + H1033C1051022

BKU615022

Brand: MBA



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

Not in current production

## SS White

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Cap End</b> High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-TO</b>	
<b>Bearing for Cap End</b> High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)  Brand: Myonic		<b>SFR144A-D27-ZSO-KC</b>	
<b>Bearing for Chuck End</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-CSN-KC</b>	
<b>Bearing for Chuck End</b> High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)  Brand: MBA		<b>SR144-D28-KC</b>	





Key **Usually stocked or Short Delivery**  
 \*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*


## 430

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic LUBEFREE  Brand: myonic	<b>263739</b>	<b>BKURMU290108</b>	
<b>Turbine Sub Assembly Without Bearings</b> For Handpiece with Wedge Type Chuck  Brand: MBA	<b>263739</b>	<b>BKURM290108</b>	

## 430 Vista

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>263968</b>	<b>BKUR290109</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>263968</b>	<b>BKUR2910945</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>263968</b>	<b>BKUR2901094</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Phenolic retainer - myonic LUBEFREE  Brand: myonic	<b>263968</b>	<b>BKUR2901083</b>	





## 430 Vista Lubefree

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)  Brand: MBA	<b>SG1303X-D</b>	<b>SFR144A-D1303-CSN-ZSO-KC</b>	





<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## 430 Vista Manual Chuck

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> Daily Lube - High Speed Bearing - Phenolic Retainer Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIFU</b>	<b>SFR144A-D37-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>418WIFU-CER</b>	<b>SFR144A-D37-CSN-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> Daily Lube - High Speed Bearing - Phenolic Retainer Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>SG1186X</b>	<b>SFR144A-D1186-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>SG1186X-D</b>	<b>SFR144A-D1186-CSN-ZSO-KC</b>	

## 430SWL

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		<b>OR-00600-080-N70</b>	
<b>Replacement Washer - Thick</b> Precision Washer 6.5 x 8.30 x 0.5 Stainless (6x8x0.5)		<b>W0060-FP-008-0050-S4</b>	
<b>Replacement Washer - Thin</b> Precision Washer 6.5 x 8.36 x 0.2 Stainless (6x8x0.2)		<b>W0060-FP-008-0020-S4</b>	
<b>Replacement Wave Washer</b> 6.4mm Shaft Size (6.35x9.32x0.76)		<b>W0064-T-010-008-W3-S2</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## C6

Position / Description / Application

Manufacturers Part

MBA Part

**Xenon Lamp**

**262522**

**BKU62522Xenon**

Brand: MBA



## COLORight System

Position / Description / Application

Manufacturers Part

MBA Part

**Xenon Lamp**

**260986**

**BKU60986Xenon**

Brand: MBA



## Concentrix

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**SR2-5A-D17-ZS-TO**

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



## Futura I and II

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**SR2-5A-D17-ZS-TO**

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


\*\* Available while stocks last

*Not in current production*


## Solara

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>263968</b>	<b>BKUR290109</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with singular contact ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>263968</b>	<b>BKUR2901083</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>263968</b>	<b>BKUR2910945</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>263968</b>	<b>BKUR2901094</b>	


## Solara Lubefree

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)  Brand: MBA	<b>SG1303X-D</b>	<b>SFR144A-D1303-CSN-ZSO-KC</b>	

## Star 430 Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement O-Ring</b> Special Dental Application O-Ring  Brand: MBA	<b>O-Ring+Washer+Disk</b>	<b>BKU69997</b>	

## StarBright ISO

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Xenon Lamp</b>  Brand: MBA	<b>262522</b>	<b>BKU62522Xenon</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Starflex Lubefree

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge

**SG1303X-D**

**SFR144A-D1303-CSN-ZSO-KC**

Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)

Brand: MBA



## Ultra

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic



## Vista PB

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor Cartridge

**418WIFU**

**SFR144A-D37-ZSO-KC**

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)

Brand: MBA



### Bearing for Rotor Cartridge

**418WIFU-CER**

**SFR144A-D37-CSN-ZSO-KC**

Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)

Brand: MBA



### Bearing for Rotor Cartridge

**SG1186X**

**SFR144A-D1186-ZSO-KC**

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)

Brand: MBA



### Bearing for Rotor Cartridge

**SG1186X-D**

**SFR144A-D1186-CSN-ZSO-KC**

Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)

Brand: MBA



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## Starflite

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## EVO M5 L

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete**

**97230612**

**BKU60612M**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



**Head Gear Complete**

**97230612**

**BKU60612**

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

Brand: MBA



Key

**Usually stocked or Short Delivery**

Long Delivery or Large Minimum Quantities

**\*\* Available while stocks last**

*Not in current production*

## Miera

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## Thera

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## Titan

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## Turbina Master

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR2-5A-D17-ZS-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)

2 Required



## Turbina Mini Master

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## X 6000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x9x4)

2 Required

**SMR93A-ZZ-MC3-JP6**



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## 115K 14000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

00542400

S625-ZZ-MC3-JP6

Stainless Bearing  
Double shielded P6 precision stainless steel bearing (5x16x5)  
2 Required



## 260 Blue and Black

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Shaft 7-14-5mm

03202200

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer  
Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)



Brand: MBA

### Bearing for Shaft 8-16-6mm

03202300

S688B-ZZ-MC3-GR3-T9H-G36

High Speed Bearing - Polyamide Retainer  
Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



## 300 Series

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine Excluding Bearings

07234102

BKUR230050

For Handpiece with Wedge Type Chuck  
Turbine sub assembly excluding bearings

Brand: MBA




### Replacement Turbine With Bearings

07234102

BKUR2300504

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

Brand: myonic



### Turbine Sub Assembly With Bearings

07234102

BKUR2300502

For Handpiece with Wedge Type Chuck  
With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic

Brand: myonic



## 400 Series

Position / Description / Application

Manufacturers Part

MBA Part


### Replacement Turbine With Bearings

07005900

BKUR230039T

For Handpiece with Wedge Type Chuck  
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



Key



Usually stocked or Short Delivery

\*\* Available while stocks last





Long Delivery or Large Minimum Quantities

*Not in current production*

## 56A

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	<b>BB 4-7-2</b>	<b>SMR74-MC3-JP6-AF12</b>	
<b>Replacement Bearing</b> Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	<b>BB 5-8-2</b>	<b>SMR85-MC3-JP6-AF12</b>	

## 796

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>02765500</b>	<b>BKUR230025</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02765500</b>	<b>BKUR2300253</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02765500</b>	<b>BKUR2300254</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>02765500</b>	<b>BKUR2300255</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 798

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>02760000</b>	<b>BKUR230026</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02760000</b>	<b>BKUR2300263</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>02760000</b>	<b>BKUR2300265</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02760000</b>	<b>BKUR230026T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02760000</b>	<b>BKUR2300264</b>	

## 896

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>02765500</b>	<b>BKUR230025</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02765500</b>	<b>BKUR2300254</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>02765500</b>	<b>BKUR2300253</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>02765500</b>	<b>BKUR2300255</b>	

Key

**Usually stocked or Short Delivery**

**\*\* Available while stocks last**



Long Delivery or Large Minimum Quantities

*Not in current production*


## 898

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	02760000	BKUR230026	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	02760000	BKUR2300265	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	02760000	BKUR230026T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	02760000	BKUR2300263	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	02760000	BKUR2300264	

## 969

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)  Brand: MBA	01579600	SMR55-0235-T4HB-MC3-AF12	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	

## 99

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	01565800	SMR85-MC3-JP6-AF12	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## 999

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Lower Shaft

01589400

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## 999E

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Headgear

04486000

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



## A20

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Rotor

01850600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



## A25

Position / Description / Application

Manufacturers Part

MBA Part


### Bearing for Rotor

01850600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



## ADEC

Position / Description / Application

Manufacturers Part

MBA Part

### Bearing for Lower Shaft

01565800

SMR85-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



## Airmotor 13000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

01850600

SR2-5A-ZZ-MC3-JP6

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)

2 Required



## AM 25

Position / Description / Application

Manufacturers Part


MBA Part

### Xenon Lamp

04300400

BKU60400Xenon

Brand: MBA



Key **Usually stocked or Short Delivery**  
**\*\* Available while stocks last**

Long Delivery or Large Minimum Quantities  
*Not in current production*

## Contra Angle Handpiece New

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**02649400**

**SMR74-ZS-T4HB-MC3-GR5-AF12**

High Precision Bearing

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)

2 Required



## Contra Angle Handpiece Old

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

**01579300**

**SMR55-0235-MC3-TO-AF12**

High Speed Bearing - Torlon Retainer

Open precision dental handpiece bearing with Torlon retainer (2.35x5.5x2)

2 Required



## EB40

Position / Description / Application

Manufacturers Part

MBA Part


**Dental Motor Set with 21 parts**

**1502217-001**

**BKU62217**

Includes parts - 3x 01862300 - 1x 04062500 - 1x 04061800 - 1x 04062300 - 1x 04058800 - 2x 04062100 - 2x 04062200 - 2x 04062400 - 2x 04058800 - 3x 04059200 - 2x 04059300 - 1x 04069600

Brand: MBA




Armature, Rotors, Electric Dental Motors & Laboratory Drives

**Replacement Motor Brush**

**04059300**

**BKU69300**

Brand: MBA



Brushes, Electric Dental Motors & Laboratory Drives

## EB40A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**04059000**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



**Bearing for Rotor**

**04062000**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## EB40E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**04059000**

**SR1-5A-ZZS-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



**Bearing for Rotor**

**04062000**

**SMR104A-ZZ-MC3-JP6**

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



## HA43A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01672300**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



## HD43

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01381900**

**SMR85-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



## HD43EM

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01672300**

**SMR85C-ZZS-MC34-TH5-L23**

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

Brand: MBA



## LA66

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft 7-14-5mm**

**03202200**

**S687A-ZZS-MC3-GR3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)

Brand: MBA



**Bearing for Shaft 8-16-6mm**

**03202300**

**S688B-ZZ-MC3-GR3-T9H-G36**

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## LT25

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor**

**01850600**

**SR2-5A-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



## Motor MF Perfecta

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft 3-8-4mm**

**01206200**

**S693A-ZZS-MC3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)



**Bearing for Shaft**

**02685500**

**S685A-ZZS-MC34-TH7-R-G48**

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)

Brand: MBA



## Motor MF Tectorque

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Shaft 3-8-4mm**

**01206200**

**S693A-ZZS-MC3-T9H**

High Speed Bearing - Polyamide Retainer

Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)



**Bearing for Shaft**


**02685500**

**S685A-ZZS-MC34-TH7-R-G48**

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)

Brand: MBA



## RA 24

Position / Description / Application

Manufacturers Part


MBA Part

**Xenon Lamp**

**04208000**

**BKU68000Xenon**

Brand: MBA



Key




Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities




\*\* Available while stocks last

*Not in current production*


## RC95BC

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05421700</b>	<b>BKUR230046</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>05421700</b>	<b>BKUR2300464</b>	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic  Brand: myonic	<b>05421700</b>	<b>BKUR2300462</b>	


## RC95RM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05421700</b>	<b>BKUR230046</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	<b>05421700</b>	<b>BKUR2300464</b>	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic  Brand: myonic	<b>05421700</b>	<b>BKUR2300462</b>	

## RQ24

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Roto Quick Coupler O-Ring</b>  Special dental applications O-Ring  Brand: MBA	<b>02060200+02060100</b>	<b>BKU60200</b>	

## RQ34

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Roto Quick Coupler O-Ring</b>  Special dental applications O-Ring  Brand: MBA	<b>02060200+02060100</b>	<b>BKU60200</b>	


**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production


## S10

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	<b>01577000</b>	<b>S623-ZZ-MC3-JP6</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>01672300</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	



## S11

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3x10x4)	<b>01577000</b>	<b>S623-ZZ-MC3-JP6</b>	

## S9

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	<b>01577000</b>	<b>S623-ZZ-MC3-JP6</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)  Brand: MBA	<b>01672300</b>	<b>SMR85C-ZZS-MC34-TH5-L23</b>	

## Straight Handpiece 6000 RPM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Chuck End</b> Stainless Bearing Open P6 precision bearing with low-torque retainer and light oil lubricant (3.18x6.35x2.38)	<b>00352800</b>	<b>SR144-MC3-JP6-AF12</b>	
<b>Bearing for Centre</b> Stainless Bearing Double shielded P6 precision stainless steel bearing (5x11x4)	<b>02685500 ?</b>	<b>SMR115-ZZ-MC3-JP6</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Synea

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Spring Washer**

**02749408**

**BKU69408**

Brand: MBA



## Synea TA98

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SR144A-D27-ZS-TO**



**Replacement Wave Washer**

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)

2 Required

**W0048-T-007-008-W3-S2**



**Replacement O-Ring**

Very High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

*OR-00605-079-VBR75*



## Synea TA98L

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Bearing**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic

**SR144A-D27-ZS-TO**



**Replacement Wave Washer**

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)

2 Required

**W0048-T-007-008-W3-S2**



**Replacement O-Ring**

Very High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

*OR-00605-079-VBR75*



## Synea Vision

Position / Description / Application

Manufacturers Part


MBA Part

**Replacement Spring Washer**

**06876000**

**BKU66000**





Brand: MBA




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>




## TA96

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>04396700</b>	<b>BKUR230029</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>04396700</b>	<b>BKUR2300295</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>04396700</b>	<b>BKUR2300294</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>04396700</b>	<b>BKUR2300293</b>	

## TA96CM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>04378600</b>	<b>SR144A-D27-ZS-TO</b>	

## TA96L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic Hybrid Bearing - Polyamide Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer and encapsulated shield one side (3.18x6.35x2.78)  Brand: MBA	<b>02770712</b>	<b>SR144A-D27-CSN-ZE-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>02770712</b>	<b>SR144A-D27-ZS-TO</b>	
<b>Xenon Lamp</b>  Brand: MBA	<b>02207600</b>	<b>BKU67600Xenon</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TA96LC

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**04378600**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



## TA96LCM

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**04378600**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



## TA97C LED

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**


**05307200**

**BKUR230030**

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA



**Replacement Turbine With Bearings**


**05307200**

**BKUR2300303**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**05307200**

**BKUR2300305**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



**Replacement Turbine With Bearings**


**05307200**

**BKUR2300304**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## TA97CLED

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**05288500**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## TA97CLM

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**05288500**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



## TA97LW

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**06307800**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



## TA98

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**


**04396600**

**BKUR230027**

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA



**Replacement Turbine With Bearings**


**04396600**

**BKUR2300275**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



**Replacement Turbine With Bearings**


**04396600**

**BKUR2300274**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**04396600**

**BKUR2300273**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic




**Replacement Turbine With Bearings**

**04396600**

**BKUR230027T**

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## TA98C

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**


**04378600**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## TA98C LED

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05307000</b>	<b>BKUR230033</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05307000</b>	<b>BKUR2300333</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>05307000</b>	<b>BKUR2300335</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05307000</b>	<b>BKUR230033T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05307000</b>	<b>BKUR2300334</b>	





## TA98CLED

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Press Button</b>  Brand: MBA	<b>05055800</b>	<b>BKU65800</b>	





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TA98L

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic Hybrid Bearing - Polyamide Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer and encapsulated shield one side (3.18x6.35x2.78)  Brand: MBA	02770712	SR144A-D27-CSN-ZE-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	02770712	SR144A-D27-ZS-TO	
<b>Replacement Press Button</b>  Brand: MBA	03503000	BKU63000	
<b>Xenon Lamp</b>  Brand: MBA	02207600	BKU67600Xenon	





## TA98L New Version

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	05306900	BKUR230033	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	05306900	BKUR2300335	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	05306900	BKUR2300334	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	05306900	BKUR2300333	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## TC95

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>03305000</b>	<b>BKUR230028</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>03305000</b>	<b>BKUR2300283</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>03305000</b>	<b>BKUR2300285</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>03305000</b>	<b>BKUR2300284</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*




## TE95

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: MBA	<b>05541600</b>	<b>SR144A-D27-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>SG1295X</b>	<b>SR144A-D27-ZS-TO</b>	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05958900</b>	<b>BKUR230036</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>05958900</b>	<b>BKUR2300365</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05958900</b>	<b>BKUR2300364</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05958900</b>	<b>BKUR2300363</b>	



Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TE95 production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Torlon retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-CSN-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-CSN-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-ZSO-KC</b>	





## TE95BC

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: MBA	<b>05541600</b>	<b>SR144A-D27-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>SG1295X</b>	<b>SR144A-D27-ZS-TO</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TE97

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05959100-05959200</b>	<b>BKUR230045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>05959100-05959200</b>	<b>BKUR2300455</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05959100-05959200</b>	<b>BKUR2300454</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05959100-05959200</b>	<b>BKUR2300453</b>	



**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TE97 production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Torlon retainer (3.18x6.35x2.77)  Brand: MBA	07003500	SFR144A-D74F-CSN-ZSO-TO	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	07003500	SFR144A-D74F-ZSO-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	07003500	SFR144A-D74F-CSN-ZSO-KC	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	07234102	BKUR230050	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic  Brand: myonic	07234102	BKUR2300504	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic  Brand: myonic	07234102	BKUR2300502	





## TE97BC

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: MBA	05541600	SR144A-D27-ZS-KC	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	SG1295X	SR144A-D27-ZS-TO	



**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production



## TE97CLEED

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05959100-05959200</b>	<b>BKUR230045</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>05959100-05959200</b>	<b>BKUR2300455</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05959100-05959200</b>	<b>BKUR2300454</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05959100-05959200</b>	<b>BKUR2300453</b>	

## TE97CM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>05958600</b>	<b>SR144A-D27-CSN-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.78)  Brand: MBA	<b>05958600</b>	<b>SR144A-D27-CSN-ZS-TO</b>	



## TE97CN

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.78)  Brand: MBA	<b>05958600</b>	<b>SR144A-D27-CSN-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded, high precision stainless steel bearing with ceramic balls and Torlon retainer (3.18x6.35x2.78)  Brand: MBA	<b>05958600</b>	<b>SR144A-D27-CSN-ZS-TO</b>	





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TE97M

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: MBA	<b>05541600</b>	<b>SR144A-D27-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>SG1295X</b>	<b>SR144A-D27-ZS-TO</b>	




## TE98

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>05958900</b>	<b>BKUR230036</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>05958900</b>	<b>BKUR2300365</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05958900</b>	<b>BKUR2300363</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>05958900</b>	<b>BKUR2300364</b>	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts &amp; WCC coated actuators</p>




**Key** Usually stocked or Short Delivery \*\* Available while stocks last

Long Delivery or Large Minimum Quantities *Not in current production*



## TE98 300 Series production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>06787501-06787502</b>	<b>BKUR230049</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings  Brand: myonic	<b>06787501-06787502</b>	<b>BKUR2300492</b>	
<b>Turbine Sub Assembly With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings with Full stainless bearings Deep Groove  Brand: myonic	<b>06787501-06787502</b>	<b>BKUR2300492</b>	

## TE98 production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Torton retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-CSN-ZSO-TO</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-CSN-ZSO-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Single shielded precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)  Brand: MBA	<b>07003500</b>	<b>SFR144A-D74F-ZSO-KC</b>	

## TE98B

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: MBA	<b>05541600</b>	<b>SR144A-D27-ZS-KC</b>	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Phenolic Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>SG1295X</b>	<b>SR144A-D27-ZS-TO</b>	





**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## TE98CLEd

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>04396600</b>	<b>BKUR230027</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>04396600</b>	<b>BKUR2300273</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>04396600</b>	<b>BKUR2300275</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>04396600</b>	<b>BKUR230027T</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>04396600</b>	<b>BKUR2300274</b>	

## TG97 400 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>07006000</b>	<b>BKUR230038</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>07006000</b>	<b>BKUR2300385</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>07006000</b>	<b>BKUR2300383</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>07006000</b>	<b>BKUR2300384</b>	

Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## TG97L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**05958600**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

Brand: Myonic



## TG98

Position / Description / Application

Manufacturers Part

MBA Part


**Replacement Turbine With Bearings**

**07005900**

**BKUR230039T**

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



## TG98 400 Series

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**


**07005900**

**BKUR230039**

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA



**Replacement Turbine With Bearings**


**07005900**

**BKUR2300393**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**07005900**

**BKUR2300394**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**07005900**

**BKUR2300395**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



## TG98L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**05958600**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)





Brand: Myonic







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## TK100 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	07272900	BKUR230051	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	07272900	BKUR2300514	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	07272900	BKUR2300515	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	07272900	BKUR2300513	

## TK94 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	07273100	BKUR230052	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	07273100	BKUR2300525	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	07273100	BKUR2300523	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	07273100	BKUR2300524	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

## TK94L

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**07106200**

**SR144-D28-CSN-ZE-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)

Brand: MBA



## TK94LM

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**07106200**

**SR144-D28-CSN-ZE-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)

Brand: MBA



## TK97 500 Series

Position / Description / Application

Manufacturers Part

MBA Part

**Replacement Turbine Excluding Bearings**


**07272600**

**BKUR230043**

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Brand: MBA



**Replacement Turbine With Bearings**


**07272600**

**BKUR2300434**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**07272600**

**BKUR2300433**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Brand: myonic



**Replacement Turbine With Bearings**


**07272600**

**BKUR2300435**

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Brand: MBA



## TK97LM

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Rotor Cartridge**

**06307900**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)





Brand: Myonic




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


## TK98 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Turbine Excluding Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings  Brand: MBA	<b>07272800</b>	<b>BKUR230044</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>07272800</b>	<b>BKUR2300444</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic  Brand: myonic	<b>07272800</b>	<b>BKUR2300443</b>	
<b>Replacement Turbine With Bearings</b> For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer  Brand: MBA	<b>07272800</b>	<b>BKUR2300445</b>	


## TK98LM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Rotor Cartridge</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)  Brand: Myonic	<b>06307900</b>	<b>SR144A-D27-ZS-TO</b>	

## Top Air 100 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required	<b>02722300</b>	<b>SR144-D28-ZS-KC</b>	

## Top Air 200 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Bearing</b> High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required	<b>02106302</b>	<b>SR144-D28-ZS-KC</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## Top Air 300 Series

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

**02106302**

**SR144-D28-ZS-KC**

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## Top Air 700 Series

Position / Description / Application

Manufacturers Part

MBA Part

### PTFE Washer-Spacer

0.13mm - half thou - thick

Deleted Item. Available while stocks last. (3x6.56x0.13)

2 Required

**W0030-F-007-001-PT**



### Replacement Bearing

**02722300**

**SR144A-D27-ZS-TO**

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic



### Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)

**W0048-T-007-008-W3-S2**



### Replacement O-Ring

Very High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

*OR-00605-079-VBR75*



Key





Usually stocked or Short Delivery

\*\* Available while stocks last





Long Delivery or Large Minimum Quantities

*Not in current production*

## Top Air 800 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required Brand: Myonic	<b>02722300</b>	<b>SR144A-D27-ZS-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00605-079-VBR75</b>	





## Top Lite 700 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required Brand: Myonic	<b>02722300</b>	<b>SR144A-D27-ZS-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00605-079-VBR75</b>	

**Key** Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## Top Lite 800 Series

Position / Description / Application	Manufacturers Part	MBA Part	
<b>PTFE Washer-Spacer</b> 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		<b>W0030-F-007-001-PT</b>	
<b>Replacement Bearing</b> High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required Brand: Myonic	<b>02722300</b>	<b>SR144A-D27-ZS-TO</b>	
<b>Replacement Wave Washer</b> 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 2719 (4.76x6.35x0.76)		<b>W0048-T-007-008-W3-S2</b>	
<b>Replacement O-Ring</b> Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<b>OR-00605-079-VBR75</b>	

## Toplight

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Spring Washer</b>  Brand: MBA	<b>02749408</b>	<b>BKU69408</b>	

## Trend

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Replacement Spring Washer</b>  Brand: MBA	<b>02749408</b>	<b>BKU69408</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## WA56

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	<b>01589400</b>	<b>S683-MC3-JP6-AF12</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	<b>01577000</b>	<b>S623-ZZ-MC3-JP6</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)  Brand: Myonic	<b>01589600</b>	<b>SMR74-MC34-GR7-R-L23</b>	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>03819000</b>	<b>BKU69000</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>05431100</b>	<b>BKU61100</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>03821000</b>	<b>BKU61000</b>	 <small>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</small>
<b>Replacement Press Button</b>  Brand: MBA	<b>03819500</b>	<b>BKU69500</b>	 <small>Press Button cpl. for Heads &amp; Air Driven Handpieces</small>
<b>WA56LT NEW</b>			
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)  Brand: Myonic	<b>05039400</b>	<b>SMR74D-MC3-GR5-L23</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## WA56LT OLD

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**01579300**

**SMR55-0235-MC3-GR3-AF12**

High Precision Bearing  
Open Abec 3 precision bearing (2.35x5.5x2)



## WA62LTA

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**01579300**

**SMR55-0235-MC3-GR3-AF12**

High Precision Bearing  
Open Abec 3 precision bearing (2.35x5.5x2)



## WA66

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing  
Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



**Head Gear Complete 1:1**

**03819000**

**BKU69000**

Brand: MBA



## WA66LTA

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**01579300**

**SMR55-0235-MC3-GR3-AF12**

High Precision Bearing  
Open Abec 3 precision bearing (2.35x5.5x2)



## WA67LTA

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**01579300**

**SMR55-0235-MC3-GR3-AF12**

High Precision Bearing  
Open Abec 3 precision bearing (2.35x5.5x2)



## WA86

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing  
Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

## WA956

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## WA966

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## WA969

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Middle Gear**

**01579600**

**SMR55-0235-T4HB-MC3-AF12**

High Speed Bearing - Torlon Retainer

Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)

Brand: MBA



## WA975

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## WA986

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)

Brand: Myonic



## WA988

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**01589600**

**SMR74-MC34-GR7-R-L23**

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



Brand: Myonic






Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



## WA99

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Speed Bearing - Torlon Retainer Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)  Brand: MBA	<b>01579600</b>	<b>SMR55-0235-T4HB-MC3-AF12</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	<b>01565800</b>	<b>SMR85-MC3-JP6-AF12</b>	

## WA99A

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)  Brand: MBA	<b>05257200</b>	<b>SMR6.50-3.40-CSN-ZZS-MC3-TO</b>	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)  Brand: MBA	<b>04486000</b>	<b>SMR74-CSN-TH7-TO</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	<b>01565800</b>	<b>SMR85-MC3-JP6-AF12</b>	

## WA99LT

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)  Brand: MBA	<b>05257200</b>	<b>SMR6.50-3.40-CSN-ZZS-MC3-TO</b>	
<b>Head Gear Complete</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: MBA	<b>05220300</b>	<b>BKU60300</b>	

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
*Not in current production*

## WD56

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01577000**

**S623-ZZ-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3x10x4)



**Head Gear Complete 1:1**

**03425000**

**BKU65000**

Brand: MBA



**Replacement Middle Gear Complete 1:1**

**03403500**

**BKU63500**

Brand: MBA



## WD66

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01589400**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## WD969

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01589400**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



## WD999

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01589400**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## WE56

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	<b>05765800</b>	<b>SMR63-MC3-JP6-AF12</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	<b>05060600</b>	<b>SR2-5A-ZZ-MC3-JP6</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)	<b>06164000</b>	<b>S684-MC3-JP6-AF12</b>	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>05843300</b>	<b>BKU63300</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>05852500</b>	<b>BKU62500W</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Press Button</b>  Brand: MBA	<b>05843800</b>	<b>BKU63800</b>	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>

## WE56A

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	<b>05060600</b>	<b>SR2-5A-ZZ-MC3-JP6</b>	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>05843300</b>	<b>BKU63300</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces &amp; Heads</p>

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## WE56LEDG

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 1:1**

**05844400**

**BKU64400**

Brand: MBA



## WE66

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**05765800**

**SMR63-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)




**Head Gear Complete 1:1**

**05843300**

**BKU63300**

Brand: MBA



## WE99A

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**02649400**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



## WE99LEDG

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**05765800**

**SMR63-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



**Bearing for Headgear**

**02649400**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

\*\* Available while stocks last

*Not in current production*

## WG56

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 1:1</b>  Brand: MBA	<b>07605900</b>	<b>BKU65900WG</b>	
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>05233200</b>	<b>BKU63200</b>	
<b>Replacement Middle Gear Complete 1:1</b>  Brand: MBA	<b>06796000</b>	<b>BKU69600</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>4 Parts</b>	<b>BKU62800</b>	

## WG99A

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete - Premium</b> For Handpiece with Wedge Type Chuck Premium head gear assembled with ceramic ball bearings tested in master electric handpiece  Brand: MBA	<b>07473400</b>	<b>BKU63400</b>	
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: myonic	<b>07473400</b>	<b>BKU63400M</b>	




## WG99LT

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete - Premium</b> For Handpiece with Wedge Type Chuck Premium head gear assembled with ceramic ball bearings tested in master electric handpiece  Brand: MBA	<b>07473400</b>	<b>BKU63400</b>	
<b>Head Gear Complete - Standard</b> For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece  Brand: myonic	<b>07473400</b>	<b>BKU63400M</b>	


Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production




## WI75

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Head Gear Complete 20:1</b>  Brand: MBA	<b>04939900</b>	<b>BKU69900</b>	
<b>Replacement Middle Gear Complete 20:1</b>  Brand: MBA	<b>04667600</b>	<b>BKU67600</b>	
<b>Replacement Press Button</b>  Brand: MBA	<b>04940500</b>	<b>BKU60500</b>	

## WI75 LEDG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	<b>04033500</b>	<b>SMR63-MC3-JP6-AF12</b>	

## WI75LEDG

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Precision Bearing Open Abec 3 precision bearing (2.35x5.5x2)	<b>04733900</b>	<b>SMR55-0235-MC3-GR3-AF12</b>	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)  Brand: Myonic	<b>04027300</b>	<b>SMR74-MC34-GR7-R-L23</b>	
<b>Bearing for Lower Shaft</b> High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)	<b>06164000</b>	<b>S684-MC3-JP6-AF12</b>	

<b>Key</b>	<b>Usually stocked or Short Delivery</b>
	<b>** Available while stocks last</b>

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

## WK56

Position / Description / Application

Manufacturers Part

MBA Part

**Head Gear Complete 1:1**

**05431100**

**BKU61100**

Brand: MBA



**Replacement Middle Gear Complete 1:1**

**05233200**

**BKU63200**

Brand: MBA



**Replacement Middle Gear Complete 1:1**

**06796000**

**BKU69600**

Brand: MBA



## WK93LT

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**06558000**

**SR166-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4.76x9.52x3.18)



## WK99LT

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**05257200**

**SMR6.50-3.40-CSN-ZZS-MC3-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)

Brand: MBA



**Bearing for Lower Shaft**

**06558000**

**SR166-ZZS-MC3-JP6**

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4.76x9.52x3.18)



## WP64

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Lower Shaft**

**01589400**

**S683-MC3-JP6-AF12**

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Key





Usually stocked or Short Delivery

\*\* Available while stocks last




Long Delivery or Large Minimum Quantities

*Not in current production*


## WS75

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Lower Shaft</b> High Precision Bearing Double shielded P6 precision stainless steel bearing (3x10x4)	<b>01577000</b>	<b>S623-ZZ-MC3-JP6</b>	
<b>Head Gear Complete 20:1</b>  Brand: MBA	<b>04939900</b>	<b>BKU69900</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Middle Gear Complete 20:1</b>  Brand: MBA	<b>06025100</b>	<b>BKU65100</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Press Button</b>  Brand: MBA	<b>04940500</b>	<b>BKU60500</b>	 <p>Press Button cpl. for Heads &amp; Air Driven Handpieces</p>

## WS75 EKM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Middle Gear</b> High Precision Bearing Open Abec 3 precision bearing (2.35x5.5x2)	<b>04733900</b>	<b>SMR55-0235-MC3-GR3-AF12</b>	
<b>Head Gear Complete 20:1</b>  Brand: MBA	<b>05101300</b>	<b>BKU61300</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>
<b>Replacement Middle Gear Complete 20:1</b>  Brand: MBA	<b>04016900</b>	<b>BKU66900</b>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces &amp; Heads</p>

## WS75E KM

Position / Description / Application	Manufacturers Part	MBA Part	
<b>Bearing for Headgear</b> High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	<b>04027300</b>	<b>SMR74-MC34-GR7-R-L23</b>	
Brand: Myonic			

Key Usually stocked or Short Delivery  
\*\* Available while stocks last

Long Delivery or Large Minimum Quantities  
Not in current production

## WS92E

Position / Description / Application

Manufacturers Part

MBA Part

**Bearing for Headgear**

**04025900**

**SMR74-CSN-TH7-TO**

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)

Brand: MBA



Key

Usually stocked or Short Delivery

\*\* Available while stocks last

Long Delivery or Large Minimum Quantities

*Not in current production*

## Friction Grip

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144-D28-ZS-KC

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



## Push Button

Position / Description / Application

Manufacturers Part

MBA Part

### Replacement Bearing

SR144A-D27-ZS-TO

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

Brand: Myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>





Miniature Bearings Australia  
<https://minibearings.com.au>  
[sales@minibearings.com.au](mailto:sales@minibearings.com.au)

We sell an extensive range of products with varying amounts of technical data. Please assist us to provide you with accurate information by placing your orders and enquiries online or by e-mail.

## Dental Handpiece Replacement Parts