

MBA

2026.1

Dental Handpiece Replacement Parts



Make	Pages	Make	Pages
ABU®	3	KMD	115
ACETON SATELEC®	4	Kool Torque®	116
ADEC®	5	Lares®	117-118
ALPHA AIR®	6	Marathon®	119-120
B.A. International®	7-11	Micro Mega®	121
BANNER®	12	Midwest®	122-125
BEAVER®	13	Miscellaneous	126
Bien Air®	14-22	MK-Dent®	127-136
Borden®	23	Morita®	137
Brasseler®	24	MTC	138
Bravo Bedetrice®	25	MTI	139
Buffalo®	26	Muss®	140
C. De Giorgi®	27	NSK®	141-166
Castellini®	28-29	Ritter®	167
Chirana®	30	Sable®	168-169
COXO®	31	Saeshin®	170-171
Dabi Atlante®	32-33	Schick®	172-173
Densim®	34-36	Siemens®	174-182
DME	37	Sirona®	183-190
Dynadent®	38	SS White®	191
Dynasurge®	39	Star®	192-194
Elco® - Austria	40	Starflite®	195
Falcon®	41	Sternweber®	196
Fangsco®	42-43	TKD	197
Faro®	44-45	Torx®	198
Futura®	46	W&H®	199-219
Gamborini®	47	Yoshida®	220
Henry Schein®	48-49		
Impact-Air®	50		
Jalenco®	51		
KaVo®	52-114		



To find and order parts online, simply type the part number from this catalogue into the search box on any web page.

Unless specified otherwise, all parts are generic offerings to replace the original parts. Some items are offered with alternative materials. Please check carefully before placing an order.

350000 RPM

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

ACETON SATELEC®

F12200 Airscalers

Position / Description / Application

Manufacturers Part

MBA Part

Bush / Rotor Repair Kit

707091

BKU67091



Key



Usually stocked or Short Delivery

** Available while stocks last



Long Delivery or Large Minimum Quantities

Not in current production

5000

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	

5100

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	





Key

Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

B.A. International®


13Z

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece		BKU65687	
Head Gear Complete - Standard For Handpiece with Wedge Type Chuck Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece		BKU65587X	
Head Gear Complete For Handpiece with Wedge Type Chuck Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece		BKU65587	
Head Gear Complete - Standard For Handpiece with Wedge Type Chuck Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece		BKU65687X	

BA 200 LTS

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit Dental replacement drive kit	6175587 + 6175603	BKU68703	

BA 40 LSS Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit Dental replacement drive kit	6405695 + 6000959	BKU69559	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

BA200LTS

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

BKU65687M



Head Gear Complete

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

BKU65687X



Head Gear Complete

Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece


BKU65587X



Head Gear Complete

Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece


BKU65687



Head Gear Complete

Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece

BKU65587



BA40

Position / Description / Application

Manufacturers Part


MBA Part

Drive Kit

Dental replacement drive kit

C875060 + C880022

BKU66022



BA45

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

Dental replacement drive kit

C875060 + C880022

BKU66022



BA45 LS

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

Dental replacement drive kit

6000983 + 6000959

BKU68359



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

B.A. International®

Diameter 8.7 BA200LTS

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece

BKU65687M



Optima® BA525

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine Excluding Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

BKUR310000




Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKUR3100003




Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

BKUR3100004




Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

BKUR3100005




Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKURC310000



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Optima® BA528

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine Excluding Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings

BKUR310001



Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKUR3100013



Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic

BKUR3100014



Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer


BKUR3100015



Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKURC310001



Key **Usually stocked or Short Delivery**
** Available while stocks last





Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Ultima® 650L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
Bearing for Cap End High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	
Bearing for Cap End High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	
Bearing for Cap End High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	

Ultimate 635





Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200244	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200243	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200245	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220024	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**





Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Ultimate 670

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220022	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200225	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200224	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200223	

Ultimate 675





Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200235	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200234	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200233	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220023	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last





Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Ultimate 680

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200235	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200234	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200233	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220023	

Ultimate 685





Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220023	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200233	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200235	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200234	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last





Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Ultimate 755

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220024	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200245	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200243	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200244	

Ultimate 758

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR220021	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR2200215	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200214	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR2200213	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

B.A. International®

Ultimate BA 250LT Mini

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6414457 + 6175603

BKU65703



Ultimate BA688 New Type

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKUR2200363



Ultimate BA755 New Type

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

BKUR2200353



Ultimate Power

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6175587 + 6175603

BKU68703

Dental replacement drive kit



Drive Kit

6405695 + 6000959

BKU69559

Dental replacement drive kit



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Banner

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Slow Speed 40000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x7x3)

2 Required

S683A-ZZ-MC3-JP6



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Black Pearl®

Position / Description / Application	Manufacturers Part	MBA Part	
Thrust Washer Stainless Steel	009.34.58	BKU63458	
Replacement Press Button	009.40.25+009.40.30+009.40.24+009.40.23	BKU64025	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1500164	BKUR2400363	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer	1500164	BKUR2400365	
Replacement Turbine Excluding Bearings Turbine sub assembly excluding bearings	1500164	BKUR240036	
Turbine Sub Assembly With Bearings With Full stainless bearings, Deep Groove, Phenolic retainer, myonic	1500164	BKUR2400364	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-KC	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-TO	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental ceramic-hybrid bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-CSN-ZS-KC	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

Black Pearl® Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

705.01.81

OR-00804-100-N70

High Temperature

Deleted Item. Available while stocks last.



Black Pearl® Friction Grip

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

701.01.00

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Key

Usually stocked or Short Delivery





** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

BIEN AIR®



Bora®

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Phenolic retainer, myonic	1500181	BKUR2400374	
Thrust Washer Stainless Steel	009.34.58	BKU63458	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer	1500181	BKUR2400375	
LED Bulb	204.33.07	BKU63307LED	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1500181	BKUR2400373	
Replacement Press Button	009.34.46+209.34.11+009.40.24+009.34.58	BKU63446	
Xenon Lamp	204.33.07 (1500000)	BKU63307Xenon	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1500181	BKUR240037	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.00	SR2-5A-D17-ZS-KC	


Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental ceramic-hybrid bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-CSN-ZS-KC	 Image Coming Soon
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1500181	BKURMU240037	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1500181	BKURM240037	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.00	SR2-5A-D17-ZS-KC	

Bora® Air Driven Handpiece




Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last.	705.01.81	OR-00804-100-N70	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last




Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®


Bora® lina TU

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer	1500164	BKUR2400365	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Phenolic retainer, myonic	1500164	BKUR2400364	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1500164	BKUR240036	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1500164	BKUR2400363	

CA 1:1

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1	1500549 / 1501085	BKU61085	
Head Gear Complete 1:1	1500289 / 1501083	BKU60289	
Replacement Press Button	1500550	BKU60550	

CA 1:1 L

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	1501085 + 1501083	BKU68583	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

CA 1:1 NOVA

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

1503014

BKU63014



CA 1:5

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

1500169


BKU60169



Replacement Middle Gear Complete 1:5

1500057

BKU60057




Head Gear Complete - Premium

1501933

BKU61933

For Handpiece with Wedge Type Chuck

Premium head gear assembled with ceramic ball bearings and tested in master electric handpiece




Head Gear Complete - Standard

1501933

BKU61933X

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings and tested in master electric handpiece



Bearing for Headgear

701.01.68

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



Bearing for Middle Gear

701.01.75

SMR4.762-2.35A-ZZS-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)



CA 1:5 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

1501933 + 1500057

BKU63357



Key

Usually stocked or Short Delivery




Long Delivery or Large Minimum Quantities

** Available while stocks last


Not in current production

BIEN AIR®





CA 20:1

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 20.1	1501026	BKU61026	
Replacement Press Button	1501027	BKU61027	
Head Gear Complete 20.1	1501898	BKU61898	

CA 20:1 L

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	1501898 + 1501026	BKU69826	

Croma M




Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	209.30.19	BKUR2400395	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	209.30.19	BKUR2400394	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	209.30.19	BKUR2400393	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	209.30.19	BKUR240039	

Key **Usually stocked or Short Delivery**
** Available while stocks last


Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®




HP11

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	701.01.65	S623-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	701.01.53	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	701.01.53	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon


LaBora® tory

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.00	SR2-5A-D17-ZS-KC	

MC2

Position / Description / Application	Manufacturers Part	MBA Part	
LED Bulb	211.60.16	BKU66016LED	 <small>LED Bulb, Couplers & Dental Motors</small>
Replacement Motor Brush	211.45.06	BKU64506	 <small>Brushes, Electric Dental Motors & Laboratory Drives</small>
Xenon Lamp	211.60.16 (1500007)	BKU66016Xenon	 <small>Xenon Lamp for Dental Couplers & Dental Motors</small>

MC2 Intra-Matic

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required	705.01.74	OR-00800-070-N70	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
 Not in current production

MC2 Isolite 300

Position / Description / Application

Manufacturers Part


MBA Part

Dental Motor Set with 15 parts

1502219-001

BKU62219

Includes parts - 3x 1300967 - 1x 705.03.15 - 1x 652.02.02 - 1x 701-01-14 - 1x 701.01.23 - 2x 705.02.07 - 2x 211.60.13 - 2x 211.45.06 - 1x 011.35.37 - 1x 1500391



MC2 Isolite100

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

701.01.23

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



MC2 Isolite300

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

701.01.23

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



Key

Usually stocked or Short Delivery




** Available while stocks last

Long Delivery or Large Minimum Quantities


Not in current production

BIEN AIR®


MC3

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Motor Brush	211.75.08	BKU67508	
Bearing for Rotor High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)	701.01.70	SR1-5A-ZZS-MC3-JP6	
Bearing for Rotor High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	701.01.71	SMR104A-ZZ-MC3-JP6	
Dental Motor Set with 21 parts Includes parts - 3x 01862300 - 1x 04062500 - 1x 04061800 - 1x 04062300 - 1x 04058800 - 2x 04062100 - 2x 04062200 - 2x 04062400 - 2x 04058800 - 3x 04059200 - 2x 04059300 - 1x 04069600	1502217-001	BKU62217	
LED Bulb	211.60.16	BKU66016LED	
Xenon Lamp	211.60.16 (1500007)	BKU66016Xenon	

Micromotor Aquilon 830

Position / Description / Application	Manufacturers Part	MBA Part	
Front and Brush End Bearings High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.14	SR2-5A-D17-ZS-KC	






Micromotor Aquilon 850

Position / Description / Application	Manufacturers Part	MBA Part	
Front and Brush End Bearings High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57) 2 Required	701.01.14	SR2-5A-D17-ZS-KC	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production





Micromotor Isolite 100 and 300

Position / Description / Application	Manufacturers Part	MBA Part	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	
Micromotor MC2 20			
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	





Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Micromotor MC2 20GT

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	





Micromotor MC2 20GTAV

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	





Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Micromotor MC2 40E

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		<i>OR-00386-091-N70</i>	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	





Micromotor MC2 40GT

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		<i>OR-00386-091-N70</i>	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


Micromotor MC2 40GTAV

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Large High Temperature Deleted Item. Available while stocks last.		OR-00804-100-N70	
Replacement O-Ring Medium Special sized o-ring used in dental Applications		OR-00400-100-N70	
Replacement O-Ring Small High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications		OR-00386-091-N70	
Front Bearing High Speed Bearing - Phenolic Retainer Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)	701.01.14	SR2-5A-D17-ZS-KC	


Motor Isolite 300

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	

Motor Isolite100

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	

Motor MC2



Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)	701.01.14	SR2-5A-ZZ-MC3-JP6	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

MX

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3.97x7.94x3.18)	701.01.94	SR155A-ZZS-MC3-JP6	
Bearing for Rotor High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4.76x7.94x3.18)	701.01.93	SR156A-ZZS-MC3-JP6	




Ondine

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	209.30.19	BKUR2400393	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required		SR144A-D27-ZS-T9H-OPT	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	209.30.19	BKUR2400394	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	209.30.19	BKUR2400395	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	209.30.19	BKUR240039	


Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


PA11

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	701.01.65	S623-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	701.01.53	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	701.01.53	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon

PM AM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Middle Gear High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)	701.01.77	SR133A-ZZS-MC3-JP6	

PM RM 1:1

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Middle Gear High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)	701.01.77	SR133A-ZZS-MC3-JP6	

Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

BIEN AIR®

Prestige





Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	009.30.82+209.31.34+009.30.63+009.31.45	BKU63082	
Xenon Lamp	204.33.07 (1500000)	BKU63307Xenon	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required		SR144A-D27-ZS-T9H-OPT	
LED Bulb	204.33.07	BKU63307LED	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)	SG1329X-D	SR144A-D31-CSN-WZ-TO	

Key **Usually stocked or Short Delivery**
** Available while stocks last




Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

Prestige S32L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1500485	BKUR2400404	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1500485	BKUR240040	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1500485	BKUR2400403	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1500485	BKUR2400405	

TD783




Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Rotor / Cartridge	208.55.01	BKUR240041	
Replacement Rotor / Cartridge	208.55.01	BKUR2400412	
Replacement Rotor / Cartridge	208.55.01	BKUR2400411	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BIEN AIR®

TDS890

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Rotor / Cartridge	208.55.01	BKUR2400411	
Replacement Rotor / Cartridge	208.55.01	BKUR2400412	
Replacement Rotor / Cartridge	208.55.01	BKUR240041	

Tornado

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1501968	BKUR2400424	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	SG1302X-D	SR144A-D27-ZS-T9H-OPT	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1501968	BKUR2400423	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1501968	BKUR2400425	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1501968	BKUR240042	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Unifix 4V Coupler

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

1300380+1300391

BKU60391

Special O-Ring for Dental Applications

Special dental applications O-Ring



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

60

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



70

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Airotor Silent Cartridge

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Silencer

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		<i>SR144-D106-ZS</i>	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Brio 440 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer		BKUR3100005	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings		BKUR310000	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKUR3100003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic		BKUR3100004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic		BKURC310000	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Bravo Bedetrice®

Drill 8 12000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

Stainless Bearing

Double shielded P6 precision stainless steel bearing (5x11x5)

2 Required

S685A-ZZ-MC3-JP6



Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

EM3000 Nose Cone

Position / Description / Application

Manufacturers Part

MBA Part

Bearing Motor End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (4x13x5)

S624-ZZ-MC3-JP6



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Goldspeed EVO M5 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

97230612 + 97230601

BKU61201



Goldspeed EVO R20 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

97230522 + 97230527

BKU62227



Goldspeed EVO S1 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

97230630 + 97230623

BKU63023



Goldspeed S1 L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

FR400N103 + FR400N104





BKU60304







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

High Power M

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	L0000874	BKUR2701033	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	L0000874	BKUR270103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	L0000874	BKUR2701034	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	L0000874	BKUR2701035	

Silent Power

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	L0001857	BKUR270105	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	L0001857	BKUR2701053	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	L0001857	BKUR2701054	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	L0001857	BKUR2701055	

Key

Usually stocked or Short Delivery





** Available while stocks last

Long Delivery or Large Minimum Quantities





Not in current production

Castellini®

Silent Power Evo

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	95529230	BKUR2701074	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	95529230	BKUR270107	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	95529230	BKUR2701073	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	95529230	BKUR2701075	

Titanium Gold 2

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	L0001370	BKUR270104	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	L0001370	BKUR2701043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	L0001370	BKUR2701044	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	L0001370	BKUR2701045	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TG 656

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	T031 2560 4/1	BKUR320003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	T031 2560 4/1	BKUR320004	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	T031 2560 4/1	BKUR320000	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	T031 2560 4/1	BKUR320005	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	T031 2560 4/1	BKURC320001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	T031 2560 4/1	BKURC3200015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	T031 2560 4/1	BKURC3200013	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

CX207 Head

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

210-G SP

BKUR350001

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



CX207 Standard Head

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

210-G SP

BKUR350001

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



CX207 Torque Head

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

210-G TP

BKUR350002

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Dabi Atlante®

L350

Position / Description / Application	Manufacturers Part	MBA Part	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	 Image Coming Soon
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	 Image Coming Soon
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	 Image Coming Soon
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	 Image Coming Soon

Key **Usually stocked or Short Delivery**
**** Available while stocks last**



Long Delivery or Large Minimum Quantities
Not in current production

Dabi Atlante®

MS350

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	

RS350

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Dublin DSLCH11

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete Dental Head 1:1
For Handpiece with Latch Type Chuck

DSSP2005

BKU62005



Intermediate Piece Dental Head 1:1

DSSP2024

BKU62024



Monaco Lux

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

DSTC2034

BKUR3301083

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Turbine Sub Assembly Without Bearings

DSTC2034

BKUR330108

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

DSTC2034

BKUR3301085

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

DSTC2034

BKURC3301085

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



Replacement Turbine With Bearings

DSTC2034

BKURC3301083

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Monaco Lux Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2035	BKUR3301095	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	DSTC2035	BKUR330109	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2035	BKUR3301093	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2035	BKURC3301093	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	DSTC2035	BKURC330109	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2035	BKURC3301095	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Monaco Lux Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2034	BKUR3301083	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2034	BKUR3301085	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2034	BKURC3301085	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	DSTC2034	BKURC330108	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2034	BKURC3301083	




Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Montreal Lux

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2041	BKUR3301105	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	DSTC2041	BKUR330110	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2041	BKUR3301103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2041	BKURC3301105	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2041	BKURC3301103	

Montreal Lux Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2042	BKURC3301205	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2042	BKURC3301203	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	DSTC2042	BKURC330120	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Montreal Lux Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2041	BKUR3301103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2041	BKUR3301105	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	DSTC2041	BKUR330110	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	DSTC2034	BKUR330108	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	DSTC2041	BKURC330110	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2041	BKURC3301103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2041	BKURC3301105	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Vienna



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2034	BKUR3301083	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	DSTC2034	BKUR330108	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2034	BKUR3301085	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	DSTC2034	BKURC3301085	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	DSTC2034	BKURC3301083	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	DSTC2034	BKURC330108	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**
 ** Available while stocks last



Long Delivery or Large Minimum Quantities
Not in current production

DME

Aden

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	

True Speed

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	

Key

Usually stocked or Short Delivery
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
 Not in current production

Hi-Torq Surgical Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

Stainless Flanged High Speed Bearing

Single shielded flanged Abec 7 precision stainless steel bearing with high speed phenolic retainer (4.76x9.52x3.18)

SFR166-ZS-MC5-GR7-KB-J260



Bearing for Cap End

Stainless Flanged High Speed Bearing

Single shielded flanged Abec 7 precision stainless steel bearing with high speed phenolic retainer (4.76x9.52x3.18)

SFR166-ZS-MC5-GR7-KB-J260



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Elco® - Austria

Elco - Austria

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

BB 7-14-5

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

11 - 4000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Precision Bearing

Deleted Item. Available while stocks last. (6x10x3)

2 Required

SMR106A-ZZS-MC34-P4P-R-AF12



Replacement Bearing

High Precision Bearing

Deleted Item. Available while stocks last. (6x10x3)

2 Required

SMR106A-ZZS-MC34-GR7-R-G48



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

A5 Small Head

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	CCHPZ4	BKUR370001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPZ4	BKUR3700013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	CCHPZ4	BKUR3700014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPZ4	BKUR3700015	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	CCHPZ4	BKURC370001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPZ4	BKURC3700013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPZ4	BKURC3700015	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

A6 Head

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKUR3700024	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKUR3700023	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	CCHPC4	BKUR370002	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPC4	BKUR3700025	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKURC3700023	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPC4	BKURC3700025	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

A6 Standard Head

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKUR3700024	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPC4	BKUR3700025	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKUR3700023	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	CCHPC4	BKUR370002	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	CCHPC4	BKURC3700023	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	CCHPC4	BKURC3700025	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	CCHPC4	BKURC370002	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

BTC77

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-T9H-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-KC-OPT



Delta

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

SR144A-D27-ZS-T9H-OPT



LC81

Position / Description / Application

Manufacturers Part

MBA Part

Bearings - Small 3.175-9.525-3.969mm

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3.18x9.52x3.97)

2 Required

SR2-ZZ-MC3-JP6



Medium Bearing Replacement

BB 7.938-12.7-3.969

SR1810-ZZS-MC3

Stainless Bearing

Double shielded standard stainless steel bearing (7.94x12.7x3.97)



Large Bearing Replacement 8-16-6mm

BB 8-16-6

S688B-ZZ-MC3-GR3-T9H-G36

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



MC01R

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

BKU64040M




Head Gear Complete - Standard

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

BKU64040X



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MC01RL

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Standard

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

BKU64040X



Head Gear Complete - Premium

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

BKU64040M



Head Gear Complete

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece

BKU64040



S380

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-ZS-KC-OPT



S400

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

SR144A-D27-ZS-T9H-OPT



S405

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

SR144A-D27-ZS-T9H-OPT



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Star

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Key




Usually stocked or Short Delivery

** Available while stocks last




Long Delivery or Large Minimum Quantities

Not in current production

MKQ

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		S687A-ZZ-MC3-JP6	
Front Armature Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		S694-ZZ-MC3-JP6	
Bearing for Chuck End Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		S695-ZZ-MC3-JP6	




MKS

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		S687A-ZZ-MC3-JP6	
Front Armature Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		S694-ZZ-MC3-JP6	
Bearing for Chuck End Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		S695-ZZ-MC3-JP6	

Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

MKT

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End Stainless Bearing Double shielded P6 precision stainless steel bearing (7x14x5)		S687A-ZZ-MC3-JP6	
Bearing for Chuck End Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		S695-ZZ-MC3-JP6	
Front Armature Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)		S694-ZZ-MC3-JP6	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

35 II

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Armature

101-26

BKU61103



35 Z

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

101-26

BKU61103



L35 Z

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

101-26

BKU61103



Maxima® ELITE

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

570-0136

BKUR3100013

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

570-0136

BKUR310001

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

570-0136

BKUR3100015

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

570-0136

BKUR3100014

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

570-0136

BKURC310001

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Henry Schein®

Maxima® ELITE L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	570-0136	BKUR310001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	570-0136	BKUR3100014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	570-0136	BKUR3100015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	570-0136	BKUR3100013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	570-0136	BKURC310001	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Henry Schein®

Maxima® ELITE Mini

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

570-0137

BKUR3100003

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

570-0137

BKUR310000

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

570-0137

BKUR3100005

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

570-0137

BKUR3100004

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

570-0137

BKURC310000

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last




Not in current production

Henry Schein®

Maxima® ELITE Mini L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	570-0137	BKUR310003	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	570-0137	BKUR310000	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	570-0137	BKUR3100005	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	570-0137	BKUR3100004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	570-0137	BKURC310000	

Maxima® PRO 45

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	112-6773	BKURC2100223	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	112-6773	BKURC210022	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	112-6773	BKURC2100225	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Henry Schein®

Maxima® PRO 45L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

112-6773

BKURC2100225

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer

Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators



Replacement Turbine Excluding Bearings

112-6773

BKURC210022

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings

Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators



Replacement Turbine With Bearings

112-6773

BKURC2100223

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic

Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Impact-Air®

High Speed Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00607-079-N70



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-KC-OPT



Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-T9H-OPT



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Lab Handpiece Air

Position / Description / Application

Manufacturers Part

MBA Part

Both Turbine Bearings

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

0768LHC

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0139

SMR104A-ZZ-MC3-JP6

High Precision Bearing
Double shielded P6 precision stainless steel bearing (4x10x4)



10B

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0126

S685A-ZZS-MC34-TH7-R-G48

High Precision Bearing
Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)



10C

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

220-0219 / 0126

S685A-ZZ-MC3-JP6

Stainless Bearing
Double shielded P6 precision stainless steel bearing (5x11x5)
2 Required



10CHC

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

0.524.2942

BKU62942

True running accuracy less than 0.03, bore diameter 1.6mm



10CN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Cap End

0.220.1302

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option
Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)



Bearing for Cap End

0.220.1302

SMR85C-ZZS-MC34-GR5-L23





High Precision Bearing Standard Cage Option
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

10E

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-TR-L23	
Replacement Chuck True running accuracy less than 0.03, bore diameter 1.6mm	0.524.2942	BKU62942	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (4x9x4)	0.220.1301	S684A-ZZS-MC34-GR7-R-G48	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

10LH

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-TR-L23	
Replacement Roll Pin	0.524.2742	BKU62742	
Replacement Pin	0.524.3262	BKU63262	
Replacement Spindle	0.524.1421	BKU61421	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (4x9x4)	0.220.1301	S684A-ZZS-MC34-GR7-R-G48	
Replacement Chuck True running accuracy less than 0.03, bore diameter 1.6mm	0.524.2942	BKU62942	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

10LP

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-TR-L23	
Replacement Pin	0.524.3262	BKU63262	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Replacement Chuck True running accuracy less than 0.03, bore diameter 1.6mm	0.524.2942	BKU62942	
Replacement Roll Pin	0.524.2742	BKU62742	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Replacement Spindle	0.524.1421	BKU61421	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (4x9x4)	0.220.1301	S684A-ZZS-MC34-GR7-R-G48	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-L23	 <p style="text-align: center;">Image Coming Soon</p>

11B

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)	0.220.0126	S685A-ZZS-MC34-TH7-R-G48	 <p style="text-align: center;">Image Coming Soon</p>

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

121

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0021

S683-GR7-Y

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)



121A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0021

S683-GR7-Y

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)



13A

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Chuck End

220-0106

SMR85-MC3-JP6-AF12

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



Bearing for Centre

220-0024 / 0103 / 1301

S684A-ZZ-MC3-JP6

4x9x4 Stainless Steel Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



Bearing for Shaft

0.220.0103

S684A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



Bearing for Cap End

220-0068

S624-MC3-JP6-AF12

Stainless Bearing

Deleted Item. Available while stocks last. (4x13x5)



13B

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0126

S685A-ZZS-MC34-TH7-R-G48

High Precision Bearing Shielded Option

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)



Bearing for Shaft

0.220.0139

SMR104A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

14E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.2005

SMR74-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



Bearing for Headgear

0.220.2004

S683B-ZZS-MC3-GR3

High Precision Bearing Shielded Option

Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)



Bearing for Shaft

0.220.0098

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



14ES

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.2004

S683B-ZZS-MC3-GR3

High Precision Bearing

Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)



15CN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.1306

SMR9525-050-MC3-L23

High Precision Bearing

Double shielded standard bearing (5x9.52x3.18)



15G

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0126

S685A-ZZS-MC34-TH7-R-G48

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)



181CB

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0097

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

181DB

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0097

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



181L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0021

S683-GR7-Y

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)



186B

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



187

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

188

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



**Image
Coming
Soon**

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



**Image
Coming
Soon**

190

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



**Image
Coming
Soon**

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



**Image
Coming
Soon**

191

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



**Image
Coming
Soon**

192

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



**Image
Coming
Soon**

Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



**Image
Coming
Soon**

Key

Usually stocked or Short Delivery

**** Available while stocks last**

Long Delivery or Large Minimum Quantities

Not in current production

193

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



196

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

0.220.0089

SMR104A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)



Bearing for Rotor

0.220.0136

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



19E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0175

S683-GR7-Y

High Precision Bearing

Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)



Bearing for Shaft

0.220.0174

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



19ES

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0175

S683-GR7-Y

High Precision Bearing



Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



200

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	1.000.7185	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Rotor High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	1.000.7185	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon


200XD

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520X	

201

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	1.000.7185	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Rotor High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	1.000.7185	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon

2068D

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	0.540.2491 + 0.540.1151	BKU69151	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

20E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.2004

S683B-ZZS-MC3-GR3

High Precision Bearing

Double shielded Abec 3 precision stainless steel bearing (3x7x2.5)



Bearing for Shaft

0.220.0098

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



20LH

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spindle

0.534.0681

BKU60681



Bearing for Shaft

0.220.0098

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



20LP

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spindle

0.534.0681

BKU60681



Bearing for Shaft

0.220.0098

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



2304LD - 553 5121

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

220 0735

SMR74-ZS-T4HB-MC3-GR5-AF12

High Speed Bearing - Torlon Retainer

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)





2 Required




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



2307

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Chuck End Stainless Bearing Deleted Item. Available while stocks last. (2x7x3)		S602B-ZZS-MC3-JP6	
Replacement Drive Ball High Precision Stainless Ball Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball 6 Required		** BL-00210-S440-G3	
Bearing for Centre Stainless Bearing Deleted Item. Available while stocks last. (4x11x4)	220-0158	S694-MC3-JP6-AF12	
Bearing for Cap End High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)	220-0157 / 0779	S692-ZZS-MC34-GR5-R-G48	 Image Coming Soon

2307L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)	0.220.0144	S692-ZS-GR3-KC	 Image Coming Soon

2307LN

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	0.220.0161	S683-MC3-JP6-AF12	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	0.220.0174	SMR104A-ZZS-MC3-GR7-T9H	 Image Coming Soon

2313

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x10x4)	0.220.0139	SMR104A-ZZ-MC3-JP6	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

2313L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0144

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)



2313LN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0161

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



2320

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Drive Ball

**** BL-00210-S440-G3**

High Precision Stainless Ball

Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball

6 Required



Bearing for Chuck End

Stainless Bearing

Deleted Item. Available while stocks last. (2x7x3)

S602B-ZZS-MC3-JP6



Bearing for Centre

220-0158

S694-MC3-JP6-AF12

Stainless Bearing

Deleted Item. Available while stocks last. (4x11x4)



Bearing for Cap End

220-0157 / 0779

S692-ZZS-MC34-GR5-R-G48

High Precision Bearing

Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)



2320L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0144

S692-ZS-GR3-KC

High Speed Bearing - Phenolic Retainer





Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


2320LN

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Chuck End Stainless Bearing Deleted Item. Available while stocks last. (2x7x3)		S602B-ZZS-MC3-JP6	
Replacement Drive Ball High Precision Stainless Ball Deleted Item. Available while stocks last. Stainless steel high precision 2.1mm ball 6 Required		** BL-00210-S440-G3	
Bearing for Centre Stainless Bearing Deleted Item. Available while stocks last. (4x11x4)	220-0158	S694-MC3-JP6-AF12	
Bearing for Cap End High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)	220-0157 / 0779	S692-ZZS-MC34-GR5-R-G48	 Image Coming Soon


2325L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	0.220.0161	S683-MC3-JP6-AF12	
Bearing for Shaft High Speed Bearing - Phenolic Retainer Single shielded Abec 3 precision stainless steel bearing with a high speed phenolic retainer. Non ISO Width (2x6x2.3)	0.220.0144	S692-ZS-GR3-KC	 Image Coming Soon

2325LN

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)	0.220.0175	S683-GR7-Y	 Image Coming Soon

2329

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open Abec 7 precision stainless steel bearing with bronze cage (3x7x2)	0.220.0175	S683-GR7-Y	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

2329LN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0174

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



23ES

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0098

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



24LH

Position / Description / Application

Manufacturers Part

MBA Part

Front Gear 1:5 Middle Gear 0.553.7201

0.553.8282

BKU68282



24LN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.0734

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



25LCA

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Standard

0.553.7341

BKU67341X

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Head Gear Complete

0.553.7341

BKU67341

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

25LH

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	0.553.7031 + 0.553.9932	BKU67031	
Drive Kit	553.7341 + 0.553.7201	BKU64101	
Drive Kit	553.7341 + 0.553.7201 + 0.553.7251	BKU64151	
Front Gear 1:5 Middle Gear 0.553.7201	0.553.8282	BKU68282	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)	0.220.0730	SMR5-0235-T4HB-MC3-AF12	
Intermediate Gear 1:5	0.553.7251	BKU67251	
Replacement Middle Gear Complete 1:5	0.553.7201	BKU67201	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	1.001.2836	SMR74-CSN-TH7-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

25LHA

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Middle Gear Complete 1:5	0.553.7201	BKU67201	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)	0.220.0730	SMR5-0235-T4HB-MC3-AF12	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	0.553.7341	BKU67341	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)	0.220.0171	S684-MC3-JP6-AF12	
Drive Kit	553.7341 + 0.553.7201 + 0.553.7251	BKU64151	
Head Gear Complete - Standard For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	0.553.7341	BKU67341X	
Replacement Press Button	0.553.7031 + 0.553.9932	BKU67031	
Intermediate Gear 1:5	0.553.7251	BKU67251	
Drive Kit	553.7341 + 0.553.7201	BKU64101	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Bearing for Headgear	0.220.0734	SMR74-CSN-TH7-TO
High Speed Ceramic-Hybrid Bearing - Torlon Retainer		
Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)		



25LHC

Position / Description / Application	Manufacturers Part	MBA Part
Bearing for Shaft	0.220.0171	S684-MC3-JP6-AF12
High Precision Bearing		
Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)		




25LHP

Position / Description / Application	Manufacturers Part	MBA Part
Front Gear 1:5 Middle Gear 1.001.3515	1.001.3516	BKU63516



Replacement Middle Gear Complete 1:5	1.001.3515	BKU63515
---	-------------------	-----------------



Intermediate Gear 1:5	1.001.3519	BKU63519
------------------------------	-------------------	-----------------



Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

25LP

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Dog Gear 1:5 Gear	1.001.3071	BKU63071	
Replacement Spring Washer	1.001.2987	BKU62987	
Replacement Middle Gear Complete 1:5	1.001.3515	BKU63515	
Front Gear 1:5 Middle Gear 1.001.3515	1.001.3516	BKU63516	
Intermediate Gear 1:5	1.001.3519	BKU63519	
Replacement Press Button	1.001.3284 + 0.201.8157 + 1.001.3283	BKU63284	
Drive Kit	1.001.3520 + 1.001.3515 + 1.001.3519	BKU62019	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520M	
Drive Kit	1.001.3520 + 1.001.3515	BKU62015	

Key Usually stocked or Short Delivery
** Available while stocks last










Long Delivery or Large Minimum Quantities
Not in current production

Replacement Bearing Special Integrated Bearing Special Size - to check	1.001.2833	BKU22833	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open Abec 7 precision bearing with Torlon retainer and light oil lubricant (2x5x1.5)	1.001.2834	S682-MC34-GR7-TO-L23	 Image Coming Soon
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	1.001.2836	SMR74-CSN-TH7-TO	 Image Coming Soon
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x1,905/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x1.9)	1.001.3167	SWR144C-D38-CSN-TO	 Image Coming Soon
Bearing for Idler Gear High Speed Bearing - Torlon Retainer Duplex tube bearing Abec 7 precision with high speed Torlon retainer (2x4.7x6)	1.001.2835	SMR4.762-200-ZZS-MC3-GR7-TO	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

25LPA


Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.001.2987	BKU62987	
Drive Kit	1.001.3520 + 1.001.3515	BKU62015	
Drive Kit	1.001.3520 + 1.001.3515 + 1.001.3519	BKU62019	
Intermediate Gear 1:5	1.001.3519	BKU63519	
Head Gear Complete Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520X	
Head Gear Complete Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520M	
Replacement Bearing Special Integrated Bearing Special Size - to check	1.001.2833	BKU22833	
Head Gear Complete Assembled with ceramic ball bearings, tested in master electric handpiece	1.001.3520	BKU63520	
Replacement Middle Gear Complete 1:5	1.001.3515	BKU63515	

Key Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
Not in current production

Replacement Press Button	1.001.3284 + 0.201.8157 + 1.001.3283	BKU63284	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open Abec 7 precision bearing with Torlon retainer and light oil lubricant (2x5x1.5)	1.001.2834	S682-MC34-GR7-TO-L23	
Bearing for Idler Gear High Speed Bearing - Torlon Retainer Duplex tube bearing Abec 7 precision with high speed Torlon retainer (2x4.7x6)	1.001.2835	SMR4.762-200-ZZS-MC3-GR7-TO	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x1,905/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x1.9)	1.001.3167	SWR144C-D38-CSN-TO	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	1.001.2836	SMR74-CSN-TH7-TO	


282A Slow Speed

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (5x16x5)	247-6053 / 0063	S625-ZZ-MC3-JP6	

29LH

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	0.220.0138	SMR104A-ZZS-MC3-GR7-T9H	



29LP

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	0.220.0138	SMR104A-ZZS-MC3-GR7-T9H	



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

31E

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	0.220.0104	SMR63-MC3-JP6-AF12	
Bearing for Headgear High Precision Bearing Open Abec 3 precision bearing (2x5x1.5)	1.003.5383	S682-MC3-GR3-GPR	 Image Coming Soon

31ES

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)	0.220.0098	S693-ZZS-MC34-TH5-R-L23	 Image Coming Soon
Bearing for Headgear High Precision Bearing Open Abec 3 precision bearing (2x5x1.5)	1.003.5383	S682-MC3-GR3-GPR	 Image Coming Soon

3548

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	0.220.0103	S684A-ZZ-MC3-JP6	

3549

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	0.220.0103	S684A-ZZ-MC3-JP6	

3550

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	0.220.0103	S684A-ZZ-MC3-JP6	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

3554

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0103

S684A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



3554 Surgical Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Nosecone Bearing

SR133-MC34-GR5-R-L23

Extra High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (2.38x4.76x1.59)



3555

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0103

S684A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x9x4)



3610

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

0.220.1302

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option

Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)

2 Required



Replacement Bearing

0.220.1302

SMR85C-ZZS-MC34-GR5-L23

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)

2 Required



3610N1 Staright Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

0-220-1306

SMR9525-050-MC3-L23

Stainless Bearing

Double shielded standard bearing (5x9.52x3.18)

2 Required



3640L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0020

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

38A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0106

SMR85-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



4005

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 7-14-5mm

0.220.0112

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)



4500

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

1.001.6391

BKU66391



Bearing for Rotor Cartridge

1.007.6975

SWMR6.45-3.175-CSN-ZZS-MC3-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

3.175x6.45x2.381/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.45x2.38)



Bearing for Rotor Cartridge

SG1426X

SWMR6.45-3.175-ZZS-MC3-GR7-KC

High Speed Bearing - Phenolic Retainer Shielded Option

Double shielded Abec 7 precision special dental bearing with high speed Phenolic retainer (3.18x6.45x2.38)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

4500B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.006.4564	BKURC2100105	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.006.4564	BKUR2100104	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.006.4564	BKURC210010	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.006.4564	BKUR2100103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.006.4564	BKURC2100103	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.006.4564	BKUR210010	
Replacement Press Button	1.006.1079	BKU61079	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.006.4564	BKUR2100105	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

465LRN

Position / Description / Application

Manufacturers Part


MBA Part

Multiflex Coupler O-ring

0.553.5262

BKU65262

Special dental applications O-Ring



Multiflex Coupler O-ring

0.200.6037+0.200.6113

BKU66113

Special dental applications O-Ring



465LRN Coupling

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Valve

0.553.9452

0.553.9452

high precision valve for dental applications



Replacement O-Ring

0.200.6038

**** 0.200.6038**

Deleted Item. Available while stocks last. special o-ring for dental applications




Replacement Ball

0.257.4013

**** 0.257.4013**

Deleted Item. Available while stocks last. High precision ball for dental applications



465R

Position / Description / Application

Manufacturers Part


MBA Part

Multiflex Coupler O-ring

0.200.6037+0.200.6113

BKU66113

Special dental applications O-Ring



465RN

Position / Description / Application

Manufacturers Part

MBA Part

Multiflex Coupler O-ring

0.553.1872

BKU61872

Special dental applications O-Ring



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

5000

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

1.001.6391

BKU66391



5000B

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

1.006.4564

BKUR2100103

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

1.006.4564

BKUR2100104

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



Replacement Press Button

1.006.1079

BKU61079




Replacement Turbine Excluding Bearings

1.006.4564

BKURC210010

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

1.006.4564

BKURC2100103

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

1.006.4564

BKUR2100105

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

1.006.4564

BKURC2100105

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine Excluding Bearings

1.006.4564

BKUR210010

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

6000

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-KC	 Image Coming Soon
Replacement Spring Washer	1.001.5032	BKU65032	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-Z-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	SG1304X	SWR144-D28W-Z-KC	 Image Coming Soon

Key Usually stocked or Short Delivery ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

6000B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.001.3888	BKUR2100114	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.001.3888	BKUR2100113	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.001.3888	BKUR2100115	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.001.3888	BKURC210011	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.001.3888	BKURC2100115	
Replacement Press Button	1.000.4770	BKU64770	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.001.3888	BKUR210011	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.001.3888	BKURC2100113	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.001.3888	BKURMU210011	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Turbine Sub Assembly Without Bearings
For Handpiece with Wedge Type Chuck

1.001.3888

BKURM210011



614A

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

220-0125

SR144A-D105-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer
Stepped OD bearing with high speed phenolic retainer for precision dental handpieces
2 Required



Replacement O-Ring

OD of Bearing

OR-00600-100-N70

Special Dental Application O-Ring
Special sized o-ring used in dental Applications
2 Required



Replacement Bearing

220-0125

SR144A-D105-CSN-ZS-KC-OPT

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer
Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces
2 Required



Replacement Bearing

220-0125

SR144A-D105-CSN-ZS-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer
Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces
2 Required



Replacement Bearing

220-0793K

SR144A-D107-CSN-ZS-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer
Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces
2 Required



Replacement Bearing

220-0125

SR144A-D105-ZS-TO-OPT

High Speed Bearing - Torlon Retainer
Stepped OD bearing with Torlon retainer for precision dental handpieces
2 Required



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


624

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production




625

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

625B SuperTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

625C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

625D

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

627

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

628 - 629 FG

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required	220-0128 / 0145	SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-TO	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-KC	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

628 BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required	220-0128 / 0145	SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-TO	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-KC	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

630

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last








Long Delivery or Large Minimum Quantities
Not in current production

630B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.5881	BKURC2100015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.5881	BKUR2100015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKURC2100013	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.5881	BKURC210001	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.5881	BKUR210001	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKUR2100014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKUR210001T	

Key Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
Not in current production

<p>Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required</p>	<p>OD of Bearing</p>	<p>OR-00600-100-N70</p>	
<p>Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p>	<p>0.553.5881</p>	<p>BKUR2100013</p>	
<p>Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-ZS-TO-OPT</p>	
<p>Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck</p>	<p>0.553.5881</p>	<p>BKURM210001</p>	
<p>Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic</p>	<p>0.553.5881</p>	<p>BKURMU210001</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0793K</p>	<p>SR144A-D107-CSN-ZS-TO</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-CSN-ZS-TO</p>	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

630C PB SuperTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

632D

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	

632D - 632E

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

632E

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	

634

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

635

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

0.220.2007

SR144A-D107-CSN-ZS-TO

High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer

Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces



635 Back

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1451X

BKUSG1451XSP

High Speed Bearing - Phenolic Retainer

3.175x6.35x2.381/2.778, stainless steel-THB Phenolic (3.18x6.35x2.38)



Bearing for Rotor Cartridge

0.220.2015

SWR144-D106W-ZS-TO

High Speed Bearing - Torlon Retainer

3.175x6.35x2.381/2.778, stainless steel, Torlon retainer



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production





635B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKURC2100035	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKURC2100033	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKURC210003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8581	BKUR2100035	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8581	BKUR2100034	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8581	BKUR210003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8581	BKUR2100033	
Replacement Press Button	0.553.8491	BKU68491	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

635BR

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8661	BKUR2100054	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8661	BKUR210005	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8661	BKUR2100053	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8661	BKUR2100055	





636

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

636B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKURC2100043	
Replacement Press Button	0.553.8061	BKU68061	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8911	BKURC210004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8911	BKURC2100045	

636CP

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8911	BKUR210004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR2100044	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR2100043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR210004T	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8911	BKURE2100045	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

637

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

0.220.2007

SR144A-D107-CSN-ZS-TO

High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer

Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces



637 Back

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1451X

BKUSG1451XSP

High Speed Bearing - Phenolic Retainer

3.175x6.35x2.381/2.778, stainless steel-THB Phenolic (3.18x6.35x2.38)



Bearing for Rotor Cartridge

0.220.2015

SWR144-D106W-ZS-TO

High Speed Bearing - Torlon Retainer

3.175x6.35x2.381/2.778, stainless steel, Torlon retainer



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



637B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8661	BKUR2100054	
Replacement Press Button	0.553.8491	BKU68491	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8661	BKUR210005	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKURC210005	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKURC2100055	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8661	BKUR2100053	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKURC2100053	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8661	BKUR2100055	



Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

637B Mini Lux 3 BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		OR-00600-100-N70	
Bearing for Cap End High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.		SR144A-D650-CSN-ZS-TO	

639A Friction Grip

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key	Usually stocked or Short Delivery
	** Available while stocks last







Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

639B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.6911	BKUR210006	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.6911	BKUR2100065	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKUR2100063	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKUR2100064	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-CSN-ZS-TO</p>	
639C Push Button			
<p>Position / Description / Application</p>	<p>Manufacturers Part</p>	<p>MBA Part</p>	
<p>Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-ZS-KC-OPT</p>	
<p>Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required</p>	<p>OD of Bearing</p>	<p>OR-00600-100-N70</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-CSN-ZS-KC-OPT</p>	
<p>Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-ZS-TO-OPT</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0125</p>	<p>SR144A-D105-CSN-ZS-TO</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required</p>	<p>220-0793K</p>	<p>SR144A-D107-CSN-ZS-TO</p>	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

640

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	

640A Friction Grip

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

640B



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKUR210001T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKUR2100013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.5881	BKURC2100015	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Spray Cover	0.553.8752	BKU68752	Spray Covers 
Replacement Threaded Ring	0.553.3322	BKU63322	Spray Covers 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKURC2100013	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKUR2100014	
Replacement Press Button	0.553.6511	BKU66511	Press Button cpl. for Heads & Air Driven Handpieces 

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.5881	BKUR210015	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.5881	BKURC210001	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.5881	BKUR210001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.5881	BKURMU210001	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	0.553.5881	BKURM210001	

640C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spray Cover For Handpiece with Wedge Type Chuck	0.553.8752	BKU68752	
Replacement Threaded Ring For Handpiece with Wedge Type Chuck	0.553.3322	BKU63322	

640C PB SuperTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production







640FG

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

642A FG

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

642B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.6911	BKUR210006	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKUR2100064	
Replacement Press Button	0.553.7391	BKU67391	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.6911	BKUR2100065	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.6911	BKURC2100065	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.6911	BKURC210006	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKUR2100063	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKURC2100063	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

642B - 642C PB BellaTorque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-ZS-TO-OPT	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

645

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	

645B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.6911	BKURC2100065	 <p>Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators</p>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.6911	BKURC210006	 <p>Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators</p>
Replacement Press Button	0.553.7391	BKU67391	 <p>Press Button cpl. for Heads & Air Driven Handpieces</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.6911	BKURC2100063	 <p>Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators</p>

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

646

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

646B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	0.553.8061	BKU68061	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR210004T	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8911	BKUR210004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8911	BKURC2100045	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKURC2100043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR2100043	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8911	BKURC210004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8911	BKUR2100044	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8911	BKURE2100045	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

647

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-ZS-TO-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.0125	SR144A-D105-CSN-ZS-TO	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

647B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	0.553.8061	BKU68061	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKUR2100074	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8541	BKURC2100075	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8541	BKUR210007	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKURC2100073	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8541	BKURC210007	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKUR2100073	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8541	BKUR2100075	

649

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Angular Contact Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces	0.220.2007	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

649B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8541	BKURC2100075	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKURC2100073	
Replacement Press Button	0.553.8061	BKU68061	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8541	BKUR210007	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKUR2100073	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8541	BKURC210007	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8541	BKUR2100074	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8541	BKUR2100075	

650

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.	0.220.2009	SR144A-D650-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

6500

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-KC	 Image Coming Soon
Replacement Spring Washer	1.003.7576	BKU67576	
Replacement Spring Washer	1.003.5900	BKU65900	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-Z-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	SG1304X	SWR144-D28W-Z-KC	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

6500B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.9664	BKURC2100155	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.9664	BKUR210012	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.9664	BKUR2100123	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.003.9664	BKUR2100124	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.9664	BKUR2100125	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.9664	BKURC2100153	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.9664	BKURC210015	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.003.9664	BKURM210012	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.9664	BKURMU210012	

Key Usually stocked or Short Delivery
** Available while stocks last




Long Delivery or Large Minimum Quantities
Not in current production

6500BR



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	1.003.5525	BKU65525	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.6884	BKURC210014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKURC2100143	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.6884	BKUR2100145	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.6884	BKURC2100145	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR2100143	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR210014T	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.6884	BKUR210014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR2100144	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKURMU210014	 <p style="font-size: 8px; color: red; margin-top: 5px;">Ultra precision sub-assemblies, cartridges</p>
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)	1.000.7980	SR144A-D31-CSN-WZ-TO	 Image Coming Soon
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.003.6884	BKURM210014	 <p style="font-size: 8px; color: red; margin-top: 5px;">Ultra precision sub-assemblies, cartridges</p>

650A FG

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		OR-00600-100-N70	
Bearing for Chuck End High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.	220-2009	SR144A-D650-CSN-ZS-TO	

Usually stocked or Short Delivery
** Available while stocks last



Long Delivery or Large Minimum Quantities
Not in current production

650B



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Threaded Ring	0.553.0543	BKU60543	 <small>Spray Covers</small>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8051	BKUR2100083	 <small>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</small>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	0.553.8051	BKUR2100085	 <small>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</small>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8051	BKUR210008	 <small>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</small>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8051	BKUR2100084	 <small>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</small>
Replacement Spray Cover	0.553.0573	BKU60573	 <small>Spray Covers</small>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	0.553.8051	BKURC210008	 <small>Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators</small>
Replacement Press Button	0.553.8061	BKU68061	 <small>Press Button cpl. for Heads & Air Driven Handpieces</small>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8051	BKURC2100083	 <small>Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators</small>

Key Usually stocked or Short Delivery
** Available while stocks last



Long Delivery or Large Minimum Quantities
Not in current production

Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	0.553.8051	BKURM210008	 <p style="font-size: 8px; color: red; margin-top: 5px;">Ultra precision sub-assemblies, cartridges</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	0.553.8051	BKURMU210008	 <p style="font-size: 8px; color: red; margin-top: 5px;">Ultra precision sub-assemblies, cartridges</p>



650B PB - 650C PB

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	220-2009	OR-00600-100-N70	
Bearing for Chuck End High Speed Ceramic-Hybrid Bearing - Torlon Retainer Deleted Item. Available while stocks last.	220-2009	SR144A-D650-CSN-ZS-TO	

650C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spray Cover	0.553.0573	BKU60573	 <p style="font-size: 8px; color: red; margin-top: 5px;">Spray Covers</p>
Replacement Threaded Ring	0.553.0543	BKU60543	 <p style="font-size: 8px; color: red; margin-top: 5px;">Spray Covers</p>

655

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.001.5032	BKU65032	 <p style="font-size: 8px; color: red; margin-top: 5px;">Spring Washers for Air Driven & Electric Driven Handpieces</p>
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)	1.000.7980	SR144A-D31-CSN-WZ-TO	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last





Long Delivery or Large Minimum Quantities
Not in current production

655B



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.2078	BKUR2100095	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR2100094	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.2078	BKURC210009	
Replacement Spray Cover	0.553.0573	BKU60573	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR210009T	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.2078	BKUR210009	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR2100093	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKURC2100093	
Replacement Threaded Ring	1.000.3391	BKU63391	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.2078	BKURC210095	
Replacement Press Button	1.000.2076	BKU62076	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKURMU210009	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.000.2078	BKURM210009	

655C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spray Cover	0.553.0573	BKU60573	
Replacement Threaded Ring	1.000.3391	BKU63391	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities <i>Not in current production</i>

659BR

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR2100093	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR210009T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
Replacement Spray Cover	0.553.0573	BKU60573	Spray Covers 
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.2078	BKUR210009	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKUR2100094	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.2078	BKUR2100095	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKURMU210009	Ultra precision sub-assemblies, cartridges 
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.000.2078	BKURM210009	Ultra precision sub-assemblies, cartridges 





660

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.001.5032	BKU65032	Spring Washers for Air Driven & Electric Driven Handpieces 


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


660B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.2078	BKURC2100095	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.2078	BKURC210009	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.2078	BKURC2100093	
Replacement Press Button	1.000.2076	BKU62076	


66GD

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	0.220.0105	SMR74-MC34-GR7-R-L23	

66LU

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	0.220.2005	SMR74-MC3-JP6-AF12	

67G

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	0.220.0104	SMR63-MC3-JP6-AF12	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

67GD

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.0105

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



67L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.0051

SWA683B-ZZ-C69-GR5-TO

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



67LDN

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Headgear

0.220.0105

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



Bearing for Headgear

0.220.0713

S683A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (3x7x3)





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

67LH

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	0.540.3671 + 0.540.1161	BKU67151	
Intermediate Piece Dental Head 2:1	0.540.1161	BKU61161	
Replacement Head Housing	0.540.3641	BKU63641	
Replacement Head Gear 2:1 For Handpiece with Latch Type Chuck	0.540.4022	BKU64022	
Head Gear Complete 1:1 For Handpiece with Latch Type Chuck	0.540.3671	BKU63671	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	0.220.2005	SMR74-MC3-JP6-AF12	

67RIC

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0023	SMR85-MC3-JP6-AF12	
Bearing for Headgear High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	0.220.0713	S683A-ZZ-MC3-JP6	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

68B - 68C - 68CN - 68D

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)

2 Required

SMR63-MC3-JP6-AF12



68CN

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

For Handpiece with Latch Type Chuck

0.540.1781

BKU61781



Drive Kit

0.540.1781 + 0.540.1151

BKU68151



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

68G

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Screw	0.201.0005	BKU60005	
Fitting set Includes 0.540.4072, 0.540.4102, 0.540.4522	0.540.4072, 0.540.4102, 0.540.4522	BKU64102	
Head Gear Complete 1:1 For Handpiece with Latch Type Chuck	0.540.1781	BKU61781	
Bush Intermediate Piece 0.540.1151	0.540.2392	BKU62392	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	0.220.0104	SMR63-MC3-JP6-AF12	
Intermediate Piece Dental Head 1:1	0.540.1151	BKU61151	
Replacement Spring	0.540.3162	BKU63162	
Drive Kit	0.540.1781 + 0.540.1151	BKU68151	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

68GDN

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

0.540.2491

BKU62491

For Handpiece with Latch Type Chuck



Replacement Head Gear 1:1

0.540.4092

BKU64092

For Handpiece with Latch Type Chuck



68L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.0051

SWA683B-ZZ-C69-GR5-TO

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

68LDN

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	0.540.2491 + 0.540.1151	BKU69151	
Replacement Press Button	0.540.2501	BKU62501	
Bearing for Chuck End Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	220-0105	SMR74-MC3-JP6-AF12	
Bearing for Headgear High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	0.220.0105	SMR74-MC34-GR7-R-L23	
Replacement Screws	0.201.0005	BKU60005	
Intermediate Piece Dental Head 1:1	0.540.1151	BKU61151	
Replacement Head Housing	0.540.2231	BKU62231	
Bush Intermediate Piece 0.540.1151	0.540.2392	BKU62392	
Head Gear Complete 1:1 For Handpiece with Latch Type Chuck	0.540.2491	BKU62491	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Replacement Spring 0.540.3162	BKU63162	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Bearing for Headgear High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	0.220.0713	S683A-ZZ-MC3-JP6 
Bearing for Cap End Stainless Bearing Double shielded P6 precision stainless steel bearing (3x7x3)	220-0713 / 0102	S683A-ZZ-MC3-JP6 
Replacement Head Gear 1:1 For Handpiece with Latch Type Chuck	0.540.4092	BKU64092  <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Fitting set Includes 0.540.4072, 0.540.4102, 0.540.4522	0.540.4072, 0.540.4102, 0.540.4522	BKU64102  <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>

Key **Usually stocked or Short Delivery**
 ** Available while stocks last





Long Delivery or Large Minimum Quantities
Not in current production

68LH

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1 For Handpiece with Latch Type Chuck	0.540.3661	BKU63661	
Replacement Spring	0.540.3162	BKU63162	
Bush Intermediate Piece 0.540.1151	0.540.2392	BKU62392	
Intermediate Piece Dental Head 1:1	0.540.1151	BKU61151	
Replacement Head Housing	0.540.3631	BKU63631	
Drive Kit	0.540.3661 + 0.540.1151	BKU66151	
Replacement Press Button	0.540.3741	BKU63741	
Replacement Head Gear 1:1 For Handpiece with Latch Type Chuck	0.540.4012	BKU64012	
Replacement Screws	0.201.0005	BKU60005	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

<p>Replacement Rear Bearing Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)</p>	<p>S683-MC3-JP6-AF12</p>	
<p>Fitting set Includes 0.540.4072, 0.540.4102, 0.540.4522</p>	<p>0.540.4072, 0.540.4102, 0.540.4522</p> <p>BKU64102</p>	 <p>Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces & Heads</p>
<p>Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)</p>	<p>0.220.2005</p> <p>SMR74-MC3-JP6-AF12</p>	
<p>Front Bearing High Speed Bearing - Torlon Retainer Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)</p>	<p>SWA683B-ZZ-C69-GR5-TO</p>	 <p>MBA Image Coming Soon</p>


Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

68LU

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1 For Handpiece with Latch Type Chuck	0.540.3661	BKU63661	
Bush Intermediate Piece 0.540.1151	0.540.2392	BKU62392	
Fitting set Includes 0.540.4072, 0.540.4102, 0.540.4522	0.540.4072, 0.540.4102, 0.540.4522	BKU64102	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	0.220.2005	SMR74-MC3-JP6-AF12	
Intermediate Piece Dental Head 1:1	0.540.1151	BKU61151	
Replacement Spring	0.540.3162	BKU63162	
Replacement Head Gear 1:1 For Handpiece with Latch Type Chuck	0.540.4012	BKU64012	
Replacement Screws	0.201.0005	BKU60005	

6N

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	0.220.0103	S684A-ZZ-MC3-JP6	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

7000

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.001.5032	BKU65032	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-KC	
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-Z-TO	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	SG1304X	SWR144-D28W-Z-KC	

Key Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
Not in current production

7000B

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.000.4772	BKUR2100134	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.4772	BKUR2100135	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.4772	BKUR210013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.4772	BKUR2100133	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.4772	BKURC2100133	
Replacement Press Button	1.000.4770	BKU64770	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.000.4772	BKURC210013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.000.4772	BKURC2100135	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.000.4772	BKURMU210013	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.000.4772	BKURM210013	
---	-------------------	--------------------	---



700KL


Position / Description / Application Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)	Manufacturers Part 0.220.0219	MBA Part S685A-ZZS-MC34-TH7-R-G48	
---	---	---	---



Position / Description / Application Bearing for Shaft High Precision Bearing Shielded Option Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)	Manufacturers Part 0.220.0218	MBA Part S693-ZZS-MC34-TH5-R-L23	
--	---	--	---



79G

Position / Description / Application Bearing for Headgear High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	Manufacturers Part 0.220.0105	MBA Part SMR74-MC34-GR7-R-L23	
---	---	---	--




7LH

Position / Description / Application Replacement Gear Unit	Manufacturers Part 0.534.1501	MBA Part BKU61501	
--	---	-----------------------------	---






Position / Description / Application Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	Manufacturers Part 0.220.0139	MBA Part SMR104A-ZZ-MC3-JP6	
---	---	---------------------------------------	---




Position / Description / Application Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	Manufacturers Part 0.220.0174	MBA Part SMR104A-ZZS-MC3-GR7-T9H	
---	---	--	---




7LP

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Gear Unit	0.534.1501	BKU61501	 <p style="font-size: 8px; color: red;">Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	0.220.0139	SMR104A-ZZ-MC3-JP6	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)	0.220.0174	SMR104A-ZZS-MC3-GR7-T9H	 <p style="text-align: center;">Image Coming Soon</p>

8000

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.003.5900	BKU65900	 <p style="font-size: 8px; color: red;">Spring Washers for Air Driven & Electric Driven Handpieces</p>

8000 Front

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)	1.000.7980	SR144A-D31-CSN-WZ-TO	 <p style="text-align: center;">Image Coming Soon</p>

Key Usually stocked or Short Delivery ** Available while stocks last



Long Delivery or Large Minimum Quantities *Not in current production*

8000B



Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR2100143	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR2100144	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.6884	BKUR2100145	
Replacement Press Button	1.003.5525	BKU65525	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.6884	BKUR210014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKUR210014T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.003.6884	BKURC210014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKURC2100143	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.003.6884	BKURC2100145	

Key Usually stocked or Short Delivery
** Available while stocks last


Long Delivery or Large Minimum Quantities
Not in current production

Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	1.003.6884	BKURM210014	 <p style="font-size: 8px; color: red; margin: 0;">Ultra precision sub-assemblies, cartridges</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.003.6884	BKURMU210014	 <p style="font-size: 8px; color: red; margin: 0;">Ultra precision sub-assemblies, cartridges</p>


80L

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	1.009.1782 + 0.540.1151	BKU68251	 <p style="font-size: 8px; color: red; margin: 0;">Drive Kits for Electric Driven Handpieces and Heads</p>
Bearing for Headgear High Precision Bearing Open Abec 7 precision stainless steel-bronze bearing - Checking details (4x7x2)	0.220.0735	SMR74-GR7-Y	 <p style="font-weight: bold; margin: 0;">Image Coming Soon</p>


80LD

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open Abec 7 precision stainless steel-bronze bearing - Checking details (4x7x2)	0.220.0735	SMR74-GR7-Y	 <p style="font-weight: bold; margin: 0;">Image Coming Soon</p>

81

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0106	SMR85-MC3-JP6-AF12	



C09

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	0.220.0020	S683-MC3-JP6-AF12	

Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

CL09

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0023	SMR85-MC3-JP6-AF12	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x11x3)	1.001.2581	S685-MC3-JP6-AF12	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

CL3

Position / Description / Application	Manufacturers Part	MBA Part	
Intermediate Piece Dental Head 3:1	1.003.4207	BKU64207	
Head Gear Complete 3:1 For Handpiece with Latch Type Chuck	1.001.5005	BKU65005	
Replacement Head Housing	1.001.9257	BKU69257	
Drive Kit	1.011.5595 + 1.003.4207	BKU69507	
Head Gear Complete 3:1 For Handpiece with Latch Type Chuck	1.011.5595	BKU65595	
Replacement Press Button	1.001.4486+0.540.3752+0.201.8958	BKU64486	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	0.220.0023	SMR85-MC3-JP6-AF12	
Bearing for Headgear High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	0.220.0051	SWA683B-ZZ-C69-GR5-TO	

Coupling - Low Speed

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Special sized o-ring used in dental Applications 2 Required	200-6159	OR-00830-075-N70	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

E 25C

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete

1.008.4040

BKU64040

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



E120

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.2005

SMR74-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



E15C

Position / Description / Application


Manufacturers Part

MBA Part

Drive Kit

0.540.3671 + 0.540.1161

BKU67151



E15L

Position / Description / Application

Manufacturers Part

MBA Part

Intermediate Piece Dental Head 2:1

0.540.1161

BKU61161



Replacement Head Gear 2:1

0.540.4022

BKU64022

For Handpiece with Latch Type Chuck



Head Gear Complete 1:1

0.540.3671

BKU63671

For Handpiece with Latch Type Chuck



Replacement Gear Unit

0.534.1501

BKU61501



Bearing for Headgear

0.220.0051

SWA683B-ZZ-C69-GR5-TO

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

E20C

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

0.540.3661 + 0.540.1151

BKU66151



E20L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

0.220.0051

SWA683B-ZZ-C69-GR5-TO

High Speed Bearing - Torlon Retainer

Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)



E25

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

1.008.4040 + 1.007.9916

BKU64016



Drive Kit

1.014.5555 + 1.014.5557

BKU65557



Front Gear 1:5 Middle Gear

1.007.9916

BKU69999



E25C

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Standard

1.008.4040

BKU64040X

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Replacement Press Button

1.008.1488

BKU61488



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

E25L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Dog Gear 1:5 Gear	1.007.5555	BKU65555	
Head Gear Complete - Standard For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.008.4040	BKU64040X	
Head Gear Complete - Premium For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.008.4040	BKU64040M	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.008.4040	BKU64040	
Replacement Spring Washer	1.001.2987	BKU62987	
Replacement Press Button	1.008.1488	BKU61488	
Replacement Middle Gear Complete 1:5	1.007.9916	BKU69916	
Bearing for Headgear Japanese High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 3 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	1.008.4039	SMR74-CSN-J3-TO	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x1,905/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x1.9)	1.008.4038	SWR144C-D38-CSN-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Bearing for Headgear	1.008.4039	SMR74-CSN-TH7-TO
High Speed Ceramic-Hybrid Bearing - Torlon Retainer		
Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)		



E31

Position / Description / Application	Manufacturers Part	MBA Part
Bearing for Headgear	0.220.0104	SMR63-MC3-JP6-AF12
High Precision Bearing		
Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)		



E40

Position / Description / Application	Manufacturers Part	MBA Part
Bearing for Shaft	0.220.0020	S683-MC3-JP6-AF12
High Precision Bearing		
Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)		



Bearing for Headgear	0.220.2005	SMR74-MC3-JP6-AF12
High Precision Bearing		
Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)		



E675

Position / Description / Application	Manufacturers Part	MBA Part
Bearing for Rotor Cartridge	1.000.7979	SWR144-D28W-CSN-Z-KC
High Speed Ceramic-Hybrid Bearing - Phenolic Retainer		
Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)		



Bearing for Rotor Cartridge	SG1304X	SWR144-D28W-Z-KC
High Speed Bearing - Phenolic Retainer		
Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)		



Bearing for Rotor Cartridge	1.000.7979	SWR144-D28W-Z-TO
High Speed Bearing - Torlon Retainer		
Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)		





Bearing for Rotor Cartridge	1.000.7979	SWR144-D28W-CSN-Z-TO
High Speed Ceramic-Hybrid Bearing - Torlon Retainer		
Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)		



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

E675 Front

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Shielded Option Double shielded Abec 7 precision special dental bearing with high speed Phenolic retainer (3.18x6.45x2.38)	SG1426X	SWMR6.45-3.175-ZZS-MC3-GR7-KC	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option 3.175x6.45x2.381/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.45x2.38)	1.007.6975	SWMR6.45-3.175-CSN-ZZS-MC3-TO	 Image Coming Soon

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production





E675L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9457	BKUR2100175	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9457	BKUR210017	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKUR2100173	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKUR2100174	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKURC2100173	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Press Button	1.007.7051	BKU67051	Press Button cpl. for Heads & Air Driven Handpieces 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9457	BKURC2100175	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9457	BKURC210017	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

E677

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-KC	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-CSN-Z-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)	SG1304X	SWR144-D28W-Z-KC	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Bearing - Torlon Retainer Single shielded dental special stainless steel bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)	1.000.7979	SWR144-D28W-Z-TO	 Image Coming Soon

E677 Front

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option 3.175x6.45x2.381/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.45x2.38)	1.007.6975	SWMR6.45-3.175-CSN-ZZS-MC3-TO	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Shielded Option Double shielded Abec 7 precision special dental bearing with high speed Phenolic retainer (3.18x6.45x2.38)	SG1426X	SWMR6.45-3.175-ZZS-MC3-GR7-KC	 Image Coming Soon

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



E677L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9457	BKURC210017	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKUR2100173	
Replacement Press Button	1.007.7051	BKU67051	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9457	BKUR2100175	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9457	BKUR210017	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKURC2100173	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9457	BKURC2100175	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9457	BKUR2100174	



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

E679

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	1.003.7576	BKU67576	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)	1.000.7980	SR144A-D31-CSN-WZ-TO	

E679 Back

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Shielded Option Double shielded Abec 7 precision special dental bearing with high speed Phenolic retainer (3.18x6.45x2.38)	SG1426X	SWMR6.45-3.175-ZZS-MC3-GR7-KC	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option 3.175x6.45x2.381/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.45x2.38)	1.007.6975	SWMR6.45-3.175-CSN-ZZS-MC3-TO	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

E679L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR2100154	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKURC210018	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKURC2100183	
Replacement Press Button	1.007.8903	BKU68903	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR210015T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKUR2100155	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKUR210015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR2100153	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKURC2100185	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

E680

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

1.000.7980

SR144A-D31-CSN-WZ-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision dental handpiece ceramic-hybrid bearing with Torlon retainer (3.18x6.35x2.78)



E680 Back

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1426X

SWMR6.45-3.175-ZZS-MC3-GR7-KC

High Speed Bearing - Phenolic Retainer Shielded Option

Double shielded Abec 7 precision special dental bearing with high speed Phenolic retainer (3.18x6.45x2.38)



Bearing for Rotor Cartridge

1.007.6975

SWMR6.45-3.175-CSN-ZZS-MC3-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

3.175x6.45x2.381/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.45x2.38)



Key

Usually stocked or Short Delivery

**** Available while stocks last**

Long Delivery or Large Minimum Quantities

Not in current production







E680L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKURC2100183	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKURC2100185	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKURC210018	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR210015T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR2100154	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.007.9313	BKUR2100153	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.007.9313	BKUR210015	
Replacement Press Button	1.007.8903	BKU68903	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.007.9313	BKUR2100155	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

ERGOgrip

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing	1.004.0203	BKU20203	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Set, Set - to check			
Bearing for Shaft 7-14-5mm	1.004.0207	S687A-ZZS-MC3-GR3-T9H	 <p style="text-align: center;">Image Coming Soon</p>
High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)			
Intra 24CN - 24LN - 553 5761			
Replacement Bearing	220 0735	SMR74-ZS-T4HB-MC3-GR5-AF12	
High Speed Bearing - Torlon Retainer Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2) 2 Required			
Intra 24LS - 553 5641			
Replacement Bearing	2200185	SMR55-0235-MC3-TO-AF12	 <p style="text-align: center;">Image Coming Soon</p>
High Speed Bearing - Torlon Retainer Open precision dental handpiece bearing with Torlon retainer (2.35x5.5x2) 2 Required			
Intra 24LZ - 24LN - 553 5641			
Replacement Bearing	2200010	SMR55-0235-MC3-GR3-AF12	 <p style="text-align: center;">Image Coming Soon</p>
High Precision Bearing Open Abec 3 precision bearing (2.35x5.5x2) 2 Required			
Intra 25CH Contra Angle			
Replacement Bearing	2201305 / 1307	SMR74-ZS-T4HB-MC3-GR5-AF12	
High Speed Bearing - Torlon Retainer Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2) 2 Required			

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Intra 25LN

Position / Description / Application

Manufacturers Part

MBA Part

Stem Bearing

2200730 in 553 7201 Stem

SMR5-0235-T4HB-MC3-AF12

High Speed Bearing - Torlon Retainer

Open standard bearing with high speed Torlon retainer and low-torque light oil lubricant (2.35x5x1.5)



Replacement Bearing

220 1307 in 553 7341 Cartridge

SMR74-ZS-T4HB-MC3-GR5-AF12

High Speed Bearing - Torlon Retainer

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

K10

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Motor Brush	0.674.2481	BKU62481	
Replacement Distance Bush Laboratory Drive	0.674.4162	BKU64162	
Replacement Chuck True running accuracy less than 0.03, bore diameter 2.35mm	0.674.7752	BKU67752	
Replacement Rod Laboratory Drive	0.674.1191	BKU61191	
Bearing for Shaft 3-8-4mm High Speed Bearing - Polyamide Retainer Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	0.220.0778	S693A-ZZS-MC3-T9H	
Replacement Quad Ring Laboratory Drive	0.224.2408	BKU62408	
Replacement Cap Laboratory Drive	0.674.2861	BKU62861	
Replacement Screw nut Laboratory Drive	0.674.2871	BKU62871	
Replacement Spring Washer Laboratory Drive	0.222.4045	BKU64045	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Replacement Lock Ring Laboratory Drive	0.222.4126	BKU64126	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Replacement Bush Laboratory Drive	0.674.6302	BKU66302	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Replacement Roll Pin Laboratory Drive	0.674.4152	BKU64152	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Replacement Screw Laboratory Drive	0.674.4172	BKU64172	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Replacement Driver Laboratory Drive	0.674.4192	BKU64192	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Replacement Spring Washer Laboratory Drive	0.222.4018	BKU64018	 <p style="font-size: 8px; color: red;">Spare Parts for Laboratory Drives</p>
Bearing for Shaft 8-16-6mm High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)	0.220.0113	S688B-ZZ-MC3-GR3-T9H-G36	 <p style="text-align: center;">Image Coming Soon</p>
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)	0.220.0779	SMR104A-ZZS-MC3-GR3-T9H	 <p style="text-align: center;">Image Coming Soon</p>

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

K11

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Set, Set - to check	0.674.3841	BKU23841	
Replacement Roll Pin Laboratory Drive	0.674.4152	BKU64152	
Replacement Distance Bush Laboratory Drive	0.674.4162	BKU64162	
Bearing for Shaft 3-8-4mm High Speed Bearing - Polyamide Retainer Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	0.220.0778	S693A-ZZS-MC3-T9H	
Replacement Screw Laboratory Drive	0.674.4172	BKU64172	
Replacement Driver Laboratory Drive	0.674.4192	BKU64192	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (4x11x4)	0.220.0779	SMR104A-ZZS-MC3-GR3-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 7 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x7.5)	0.220.0045	SR11010B-ZZS-MC3-GR7-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x6)	0.220.0044	SR11010A-ZZS-MC3-GR3-T9H	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

K12 PowerGrip®

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

0.220.0055

BKU20055

High Speed Bearing - Polyamide Retainer

7x14x17.5, stainless steel, Polyamide retainer - Checking Details



K186

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Motor Brush

0.674.0241

BKU60241



K190

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Dental Motor

0.676.0871


BKU60871



Replacement Motor Brush

0.676.0881

BKU60881




Replacement Motor O-Ring

0.200.6120

BKU66120

Special dental applications O-Ring



K192

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Dental Motor

0.676.0971

BKU60971



K196

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Dental Motor

0.676.0871

BKU60871



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

K200

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Motor O-Ring

0.200.6120

BKU66120


Special dental applications O-Ring



Replacement Dental Motor

1.000.6859

BKU66859



K201

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Motor Brush

0.676.0881

BKU60881



K4

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

0.674.4871

BKU64871

True running accuracy less than 0.03, bore diameter 2.35mm



Replacement Bearing

0.220.0058

BKU20058

Set, Set - to check



Bearing for Shaft 10x19x7

0.220.0056

S6800B-ZZS-MC3-GR3-KC

High Speed Bearing - Phenolic Retainer




Double shielded Abec 3 stainless steel bearing with a high speed phenolic cage (10x19x7)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

K5

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Set, Set - to check	0.220.0058	BKU20058	
Replacement Chuck True running accuracy less than 0.03, bore diameter 2.35mm	0.674.4871	BKU64871	
Bearing for Shaft 10x19x7 High Speed Bearing - Phenolic Retainer Double shielded Abec 3 stainless steel bearing with a high speed phenolic cage (10x19x7)	0.220.0056	S6800B-ZZS-MC3-GR3-KC	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

K9



Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Centre Stainless Bearing Double shielded P6 precision stainless steel bearing (4x11x4)	220-0157 / 0779	S694-ZZ-MC3-JP6	
Replacement Motor Brush	0.674.0241	BKU60241	
Replacement Screw Laboratory Drive	0.674.4172	BKU64172	
Replacement Spring Washer Laboratory Drive	0.222.4018	BKU64018	
Replacement Driver Laboratory Drive	0.674.4192	BKU64192	
Replacement Roll Pin Laboratory Drive	0.674.4152	BKU64152	
Replacement Chuck True running accuracy less than 0.03, bore diameter 2.35mm	0.674.7752	BKU67752	
Replacement Lock Ring Laboratory Drive	0.222.4126	BKU64126	
Replacement Bush Laboratory Drive	0.674.6302	BKU66302	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Replacement Quad Ring Laboratory Drive	0.224.2408	BKU62408	
Replacement Distance Bush Laboratory Drive	0.674.4162	BKU64162	
Replacement motor Laboratory Drive	0.674.1261	BKU61261	
Replacement Spring Washer Laboratory Drive	0.222.4045	BKU64045	
Bearing for Cap End High Precision Bearing Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)	220-0158	S692-ZZS-MC34-GR5-R-G48	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6) 2 Required	220-0113	S688B-ZZ-MC3-GR3-T9H-G36	





K9 Blau

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Double shielded P6 precision stainless steel bearing (5x11x4)	0.220.0223	SMR115-ZZ-MC3-JP6	
Replacement Bearing Set, Set - to check	1.000.3586	BKU23586	




Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

K9 Blue

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Dental Motor	0.674.4811	BKU61261	
Replacement Chuck True running accuracy less than 0.03, bore diameter 2.35mm	0.674.4871	BKU64871	
Bearing for Shaft High Speed Bearing - Polyamide Retainer 7x14x1x7.5, stainless steel, Polyamide retainer - Checking Details	0.220.0055	BKU20055	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (7.94x15.88x6)	0.220.0044	SR11010A-ZZS-MC3-GR3-T9H	




K9 EWL 950

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Rod Laboratory Drive	0.674.1191	BKU61191	
Replacement Cap Laboratory Drive	0.674.2861	BKU62861	
Replacement Screw nut Laboratory Drive	0.674.2871	BKU62871	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production




K9 Umbra

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	0.220.0157	SMR104A-ZZ-MC3-JP6	
Bearing for Shaft 8-16-6mm High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)	0.220.0113	S688B-ZZ-MC3-GR3-T9H-G36	 Image Coming Soon
Bearing for Rotor High Precision Bearing Shielded Option Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)	0.220.0158	S692-ZZS-MC34-GR5-R-G48	 Image Coming Soon

K9 Umbra EWL 955

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)	0.220.0148	SMR148A-ZZS-MC3-GR3-T9H	




KL702

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	1.003.6023	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Shaft High Precision Bearing Shielded Option Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)	1.003.5818	S693-ZZS-MC34-TH5-R-L23	 Image Coming Soon
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	1.003.6023	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

KL703 LED

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x7x3)	1.007.0966	S683A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	1.003.6023	SMR85C-ZZS-MC34-GR5-TR-L23	
Bearing for Shaft High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	1.003.6023	SMR85C-ZZS-MC34-GR5-L23	 Image Coming Soon





L80

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	=L890:L891+U890:U891=C874	BKU63606	
Head Gear Complete 1:1 For Handpiece with Wedge Type Chuck	1.009.1782	BKU61782	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

M05L

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.011.1374	BKU61374	
Replacement Middle Gear Complete 1:5	1.011.3026	BKU63206	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Precision ceramic-hybrid dental handpiece bearing with Torlon retainer (3.18x5.79x2270)	1.010.9766	SMR5.79-3.175-CSN-ZZS-MC3-TO	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.5x6.5x2)	1.010.9767	SMR6.50-3.50A-CSN-ZZS-MC3-TO	

M25

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	1.011.0129 + 1.001.3515 + 1.001.3519	BKU62919	
Drive Kit	1.011.0129 + 1.001.3515	BKU62915	
Replacement Head Housing	1.009.3635	BKU63635	
Replacement Sprayinset	1.009.3629	BKU63629	
Replacement Shaft	1.001.3518	BKU63518	




Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

M25L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	1.009.3634 + 0.201.8157 + 1.009.3638	BKU63634	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	1.011.0129	BKU60129	
Replacement Dog Gear 1:5 Gear	1.001.3071	BKU63071	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.5x6.5x2)	1.010.9767	SMR6.50-3.50A-CSN-ZZS-MC3-TO	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Precision ceramic-hybrid dental handpiece bearing with Torlon retainer (3.18x5.79x2270)	1.010.9766	SMR5.79-3.175-CSN-ZZS-MC3-TO	

M8700L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	3.001.9915	BKURC210020	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	3.001.9915	BKURC2100203	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	3.001.9915	BKURC2100205	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

M8700L Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Application O-Ring	3.002.0256	BKU60256	
M8900L			
Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKUR2100193	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKUR2100194	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2.000.2266	BKUR2100195	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2.000.2266	BKUR210019	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2.000.2266	BKURC2100195	
Replacement Press Button	3.000.1874	BKU61874	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKURC2100193	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2.000.2266	BKURC210019	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

M9000L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2.000.2266	BKURC2100195	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKURC2100193	
Replacement Press Button	3.000.1874	BKU61874	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2.000.2266	BKUR210019	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKUR2100194	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2.000.2266	BKUR2100195	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2.000.2266	BKURC210019	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2.000.2266	BKUR2100193	
M9000L Air Driven Handpiece			
Replacement O-Ring	3.000.2581	BKU62581	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production



Magno 634

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required	220-0128 / 0145	SR144-D106-ZS	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-KC	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required	220-0128 / 0145	SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required	220-0128 / 0145	SR144-D106-ZS-TO	
Miscellaneous Refer OEM Part Number			
Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38)	0.220.0128	SR144-D106-CSN-ZS-KC	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38)	0.220.0128	SR144-D106-CSN-ZS-TO	
Bearing for Rotor Cartridge High Precision Bearing Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer	0.220.0128	SR144-D106-ZS-TO	
Bearing for Rotor Cartridge High Precision Bearing Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer	0.220.0128	SR144-D106-ZS-KC	



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production



Multiflex® Coupling

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Small - Silicone Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required	200-6113 (Part of Set)	OR-00370-099-S60-R	
Replacement O-Ring Large Special sized o-ring used in dental Applications 3 Required	200-6113 (Part of Set)	OR-00640-080-N70	


Multiflex® Coupling Motor K190

Position / Description / Application	Manufacturers Part	MBA Part	
Xenon Lamp	0.553.3881	BKU63881Xenon	
LED Bulb	0.553.3881	BKU63881LED	

Multiflex® Coupling Motor K201

Position / Description / Application	Manufacturers Part	MBA Part	
LED Bulb	0.553.3881	BKU63881LED	
Xenon Lamp	0.553.3881	BKU63881Xenon	

S201 C

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	1.001.4490 + 1.001.4501	BKU69001	




Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production




S201XL

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 20:1 For Handpiece with Latch Type Chuck	3.000.5177	BKU65177	
Replacement Head Housing	3.000.5315	BKU65315	

S459C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.001.2518	BKURC2100223	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.001.2518	BKURC2100225	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.001.2518	BKURC210022	

S459L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.001.2518	BKURC210022	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.001.2518	BKURC2100223	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.001.2518	BKURC2100225	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

S550

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

1.000.8101

S693-CSN-ZBS-MC34-GR7-L23

High Precision Bearing

Double shielded with special bronze shield one side. Abec 7 precision stainless ceramic-hybrid bearing (3x8x3)



S605

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

1.000.7979

SWR144-D28W-CSN-Z-KC

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded dental special ceramic-hybrid bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)



Replacement Spring Washer

1.001.5032

BKU65032



Replacement Press Button

1.007.7277

BKU67277



Bearing for Rotor Cartridge

1.000.7979

SWR144-D28W-CSN-Z-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Single shielded dental special angular contact ceramic-hybrid bearing with extended inner ring and high speed Torlon retainer (3.18x6.35x2.38)



Bearing for Rotor Cartridge

SG1304X

SWR144-D28W-Z-KC

High Speed Bearing - Phenolic Retainer




Single shielded dental special stainless steel bearing with extended inner ring and high speed phenolic retainer (3.18x6.35x2.38)





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

S605C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.008.4164	BKURC210021	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.008.4164	BKURC2100215	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4164	BKURC2100213	

S609

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	1.007.0854	BKU60854	
Replacement Spring Washer	1.001.5032	BKU65032	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>





S609C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.008.4163	BKUR210016	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR2100163	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR2100164	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR210016T	 Triple pack - Unit with 3 pieces
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly without ball bearings	1.008.4163	BKURC210016	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.008.4163	BKUR2100165	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

S615L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4164	BKURC2100213	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.008.4164	BKURC2100215	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.008.4164	BKURC210021	
Replacement Press Button	1.007.7277	BKU67277	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

S619C

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	1.008.4163	BKUR2100165	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly without ball bearings	1.008.4163	BKURC210016	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	1.008.4163	BKUR210016	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR2100163	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR210016T	 <p style="font-size: small;">Triple pack - Unit with 3 pieces</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	1.008.4163	BKUR2100164	

S619L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	1.007.0854	BKU60854	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

SF4004

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

0.671.1562

BKU61562

True running accuracy less than 0.03, bore diameter 2.35mm



Bearing for Shaft

0.247.6056

S608-ZZS-MC3-GR3-KC

High Speed Bearing - Phenolic Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (8x22x7)



Bearing for Shaft 7-14-5mm

0.220.0112

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)



SF4005

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

0.671.1562

BKU61562

True running accuracy less than 0.03, bore diameter 2.35mm



Bearing for Shaft

0.247.6057

S608-ZZS-MC3-GR3-KC

High Speed Bearing - Phenolic Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (8x22x7)



SONICflex 2000 Airscaler

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Bush

0.571.5032

BKU65032S

Refer OEM Part Number



SONICflex 2003 Airscaler

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Bush

0.571.5032

BKU65032S




Refer OEM Part Number



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Surgical Instrument 3554

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Centre High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-TR-L23	
Lower Bearing 4x9x4 Stainless Steel Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	220 0103	S684A-ZZ-MC3-JP6	
Bearing for Centre High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	0.220.1302	SMR85C-ZZS-MC34-GR5-L23	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Colibri

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Kool Torque®

USA

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

BB 6-13-5

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)

2 Required



**Image
Coming
Soon**

Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

330E

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-ZS-KC-OPT



557 CX ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

10632

BKUR360001

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



557 Push Button

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x1.91)

SR144C-D26-KC



Replacement O-Ring

Special Dental Application O-Ring

Special sized o-ring used in dental Applications

2 Required

OR-00600-075-N70



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-ZS-KC-OPT



557 Turbo ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

10633

BKUR360003

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



557SR

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-ZS-KC-OPT



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

757

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x1.91)

SR144C-D26-KC



Replacement O-Ring

Special Dental Application O-Ring

Special sized o-ring used in dental Applications

2 Required

OR-00600-075-N70



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-ZS-KC-OPT



757 CX ProStyle SF

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

10629

BKUR360002

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



757 Multi Port

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

10630

BKUR360004

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



757 Single Port

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

10631

BKUR360005

For Handpiece with Clamp Pad Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Key




Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Workhorse

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x1.91)		SR144C-D26-KC	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required		OR-00600-075-N70	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-ZS-KC-OPT	 Image Coming Soon

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

BM50S1

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

SUS420J2

BKU64202

True running accuracy less than 0.03, bore diameter 2.35mm



H35LSP

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60108



H35SP1


Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61107



H37L


Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60108



Handy 700

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

SUS420J2

BKU64202

True running accuracy less than 0.03, bore diameter 2.35mm







Key	Usually stocked or Short Delivery
	** Available while stocks last



Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Marathon®


Motor DC Type

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded bearing with a high speed polyamide retainer and light oil lubricant (3x10x4)	1030ZZ	S623-ZZS-MC3-T9H-AF12	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	830ZZ	S693A-ZZS-MC3-T9H	
Bearing for Shaft 6-13-5mm High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)	1360ZZ	S686B-ZZS-MC3-GR3-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)	188ZZ	SR188A-ZZS-MC3-GR3-T9H	

Motor E Type

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	830ZZ	S693A-ZZS-MC3-T9H	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded bearing with a high speed polyamide retainer and light oil lubricant (3x10x4)	1030ZZ	S623-ZZS-MC3-T9H-AF12	

Motor ETC Type

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft 6-13-5mm High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)	1360ZZ	S686B-ZZS-MC3-GR3-T9H	

Multi 600

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Chuck True running accuracy less than 0.03, bore diameter 2.35mm	SUS420J2	BKU64202	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Marathon®

SDEH37L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61107



SDESH40C

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61108



SH40C


Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60108



SH45C

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Armature

BKU61108



Replacement Motor Brush

BKU60108



Spindle DC Type

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

1480ZZ

SMR148A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)



Bearing for Shaft

1260ZZ

SMR126A-ZZS-MC3-T9H

High Speed Bearing - Polyamide Retainer Shielded Option



Double shielded high speed stainless steel bearing with Polyamide retainer (6x12x4)



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Spindle ETC Type

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (6x12x4)	1260ZZ	SMR126A-ZZS-MC3-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)	1480ZZ	SMR148A-ZZS-MC3-GR3-T9H	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Micro Mega®

350

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SFR144A-D27-ZSO-KC-OPT

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

1000

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



400

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-T9H-OPT



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-KC-OPT



Bearing for Chuck End

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-T9H-OPT



Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-KC-OPT



Bearing for Cap End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-CSN-KC-OPT



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

5 Hole Hose

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

790257R

BKU60257Xenon



5 Hole Hoses

Position / Description / Application

Manufacturers Part

MBA Part

LED Bulb

790257R

BKU63307LED



6 Pin Coupling

Position / Description / Application


Manufacturers Part

MBA Part

Xenon Lamp

790257R

BKU60257Xenon



LED Bulb

790257R

BKU63307LED



8000 - 8000i

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-KC-OPT



Bright ISO

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

262522

BKU62522Xenon



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

C6

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

262522

BKU62522Xenon



Easycare Mini

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1223X-D

SR145-D75-CSN-ZZE-KC

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Shielded Option

Double shielded precision special ceramic-hybrid dental bearing with stepped OD and high speed Phenolic retainer



Bearing for Rotor Cartridge

SG1223X.2mq

SR145-D75-ZZE-KC

High Speed Bearing - Phenolic Retainer Shielded Option

Double shielded Abec 7 precision bearing with stepped OD and high speed Phenolic retainer



eStylus Motor

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

774150

BKU66016LED



Quiet Air

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-T9H-OPT



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-KC-OPT



Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-KC-OPT



Bearing for Cap End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)





SR144-D28-CSN-KC-OPT





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Quiet Air Compact

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	
Bearing for Cap End High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-CSN-KC-OPT	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-KC-OPT	





Stylus

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Shielded Option Double shielded Abec 7 precision bearing with stepped OD and high speed Phenolic retainer	SG1223X.2mq	SR145-D75-ZZE-KC	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Shielded Option Double shielded precision special ceramic-hybrid dental bearing with stepped OD and high speed Phenolic retainer	SG1223X-D	SR145-D75-CSN-ZZE-KC	




Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


Tradition

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-KC-OPT	
Bearing for Cap End High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-CSN-KC-OPT	

Tradition Lever

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece ceramic-hybrid bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)	SG1222X-D	SR144-D28-CSN-ZS-KC	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	780118	BKURMU260101	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	780118	BKURM260101	

Tradition Lever Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge Daily Lube - High Speed Bearing - Phenolic Retainer Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)	SG1301X-D	SFR144A-D1301-CSN-ZSO-KC	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Traditon XGT

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

790118

BKURMU260100

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Turbine Sub Assembly Without Bearings

790118

BKURM260100

For Handpiece with Wedge Type Chuck



Traditon Lever

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1222X

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)



Traditon Lever Push Button

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1301X

SFR144A-D1301-ZSO-KC

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)



Tru-Line

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-T9H-OPT



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

SFR144A-D27-ZSO-KC-OPT



Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-KC-OPT



Bearing for Cap End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-CSN-KC-OPT



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

XGT Lever Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

BKU69998

Special Dental Application O-Ring



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Miscellaneous

Dental Washers

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)

W0048-T-007-008-W3-S2



PTFE Washer-Spacer

0.13mm - half thou - thick

Deleted Item. Available while stocks last. (3x6.56x0.13)

W0030-F-007-001-PT



Key




Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

EM1190K

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)	BE5015	SMR85C-ZZS-MC34-GR5-TR-L23	
Replacement Motor Brush	SP1219	BKU60881	
Bearing for Rotor High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)	BE5015	SMR85C-ZZS-MC34-GR5-L23	

HB21

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1034	BKUR2801085	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1034	BKUR2801083	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2034	BKURC280003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2034	BKURC2800033	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1034	BKUR280108	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2034	BKURC2800035	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HC20

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2038	BKURC280007	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2038	BKUR2801205	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2038	BKUR2801203	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC2038	BKUR280120	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2038	BKURC2800075	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2038	BKURC2800073	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HC21

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2038	BKURC280007	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC2038	BKUR280120	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2038	BKUR2801203	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2038	BKUR2801205	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2038	BKURC2800073	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2038	BKURC2800075	

Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

HC4021

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2031	BKURC2800013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2031	BKURC2800015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1031 / TC2031	BKUR2801065	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801063	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2031	BKURC2800001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801064	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1031 / TC2031	BKUR280106	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HC4022

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801074	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801073	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1032 / TC2032	BKUR280107	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2032	BKURC2800023	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2032	BKURC280002	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1032 / TC2032	BKUR2801075	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2032	BKURC2800025	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

HC5021

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2031	BKURC2800015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR280106T	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1031 / TC2031	BKUR280106S	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2031	BKURC2800013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801063	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801064	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1031 / TC2031	BKUR280106	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2031	BKURC280001	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HC5022

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801074	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801073	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1032 / TC2032	BKUR280107	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2032	BKURC2800025	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR280107T	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators Triple pack - Unit with 3 pieces 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1032 / TC2032	BKUR2801075	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2032	BKURC2800023	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2032	BKURC280002	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

HC60012

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801074	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2032	BKURC2800025	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1032 / TC2032	BKUR2801075	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1032 / TC2032	BKUR280107	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2032	BKURC2800023	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2032	BKURC2800002	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR280107T	 <p>Triple pack - Unit with 3 pieces</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1032 / TC2032	BKUR2801073	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®




HC6021

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR280106T	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801063	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2031	BKURC2800015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1031 / TC2031	BKUR2801065	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2031	BKURC280001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TC1031 / TC2031	BKUR2801064	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2031	BKURC2800013	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1031 / TC2031	BKUR280106	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

HC8022

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2042	BKURC280006	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2042	BKURC2800065	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2042	BKURC2800063	

HC9021

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2041	BKUR2801105	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2041	BKURC2800055	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2041	BKURC2800053	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2041	BKURC280005	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC2041	BKUR280110	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2041	BKUR2801103	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

HE17

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2034	BKURC2800033	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1034	BKUR2801085	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2034	BKURC2800035	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2034	BKURC280003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1034	BKUR2801083	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1034	BKUR280108	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HE18

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2035	BKURC2800043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1035	BKUR2801095	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1035	BKUR2801093	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2035	BKURC2800045	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1035	BKUR280109	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2035	BKURC280004	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

HE22

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC1035	BKUR2801095	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC1035	BKUR2801093	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC1035	BKUR280109	Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2035	BKURC2800045	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2035	BKURC2800043	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2035	BKURC280004	Precision sub-assemblies, cartridges with carbide inserts & WCC coated actuators 

Key **Usually stocked or Short Delivery**
 ** Available while stocks last




Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®

HP21

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2041	BKUR2801105	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2041	BKURC280005	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2041	BKURC2800053	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TC2041	BKUR280110	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2041	BKUR2801103	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2041	BKURC2800055	

HP22

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TC2042	BKURC2800063	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TC2042	BKURC280006	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TC2042	BKURC2800065	

Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

MK-Dent®

LC01L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5015

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option
Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)



Bearing for Shaft

BE5014

S684A-ZZ-MC3-JP6

High Precision Bearing Shielded Option
Double shielded P6 precision stainless steel bearing (4x9x4)



Bearing for Shaft

BE5015

SMR85C-ZZS-MC34-GR5-L23

High Precision Bearing Standard Cage Option
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



LC11L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5013

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing
Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



LC27L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5018

SMR104A-ZZ-MC3-JP6

High Precision Bearing
Double shielded P6 precision stainless steel bearing (4x10x4)



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities




** Available while stocks last

Not in current production

LCH11

Position / Description / Application	Manufacturers Part	MBA Part	
Intermediate Piece Dental Head 1:1	SP2024	BKU62024	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	BE5012	SMR74-MC3-JP6-AF12	
Drive Kit	SP2005 + SP2024	BKU60524	
Head Gear Complete Dental Head 1:1 For Handpiece with Latch Type Chuck	SP2005	BKU62005	
Bearing for Headgear High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	BE5011	SWA683B-ZZ-C69-GR5-TO	

LCH21

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	SP2105 + SP2124	BKU60254	
Head Gear Complete Dental Head 2:1 For Handpiece with Latch Type Chuck	SP2105	BKU62105	
Intermediate Piece Dental Head 2:1	SP2124	BKU62124	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

LCH21P

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

BE6055

SMR63-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



LE11

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5013

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



LE148L

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Shaft

BE5019

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Bearing for Shaft

BE5017

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



LE15L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

BE5027

SR133A-ZZS-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)



Bearing for Headgear

BE5029

SWR144C-D38-CSN-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

3.175x6.35x1.905/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x1.9)



Bearing for Headgear

BE5021B

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



Bearing for Middle Gear

BE5027

SMR4.762-2.35A-ZZS-MC3-GR7-TO

High Speed Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

LE21

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5013

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



LE27L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5018

SMR104A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



Bearing for Shaft

BE5016

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



LE74L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5019

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Bearing for Shaft

BE5017

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



LP01L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5015

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option

Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)



Bearing for Shaft

BE5014

S684A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x9x4)



Bearing for Shaft

BE5015

SMR85C-ZZS-MC34-GR5-L23

High Precision Bearing Standard Cage Option

Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

LP11L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5013

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



LP15L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

BE5027

SR133A-ZZS-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x4.76x2.38)



Bearing for Middle Gear

BE5027

SMR4.762-2.35A-ZZS-MC3-GR7-TO

High Speed Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision bearing with high speed Torlon retainer (2.35x4.76x2.35)



Bearing for Headgear

BE5029

SWR144C-D38-CSN-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

3.175x6.35x1.905/2.778, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x1.9)



Bearing for Headgear

BE5021B

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



LP54L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5018

SMR104A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



Bearing for Shaft

BE5016

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer Shielded Option

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



LP74L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

BE5016

SMR104A-ZZS-MC3-GR7-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 7 precision stainless steel bearing with high speed Polyamide retainer (4x11x4)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

LPH11

Position / Description / Application	Manufacturers Part	MBA Part	
Intermediate Piece Dental Head 1:1	SP2024	BKU62024	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	BE5012	SMR74-MC3-JP6-AF12	
Drive Kit	SP2005 + SP2024	BKU60524	 <p>Drive Kits for Electric Driven Handpieces and Heads</p>
Head Gear Complete Dental Head 1:1 For Handpiece with Latch Type Chuck	SP2005	BKU62005	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Bearing for Headgear High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	BE5011	SWA683B-ZZ-C69-GR5-TO	 <p>MBA Image Coming Soon</p>

Key **Usually stocked or Short Delivery**
** Available while stocks last



Long Delivery or Large Minimum Quantities
Not in current production

MK-Dent®


LPH21

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	BE5012	SMR74-MC3-JP6-AF12	
Head Gear Complete Dental Head 2:1 For Handpiece with Latch Type Chuck	SP2105	BKU62105	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Intermediate Piece Dental Head 2:1	SP2124	BKU62124	 <p>Head Gears, Intermediate Pieces, Middle Gears, Parts for Electric Driven Handpieces & Heads</p>
Drive Kit	SP2105 + SP2124	BKU60254	 <p>Drive Kits for Electric Driven Handpieces and Heads</p>
Bearing for Headgear High Speed Bearing - Torlon Retainer Shielded Option Double shielded special dental bearing with extended inner ring and high speed Torlon retainer (3x7x2.5)	BE5011	SWA683B-ZZ-C69-GR5-TO	 <p>Image Coming Soon</p>

LS2115L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	BE5021A	SMR74-CSN-TH7-TO	 <p>Image Coming Soon</p>
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (2.38x4.76x1588)	BE5022	SR133-MC3-TH7-TO	 <p>Image Coming Soon</p>

Miscellaneous Refer OEM Part Number

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)	BE5010	SR1-5A-ZZS-MC3-JP6	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

QC6016K

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb and Cap

BU8012

BKU68012LED



QC6016KW

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb and Cap

BU8012

BKU68012LED



Key



Usually stocked or Short Delivery

** Available while stocks last



Long Delivery or Large Minimum Quantities

Not in current production



Falcon

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Chuck End High Precision Bearing Shielded Option Deleted Item. Available while stocks last. (6x10x3) 3 Required		SMR106A-ZZS-MC34-P4P-R-AF12	
Bearing for Chuck End High Precision Bearing Shielded Option Deleted Item. Available while stocks last. (6x10x3) 3 Required		SMR106A-ZZS-MC34-GR7-R-G48	

Micromotor

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (3x9x4)	BB 3-9-4	SMR93A-ZZ-MC3-JP6	
Replacement Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (2.5x8x4)	BB 2.5-8-4	S602XA-ZZ-MC3-JP6	

PAR4HEXO

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)	SG1357X-D	SR144-D28-CSN-ZE-TO-OPT	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-TO	

PAR4HUEXO

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-TO	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)	SG1357X-D	SR144-D28-CSN-ZE-TO-OPT	 Image Coming Soon

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Torx 1000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

Stainless Bearing Shielded Option

Special anti-corrosion stainless TBN-6H bearing designed primarily for fishing reels and salt-water applications (5x8x2.5)

2 Required

SACB-MR85A-ZZ-MC3



Replacement Bearing

Stainless Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (5x8x2.5)

2 Required

SMR85A-ZZS-MC3-JP6



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

MTC

105

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Lynx

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement O-Ring Special Dental Application O-Ring Special sized o-ring used in dental Applications 2 Required	OD of Bearing	OR-00600-100-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Micromat-Athlet Type 200 Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

0229

**** S687A-ZZ-MC3-2OR**

Stainless Bearing

Deleted Item. Available while stocks last. Special bearing with 2 o-rings in OD grooves. (7x14x5)



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

3 405

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



450L

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NCH-SU03

BKUR250059

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



455L Phatelus 45S

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NCH-TU03

BKUR250060

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



A200 N

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

NMC-MU03

BKUR250052

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



A300 N

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

NMC-SU03

BKUR250053

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



A400 N

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NMC-TU03

BKUR250054

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

A500

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

Ti-MU03

BKUR250045

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



A600

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

Ti-SU03

BKUR250046

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



A700

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

Ti-TU03

BKUR250047

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Key






Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production










Alpha Air

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Alpha Light

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)	418WIF	SFR144A-D27-ZSO-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)	418WIF	SFR144A-D27-ZSO-T9H-OPT	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)	418WIF-CER	SFR144A-D27-CSN-ZSO-KC-OPT	 Image Coming Soon

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
 Not in current production

Bearing for Rotor Cartridge

418WIF-CER

SFR144A-D27-CSN-ZSO-T9H-OPT

High Speed Ceramic-Hybrid Bearing - Polyamide Retainer

Precision flanged dental handpiece ceramic-hybrid bearing with polyamide retainer (3.18x6.35x2.78)



BB Y

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C032012

BKU62012



Head Gear Complete 1:1

C037015

BKU67015



Drive Kit

C037015 + C032012

BKU61512



CHQD

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-CSN-KC



Bearing for Chuck End

High Speed Bearings - Polyamide Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-ZSO



Bearing for Chuck End

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

SFR144-D28-KC



Bearing for Cap End

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

SR144-D28-ZS-KC-OPT



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

Coupling

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

Small
Special sized o-ring used in dental Applications
2 Required

OR-00460-070-N70



Replacement O-Ring

Large
Special sized o-ring used in dental Applications
4 Required

OR-00635-089-N70



Coupling KCL

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

10000198

BKU60198LED



Coupling KCL SCL

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

10000198

BKU60198LED



DYNAled

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic

SX-SU03

BKUR250073



E4R Slow Speed - Contra Angle Green Band

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Cap End

Stainless Bearing
Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)

SMR63-MC3-JP6-AF12



Bearing for Chuck End

Stainless Bearing
Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)

SMR105-MC3-JP6-AF12



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

EX100

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Shaft

B22227-7923

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



EX103

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B22227-7923

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



EX12

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B10321-1005

SMR105-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



EX200

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B22227-7923

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



EX5

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B10321-1005

SMR105-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



Bearing for Shaft

B13213-0804

SMR84A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x8x3)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



** Available while stocks last

Not in current production

EX6

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	B10321-1005	SMR105-MC3-JP6-AF12	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	B13213-0804	SMR84A-ZZ-MC3-JP6	



EX6 Blue Band

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Stainless Bearing Double shielded P6 precision stainless steel bearing (4x8x3)	BB 4-8-3	SMR84A-ZZ-MC3-JP6	
Replacement Bearing Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	BB 5-10-3	SMR105-MC3-JP6-AF12	

EX7

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	B10321-1005	SMR105-MC3-JP6-AF12	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	B13213-0804	SMR84A-ZZ-MC3-JP6	

EX8

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	B13213-0804	SMR84A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	B10321-1005	SMR105-MC3-JP6-AF12	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

EX9

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B10321-1005

SMR105-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



Bearing for Shaft

B13213-0804

SMR84A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x8x3)



FPB

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

C875050

BKU65050



FPB Y

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C875012

BKU65012



Drive Kit

C875060 + C875012

BKU66012



FPB03

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

C875060

BKU65060



FX15

Position / Description / Application


Manufacturers Part

MBA Part

Drive Kit

C1052700 + C10000946

BKU60046



Key




Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities





** Available while stocks last

Not in current production

FX23

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Middle Gear Complete 1:1	10000946	BKU60946	
Replacement Press Button	C1052500 - FX25-04	BKU62500X	
Head Gear Complete 1:1	C1052700	BKU62700	





FX25

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	C1052500 - FX25-04	BKU62500X	
Drive Kit	C1052700 + C10000946	BKU60046	
Replacement Middle Gear Complete 1:1	10000946	BKU60946	
Head Gear Complete 1:1	C1056700	BKU62700	


Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production



FX57

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement End Cap 4:1	20002065	BKU62065	
Head Gear Complete 4:1	10001571	BKU61571	
Replacement End Cap 4:1	20002066	BKU62066	
Replacement Middle Gear Complete 4:1	10000996	BKU60996	

IS205

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)	B22227-7923	SR1-5A-ZZS-MC3-JP6	

IS65

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)	B10321-1005	SMR105-MC3-JP6-AF12	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x8x3)	B13213-0804	SMR84A-ZZ-MC3-JP6	



M25

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	C601050	BKU61050	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


M25L New

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Middle Gear Complete 1:1	C1142022	BKU62022	 Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces & Heads
Head Gear Complete 1:1	C1038015	BKU68015	 Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces & Heads


M25L Old

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Middle Gear Complete 1:1	10000574	BKU60574	 Head Gears, Intermediate Pieces, Middle Gears Parts for Electric Driven Handpieces & Heads


M500

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	P1001500 - SX-MU04	BKU61500	 Press Button cpl. for Heads & Air Driven Handpieces


M500 Pana Max M Plus

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	SX-MU03	BKUR250048	 Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators

M600

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	P1002500 - SX-SU04	BKU62500M	 Press Button cpl. for Heads & Air Driven Handpieces

M65

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Chuck True running accuracy less than 0.03	H190180	BKU60180	 Chucks for Laboratory Drives, Electric Handpieces




Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

M95

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	C600015	BKU60015	
Replacement Press Button	10000664	BKU60664	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427062	BKU67062	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427168	BKU67168	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	740-CER	SMR74-CSN-TH7-TO	 Image Coming Soon

M95L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	10000664	BKU60664	
Drive Gear Assembly - Dog Clutch Complete 1:5 Complete Assembly	10000560	BKU60560	
Dog Gear 1:5 Dog Clutch Clutch	C600C013	BKU69993	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

M95L New

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:5

C1141022

BKU61022N



M95L Old

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:5

10000559

BKU60559



Mach - Mach-LITE P045SU03 - P046TU03 Green-Black

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00620-100-N70



Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)

W0048-T-007-008-W3-S2



Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Mach Lite 65 T

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NCH-TU03

BKUR250060

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Mach Lite XT M

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NMC-MU03

BKUR250052

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Mach Lite XT S

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NMC-SU03

BKUR250053

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Mach Lite XT T

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

NMC-TU03

BKUR250054

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Mach Mini- Mach-LITE Mini P047MU03 Red Band

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

SFR144-D28-CSN-KC



Replacement Bearing

High Speed Bearing - Polyamide Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

SFR144-D28-ZSO



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Deleted Item. Available while stocks last. (3.18x6.35x2.38)

2 Required

SFR144-D28-KC



Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00620-100-N70



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Mach N75

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

418NA-CER

SR144-D28-CSN-KC-OPT

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)



Bearing for Rotor Cartridge

418NA

SR144-D28-T9H-OPT

High Speed Bearing - Polyamide Retainer

Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)



Bearing for Rotor Cartridge

418NA

SR144-D28-KC-OPT

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)



Bearing for Rotor Cartridge

418NA-CER

SR144-D28-CSN-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)



Mach QD 45

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

NCH-SU03

BKUR250059

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer



Motor KCL Coupler

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

Y900132

BKU60132Xenon



Motor Ultimate XL 450 UHR

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Motor Shaft

B11253-0904

S684A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (4x9x4)



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



** Available while stocks last

Not in current production



MP F

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C521015 + C547022	BKU61522	
Head Gear Complete 1:1	C521015	BKU61015M	

MP Y

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Middle Gear Complete 1:1	C547022	BKU67022	
Head Gear Complete 1:1	C521015	BKU62101	




MPA Y

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1	C521015	BKU61015M	
Drive Kit	C521015 + C547022	BKU61522	


Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production



NAC Y

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C030011 + C032012	BKU61112	
Head Gear Complete 1:1	C030011	BKU60011	
Replacement Middle Gear Complete 1:1	C032012	BKU62012	

Nano 25LS

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C1038015 + C1038022	BKU61580	

Nano 95LS

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	C1034017	BKU64017	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	C1034015	BKU64015	

Key





Usually stocked or Short Delivery

** Available while stocks last





Long Delivery or Large Minimum Quantities

Not in current production

NL45

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-TO	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)	418NA	SR144-D28-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)	418NA	SR144-D28-T9H-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-KC-OPT	





NL65 CH QD

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Open high precision stainless bearing for dental applications with high speed polyamide retainer (3.18x6.35x2.38)	418NA	SR144-D28-T9H-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)	418NA	SR144-D28-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer Precision ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)	418NA-CER	SR144-D28-CSN-TO	




Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

NL9000L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TiX-TU03	BKUR2500445	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TiX-TU03	BKUR250044	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR2500444	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR2500443	

NL9000M

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-MU03	BKUR2500423	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TiX-MU03	BKUR2500425	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TiX-MU03	BKUR250042	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-MU03	BKUR2500424	

Key

Usually stocked or Short Delivery
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


NL9000S

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR2500434	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TiX-SU03	BKUR250043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TiX-SU03	BKUR2500435	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR250043T	 Triple pack - Unit with 3 pieces
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR2500433	

NL9000S Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring	0312-063100	BKU63100	

NL9000T

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR250044T	 Triple pack - Unit with 3 pieces

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

NL95S

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

Ti-SU03

BKUR250046

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic




Replacement Turbine With Bearings

Ti-MU03

BKUR250045

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



NL95T

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

Ti-TU03

BKUR250047

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



NPB Slow Speed

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

BB 4-8-2

SMR84-MC3-JP6-AF12

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)



Replacement Bearing

BB 4-7-2

SMR74-MC3-JP6-AF12

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)



NSK Presto Rotor

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

418WIF

SFR144A-D27-ZSO-KC-OPT

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)



Bearing for Rotor Cartridge

418WIF-CER

SFR144A-D27-CSN-ZSO-KC-OPT

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

NTF TU03

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production







Optica

Position / Description / Application	Manufacturers Part	MBA Part	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)		SR144A-D27-ZS-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	


Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Pana Air Mini MU03 Red Band

Position / Description / Application	Manufacturers Part	MBA Part	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)		SR144A-D27-ZS-T9H-OPT	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	





Pana Air Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer	NPA-SU03	BKUR250056	 <p>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>


Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Pana Air SU03 - TU03 Green-Black

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)		SR144A-D27-ZS-T9H-OPT	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)		SFR144-D28-KC	
Bearing for Chuck End High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)		SFR144-D28-CSN-KC	
Bearing for Chuck End High Speed Bearings - Polyamide Retainer Deleted Item. Available while stocks last. (3.18x6.35x2.38)		SFR144-D28-ZSO	




Pana Max 2

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer	PAX2-SU03	BKUR250068	 <p>Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Press Button	P1190050 - PAX2-SU03	BKU60050	 <p>Press Button cal. for Heads & Air Driven Handpieces</p>


Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production




Pana Max S

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings	PAX-SU	BKUR250051	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	PAX-SU	BKUR2500513	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	PAX-SU	BKUR2500514	

Pana Max S Plus

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	SX-SU03	BKUR250073	





Pana Max T

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	PAX-TU	BKUR2500503	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	PAX-TU	BKUR2500504	
Replacement Turbine Excluding Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings	PAX-TU	BKUR250050	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>







Phatelus

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)	418WIF	SFR144A-D27-ZSO-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)	418WIF	SFR144A-D27-ZSO-T9H-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Polyamide Retainer Precision flanged dental handpiece ceramic-hybrid bearing with polyamide retainer (3.18x6.35x2.78)	418WIF-CER	SFR144A-D27-CSN-ZSO-T9H-OPT	 Image Coming Soon
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision dental handpiece bearing ceramic-hybrid and phenolic retainer (3.18x6.35x2.78)	418WIF-CER	SFR144A-D27-CSN-ZSO-KC-OPT	 Image Coming Soon

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Phatelus II

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)		SR144A-D27-ZS-T9H-OPT	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production


Phatelus III

Position / Description / Application	Manufacturers Part	MBA Part	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Bearing for Chuck End High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)		SR144A-D27-ZS-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Phatelus MU03

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00607-079-N70	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	

Phatelus Set Air Driven Handpiece

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring	Y900580	BKU60580	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Phatelus SU03

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)

W0048-T-007-008-W3-S2




PTFE Washer-Spacer

0.13mm - half thou - thick

Deleted Item. Available while stocks last. (3x6.56x0.13)

2 Required

W0030-F-007-001-PT



Replacement O-Ring

Quad Ring

High Temperature

Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications

2 Required

**** QR-00600-080-N70-NSK**



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-KC-OPT



Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00607-079-N70



Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-T9H-OPT



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Phatelus TU03

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00607-079-N70



Replacement Bearing

High Speed Bearing - Polyamide Retainer

Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-T9H-OPT



Replacement O-Ring

High Temperature

Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications

2 Required

Quad Ring

** QR-00600-080-N70-NSK



Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)

W0048-T-007-008-W3-S2



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)

2 Required

SFR144A-D27-ZSO-KC-OPT




PTFE Washer-Spacer

0.13mm - half thou - thick

Deleted Item. Available while stocks last. (3x6.56x0.13)

2 Required

W0030-F-007-001-PT



Presto AQUA Cartridge

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Rotor / Cartridge

T816-016

BKUR250067



Presto Aqua II

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

T815019

BKU65019

True running accuracy less than 0.03



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

Presto Aqua Lux

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

T815019

BKU65019

True running accuracy less than 0.03



Presto Cartridge

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Rotor / Cartridge

T815-016

BKUR250066



PTL

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

Y900529


BKU60529LED



Xenon Lamp

Y900529

BKU60529Xenon



PTL 950TI

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B22227-7923

SR1-5A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



QD J

Position / Description / Application

Manufacturers Part


MBA Part

Replacement O-Ring




0312-021083/063080

BKU61083




Special dental applications O-Ring




QD S

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	PAX-SU	BKUR2500514	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	PAX-SU	BKUR2500513	
Replacement Turbine Excluding Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings	PAX-SU	BKUR250051	


QD T

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	PAX-TU	BKUR2500503	
Replacement Turbine Excluding Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly excluding bearings	PAX-TU	BKUR250050	
Replacement Turbine With Bearings For Handpiece with Clamp Pad Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	PAX-TU	BKUR2500504	

S Max M25 NEW

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C1038015 + C1142022	BKU61502	

S Max M25 OLD

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C601015 + 10000574	BKU61574	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

S Max M95L NEW

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C600015 + C1141022

BKU61522



S Max pico

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

SX-PU03

BKUR250065

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



S Max SG20

Position / Description / Application

Manufacturers Part


MBA Part

Drive Gear Assembly - Dog Clutch Complete 20:1

6 Parts

BKU63001

Complete Assembly



Drive Kit

C1010150 + C1010022


BKU65021



Replacement Middle Gear Complete 20:1

C1010022

BKU60022SG



S Max-Reg tm- M500

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

SX-MU03

BKUR250074

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



S Max-Reg tm- M600

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

SX-SU03

BKUR250073

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

S Max-Reg tm- M800

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

SX-MU03

BKUR250074

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



S Max-Reg tm- M900

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

SX-SU03

BKUR250073

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic



S900 Airscaler

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Bush

T695104B

BKU65104

Refer OEM Part Number



S950 Airscaler

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bush

T695104B

BKU65104

Refer OEM Part Number



SG20

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

10000228

BKU60228



SG20L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C1010150 + C1003022


BKU65022



Head Gear Complete 20:1

C1010150002

BKU60002



Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

SG20X

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

10000228

BKU60228



SGAE

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B12213-0704

SMR74A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



SGAE2G

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B12213-0704

SMR74A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



SGM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C293B031

BKU63031



SGMS

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C293B031

BKU63031



SGR2E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B12213-0704

SMR74A-ZZS-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x7x2.5)



Bearing for Shaft

B10313-0703

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

SGSE

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B12213-0704

SMR74A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



SGT2E0

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B12213-0704

SMR74A-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4x7x2.5)



SGT2E1

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B10313-0703

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



TEP Y

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C837012

BKU67012



Head Gear Complete 1:1

C837015

BKU67015TEP



Ti Max® M95L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C600015 + 10000664

BKU61564



Ti Max® Ti25L

Position / Description / Application

Manufacturers Part


MBA Part

Drive Kit

C875060 + C880022

BKU66022

Dental replacement drive kit



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

Ti Max® X

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C1010150 + C1003022

BKU65022



Ti Max® X25L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C601015 + C601022

BKU61522X



Ti Max® X500

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

TiX-MU03

BKUR2500423

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

TiX-MU03

BKUR2500425

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

TiX-MU03

BKUR2500424

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

TiX-MU03

BKUR250042

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



** Available while stocks last

Not in current production

Ti Max® X600

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TiX-SU03	BKUR250043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TiX-SU03	BKUR2500435	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR2500434	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR250043T	 <small>Triple pack - Unit with 3 pieces</small>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKUR2500433	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TiX-SU03	BKURM250043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKURMU250043	

Ti Max® X650

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	TiX-SU03	BKURM250043	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-SU03	BKURMU250043	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Ti Max® X700

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR250044T	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR2500444	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	TiX-TU03	BKUR2500445	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	TiX-TU03	BKUR2500443	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	TiX-TU03	BKUR250044	

Ti Max® X95L

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C600015 + C600022	BKU61522	

Ti Max® Z25

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C1038015 + C1038022	BKU61580	

Ti Max® Z45L

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	C1064015 + C1064022	BKU61522Z	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Ti Max® Z800L

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine Excluding Bearings

TiZ-SU03

BKUR250072

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

TiZ-SU03

BKUR2500723

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Ti Max® Z900L

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

TiZ-SU03

BKUR2500693

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

TiZ-SU03

BKUR250069

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings



Ti Max® Z95

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

C1034D017 + C1034022

BKU61722



Ti15

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C880022

BKU68022



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Ti25

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

C880022

BKU68022



Replacement Press Button

C875050

BKU65050



Ti25L

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

C875060

BKU65060



Ti85

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

B19157-4723

SR133-MC3-TH7-TO

High Speed Bearing - Torlon Retainer

Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (2.38x4.76x1588)



Ti95

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Dog Clutch

SMR84-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)



Bearing for Middle Gear

B19157-4723

SR133-MC3-TH7-TO

High Speed Bearing - Torlon Retainer

Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (2.38x4.76x1588)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

UHR 50T C

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



UHR45C

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 6-13-5mm

B11243-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Bearing for Shaft

B20245-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



UHR50C

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spindle Washer

H063B109A

BKU63109

Laboratory Drive



Replacement Dust Proof Cover

H063B106A

BKU63106

Laboratory Drive



UHR50T

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spindle Washer

H063B109A

BKU63109

Laboratory Drive



Replacement Dust Proof Cover

H063B106A

BKU63106

Laboratory Drive



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

UHRXL T

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B20245-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Bearing for Shaft 6-13-5mm

B11243-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Ultimate 450

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Ultimate 450-VH 2

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Ultimate 450-VH 3

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Ultimate XL

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B20245-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Bearing for Shaft 6-13-5mm

B11243-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Ultimate XL 450

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

H203180

BKU63180

True running accuracy less than 0.03, bore diameter 2.35mm



Ultimate XL 500

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

H203180

BKU63180

True running accuracy less than 0.03, bore diameter 2.35mm



VH 2

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

B20215-2763

SR188A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)



Bearing for Shaft 6-13-5mm

B11223-1306

S686B-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer









Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Vip-Air

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<i>OR-00607-079-N70</i>	
Replacement Bearing High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-KC-OPT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78) 2 Required		SFR144A-D27-ZSO-T9H-OPT	
Replacement O-Ring High Temperature Deleted Item. Available while stocks last. Special quad-ring o-ring for dental applications 2 Required	Quad Ring	** QR-00600-080-N70-NSK	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
VOLVERE 8 30000 RPM			
Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End Stainless Bearing Double shielded P6 precision stainless steel bearing (6.35x12.7x4.76)		SR188A-ZZ-MC3-JP6	
Bearing for Chuck End Stainless Bearing Double shielded P6 precision stainless steel bearing (5x13x4)		S695-ZZ-MC3-JP6	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

VOLVERE GX 40M

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)

SMR84-MC3-JP6-AF12



Bearing for Cap End

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)

SMR74-MC3-JP6-AF12



VOLVERE GX 6000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x6x2.5)

SMR63A-ZZ-MC3-JP6



Bearing for Cap End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x7x3)

S683A-ZZ-MC3-JP6



VOLVERE GX35M VH

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Cap End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (6.35x12.7x4.76)

SR188A-ZZ-MC3-JP6



Bearing for Chuck End

Stainless Bearing

Double shielded P6 precision stainless steel bearing (6x13x5)

S686B-ZZ-MC3-JP6



Volvere i7

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

H203180

BKU63180

True running accuracy less than 0.03, bore diameter 2.35mm



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

VR E

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft 6-13-5mm High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (6x13x5)	B11243-1306	S686B-ZZS-MC3-GR3-T9H	 Image Coming Soon
Bearing for Shaft High Speed Bearing - Polyamide Retainer Double shielded Abec 3 precision stainless steel bearing with high speed polyamide retainer (6.35x12.7x4.76)	B20245-2763	SR188A-ZZS-MC3-GR3-T9H	 Image Coming Soon

X

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	10000228	BKU60228	

X SG20L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Head Housing	10000254001	BKU64001	
Replacement Middle Gear Complete 20:1	C1003022	BKU63022	

X15L

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1	C601015	BKU61015	
Replacement Middle Gear Complete 1:1	C601022	BKU61022	




Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

X25

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427062	BKU67062	
Replacement Head Housing	C601011	BKU61011	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427168	BKU67168	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427064	BKU67064	
Replacement Press Button	C601050	BKU61050	

X25L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Head Housing	C601011	BKU61011	
Head Gear Complete 1:1	C601015	BKU61015	
Replacement Middle Gear Complete 1:1	C601022	BKU61022	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

X450

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

P1078500 - TiX45-SU04

BKU68500



X500

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

P067B050 - TiX-MU04

BKU67050



X600

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

P068B050 - TiX-SU04

BKU68050



Bearing for Rotor Cartridge

SG1441X-D

SR144-D29-CSN-TO-OPT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision angular contact ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)



X600 Air Driven Handpiece

Position / Description / Application


Manufacturers Part

MBA Part

Replacement O-Ring

0312-063100

BKU63100



X65

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Chuck

H190180

BKU60180

True running accuracy less than 0.03



X650

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1441X-D

SR144-D29-CSN-TO-OPT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision angular contact ceramic-hybrid dental handpiece bearing with high speed Torlon retainer. ceramic balls (3.18x6.35x2.38)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

X700

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Press Button

P069B050 - TiX-TU04

BKU69050



X95

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Dog Clutch

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x8x2)

SMR84-MC3-JP6-AF12




Replacement O-Ring

C427064

BKU67064

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring



Replacement O-Ring

C427168

BKU67168

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring



Head Gear Complete

C600015

BKU60015

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Replacement O-Ring

C427062

BKU67062

Special Dental Applications O-Ring for Electric handpiece

Special dental applications O-Ring



Bearing for Headgear

740-CER

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

X95L

Position / Description / Application	Manufacturers Part	MBA Part	
Dog Gear 1:5 Dog Clutch C600C013 Clutch		BKU69993	
Replacement Middle Gear Complete 1:5	C600022	BKU60022	
Replacement Head Housing	C600011	BKU60011A	

XSG20L

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 20:1	C1010150002	BKU60002	

Z25

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427168	BKU67168	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	C427062	BKU67062	
Replacement Press Button	C1038050001	BKU60001	
Replacement Middle Gear Complete 1:1	C1038022	BKU68022N	

Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Z25L

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

C1038015

BKU68015



Z45L

Position / Description / Application


Manufacturers Part

MBA Part

Head Gear Complete 4:1.2

C1064015


BKU64015Z



Replacement Middle Gear Complete 4:1.2

C1064022

BKU64022Z



Z800

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Press Button

P1051050 - TiZ-MU04

BKU61050Z



Z800 Back

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Rotor Cartridge

SG1590X-D

BKUSG1590XCT

High Speed Ceramic-Hybrid Bearing - Torton Retainer

Special Size, ceramic-hybrid, Torton retainer - Checking Details



Z800 Front

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Rotor Cartridge

SG1569X-D

BKUSG1569XCT

High Speed Ceramic-Hybrid Bearing - Torton Retainer

Special Size, ceramic-hybrid, Torton retainer - Checking Details



Z900

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Press Button

P1052050 - TiZ-SU04

BKU62050



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Z900 Back

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Rotor Cartridge

SG1590X-D

BKUSG1590XCT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Special Size, ceramic-hybrid, Torlon retainer - Checking Details



Z900 Front

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Rotor Cartridge

SG1569X-D

BKUSG1569XCT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Special Size, ceramic-hybrid, Torlon retainer - Checking Details



Z95 Front

Position / Description / Application

Manufacturers Part

MBA Part


Bearing for Headgear

SG1614-D

BKUSG1614XCT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

6, ceramic-hybrid, Torlon retainer - to check



Z95L

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete

C1034017

BKU64017

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Replacement Middle Gear Complete 1:5

C1034022

BKU64022N

Drive Gear Assembly - Dog Clutch Complete 1:5

C1034013

BKU64013

Complete Assembly



Head Gear Complete

C1034015

BKU64015

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

Accutorg

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

SR144A-D27-ZS-T9H-OPT



Borden

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-CSN-KC-OPT



Replacement Bearing

High Speed Bearing - Phenolic Retainer

Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)

2 Required

SR144-D28-KC-OPT



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

ACCESS PRO NL

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKUR3100003	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2009914	BKUR3100005	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2009914	BKUR310000	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKUR3100004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKURC310000	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

KL Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKUR3100004	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKUR3100003	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2009914	BKUR310000	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2009914	BKUR3100005	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009914	BKURC310000	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

ROTAMAX PRO KL Torque

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKUR3100014	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKUR3100013	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2009915	BKUR310001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2009915	BKUR3100015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKURC310001	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

ROTAMAX PRO NL

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	2009915	BKUR3100015	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKUR3100014	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	2009915	BKUR310001	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKUR3100013	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	2009915	BKURC310001	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

102L Strong

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61102



103

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61103



106

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61106



204 Strong 90

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Armature

BKU61102



Forte 100

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

BKU68882

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 100 II

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU61073

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 100a III

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

BKU68882

True running accuracy less than 0.03, bore diameter 2.35mm



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Forte 102

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU68883

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 103L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU61073

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 106

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU68882

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 106ZO

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU68882

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 300 II

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU68883

True running accuracy less than 0.03, bore diameter 2.35mm



Forte 300 III

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Chuck

107-3

BKU68883

True running accuracy less than 0.03, bore diameter 2.35mm



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities



** Available while stocks last

Not in current production

Lab Drive

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	830ZZ	S693A-ZZS-MC3-T9H	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x10x4)	1040ZZ	SMR104A-ZZ-MC3-JP6	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (8x14x4)	1480ZZ	SMR148A-ZZS-MC3-GR3-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with high speed polyamide cage (5x10x4)	1050ZZ	SMR105A-ZZS-MC3-GR3-T9H	
Bearing for Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide Cage (7x13x4)	1370ZZ	SMR137A-ZZS-MC3-GR3-T9H	
Bearing for Motor Shaft High Speed Bearing - Polyamide Retainer Shielded Option Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (4x9x4)	940ZZ	S684A-ZZS-MC3-GR3-T9H	

Micro-Motor

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Shielded Option Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)	SAESHIN003	S693A-ZZS-MC3-T9H	
Replacement Bearing 4x9x4 Stainless Steel Bearing Shielded Option Double shielded P6 precision stainless steel bearing (4x9x4)	SAESHIN004	S684A-ZZ-MC3-JP6	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

Strong 102


Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60107



Strong 105

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60107



Strong 105L


Position / Description / Application

Manufacturers Part

MBA Part

Replacement Motor Brush

BKU60107



ZO

Position / Description / Application

Manufacturers Part

MBA Part





Replacement Chuck

BKU68882





True running accuracy less than 0.03, bore diameter 2.35mm



C1

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	
Replacement Felt Washer Laboratory Drive	4129	BKU64129	
Bearing for Motor Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)	4120	S688B-ZZS-MC3-GR3	




C2

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	
Replacement Felt Washer Laboratory Drive	4129	BKU64129	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)	4120	S688B-ZZS-MC3-GR3	

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production



C3

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Precision Bearing Special design for Schick dental handpiece replacement bearings part number 7120 (5x9.52x3.75)	7212	SWMR952505A-ZZS-MC3-GR7	
Bearing for Shaft 8x15.875x6 High Speed Bearing - Phenolic Retainer Double shielded precision stainless steel bearing with high speed polyamide retainer (8x15.88x6)	7111	SMR15875-08-ZZ-MC3-KC	
Bearing for Shaft High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with high speed phenolic cage (10.15x19x7)	7120	SMR19-10.15-ZZS-MC3-GR3-KC	

C3 7050CN

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft High Precision Bearing Double shielded Abec 3 precision stainless steel bearing (8x16x6)	4120	S688B-ZZS-MC3-GR3	




CN

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	



Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Q Profi

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Motor Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4x9x4)	5212	S684A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (7x14x5)	611/1	S687A-ZZ-MC3-JP6	
Bearing for Shaft High Precision Bearing Shielded Option Double shielded Abec 3 precision stainless steel bearing (8x16x6)	4120	S688B-ZZS-MC3-GR3	 Image Coming Soon




Qube

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Shaft 5x11x5 High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (5x11x5)	9160	S685A-ZZS-MC3-GR3-KC	 Image Coming Soon
Bearing for Shaft High Speed Bearing - Phenolic Retainer Double shielded Abec 3 precision stainless steel bearing with a high speed phenolic cage (6x13x5)	9162	S686B-ZZS-MC3-GR3-KC	 Image Coming Soon


Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

001

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6710441	BKUR2200343	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6710425	BKUR2200313	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6710433	BKUR2200283	

3000

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required		SR144-D28-ZS-KC-OPT	

Key Usually stocked or Short Delivery ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

4000 MS-ML

Position / Description / Application

Manufacturers Part

MBA Part

<p>Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required</p>	<p>SR144-D106-ZS</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required</p>	<p>SR144-D106-CSN-ZS-TO</p>	
<p>Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required</p>	<p>SR144-D106-ZS-TO</p>	
<p>Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required</p>	<p>SR144-D106-ZS-KC</p>	
<p>Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required</p>	<p>SR144-D106-CSN-ZS-KC</p>	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

4000 S-SL Push Button

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	 Image Coming Soon
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	 Image Coming Soon
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	 Image Coming Soon
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	 Image Coming Soon





P Control

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6084979	BKUR2200223	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	6084979	BKUR2200224	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	6084979	BKUR2200225	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly excluding bearings	6084979	BKUR220022	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

P Racer

Position / Description / Application	Manufacturers Part	MBA Part	
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly excluding bearings	6084987	BKUR220023	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6084987	BKUR2200233	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	6084987	BKUR2200234	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck - Use only with Oil Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	6084987	BKUR2200235	

Siroboost

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	6334598	BKUR2200244	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	6334598	BKUR220024	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	6334598	BKUR2200243	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	6334598	BKUR2200245	

Key



Usually stocked or Short Delivery

** Available while stocks last





Long Delivery or Large Minimum Quantities

Not in current production

T1 Control

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	4175217	BKUR2200204	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	4175217	BKUR2200205	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	4175217	BKUR220020	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	4175217	BKUR2200203	

T1 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	4178591	BKUR2200183	
Replacement Turbine With Bearings Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	4178591	BKUR2200185	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	4178591	BKUR2200184	
Replacement Turbine Excluding Bearings Turbine sub assembly excluding bearings	4178591	BKUR220018	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

T1 Mini TMC3LK

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	 Image Coming Soon
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	 Image Coming Soon
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	 Image Coming Soon
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	 Image Coming Soon

T2 Control

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	5868299	BKUR220022	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	5868299	BKUR2200223	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	5868299	BKUR2200224	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	5868299	BKUR2200225	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

T2 Control Serial Numbers 602 001 to 802 000 with

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

6432590

BKUR2200303

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



T2 Control Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

6710425

BKUR2200313

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Mini

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

5957142

BKUR2200213

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Replacement Turbine With Bearings

5957142

BKUR2200214

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

5957142

BKUR2200215

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



Replacement Turbine Excluding Bearings

5957142

BKUR220021

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings



T2 Mini Serial Numbers 602 001 to 702 000 with Wed

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6432632

BKUR2200323

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

T2 Mini Serial Numbers 702 001 to 802 000 with Wed

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6548130

BKUR2200333

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Mini Serial Numbers 802

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

6710441

BKUR2200343

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Mini Serial Numbers 802 001

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710441

BKUR2200343

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Racer

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

5868307

BKUR2200234

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

5868307

BKUR220023

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings



Replacement Turbine With Bearings

5868307

BKUR2200235

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

5868307

BKUR2200233

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

**** Available while stocks last**

Not in current production

T2 Serial Numbers 602 001 to 702 000

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

6432616

BKUR2200263

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

6548122

BKUR2200273

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Serial Numbers 802

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710433

BKUR2200283

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T2 Serial Numbers 802 001

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710433

BKUR2200283

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Boost Serial Numbers 602 001 to 702 000 with We

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6548122

BKUR2200273

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

6432616

BKUR2200263

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Boost Serial Numbers 802

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710433

BKUR2200283

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

T3 Boost Serial Numbers 802 001

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Turbine With Bearings

6710433

BKUR2200283

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Mini

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine Excluding Bearings

5957142

BKUR220021

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

5957142

BKUR2200213

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

5957142

BKUR2200214

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



Replacement Turbine With Bearings

5957142

BKUR2200215

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



T3 Mini Serial Numbers 602 001 to 702 000 with Wed

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6432632

BKUR2200323

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Mini Serial Numbers 702 001 to 802 000 with Wed

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6548130

BKUR2200333

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

T3 Mini Serial Numbers 802

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710441

BKUR2200343

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Mini Serial Numbers 802 001

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

6710441

BKUR2200343

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



T3 Racer

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine Excluding Bearings

5868307

BKUR220023

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

5868307

BKUR2200235

For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

5868307

BKUR2200234

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

5868307

BKUR2200233





For Handpiece with Wedge Type Chuck
Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

T3 Racer Serial Numbers 600 000

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	6486240	BKUR220025	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	6486240	BKUR2200253	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer	6486240	BKUR2200254	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	6486240	BKUR2200255	

TC3

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	4175217	BKUR2200203	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	4175217	BKUR220020	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	4175217	BKUR2200204	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	4175217	BKUR2200205	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TM1

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Deleted Item. Available while stocks last. Narrow stepped OD precision dental handpiece bearing with phenolic 2 Required		SR144-D106-ZS	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.381, ceramic-hybrid, Torlon retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Narrow stepped OD precision dental handpiece bearing with high speed Torlon retainer 2 Required		SR144-D106-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.381, ceramic-hybrid, Phenolic retainer - to check (3.18x6.35x2.38) 2 Required		SR144-D106-CSN-ZS-KC	
Replacement Bearing High Speed Bearing - Phenolic Retainer Narrow stepped OD precision dental handpiece bearing with high speed phenolic retainer 2 Required		SR144-D106-ZS-KC	

TMC

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	4178591	BKUR2200184	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	4178591	BKUR2200183	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	4178591	BKUR220018	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	4178591	BKUR2200185	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TS1

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-KC-OPT	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	5653054	BKUR2200174	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	5653054	BKUR2200175	
Replacement Turbine With Bearings Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	5953054	BKUR2200173	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	5953054	BKUR220017	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TS2

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	5953054	BKUR220017	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	5653054	BKUR2200174	
Replacement Bearing High Speed Bearing - Phenolic Retainer Stepped OD bearing with high speed phenolic retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-KC-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Stepped OD ceramic-hybrid bearing with high speed phenolic retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-KC-OPT	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	5953054	BKUR2200173	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	5653054	BKUR2200175	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0125	SR144A-D105-CSN-ZS-TO	
Replacement Bearing High Speed Bearing - Torlon Retainer Stepped OD bearing with Torlon retainer for precision dental handpieces 2 Required		SR144A-D105-ZS-TO-OPT	
Replacement Bearing High Speed Ceramic-Hybrid Bearing - Torlon Retainer Stepped OD angular contact ceramic-hybrid bearing with Torlon retainer for precision dental handpieces 2 Required	220-0793K	SR144A-D107-CSN-ZS-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

000

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

5828889

BKU68889



001

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)



001X

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

6367648

BKU67648



13Z

Position / Description / Application

Manufacturers Part

MBA Part


Head Gear Complete

6175587 / 6405687

BKU65587

For Handpiece with Wedge Type Chuck

Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece




Head Gear Complete - Standard

6175587 / 6405687

BKU65587X

For Handpiece with Wedge Type Chuck

Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece



13Z Mini

Position / Description / Application

Manufacturers Part

MBA Part


Head Gear Complete

6414457 / 6405679

BKU65687

For Handpiece with Wedge Type Chuck

Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece



13Z Mini Standard

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Standard

6414457 / 6405679

BKU65687X

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

40

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)



6

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)



6:1 VDW

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

5887174

SMR74D-MC3-GR5-L23

High Precision Bearing

Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)



A170

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete

6414457 / 6405679

BKU65687M

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece



A200

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:5

6175603


BKU65603WB



Replacement Middle Gear Complete 1:5

6175603

BKU65603



A200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

6049337

BKU69337



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

A40L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

6000959

BKU60959



Replacement Head Gear 1:1

6000983

BKU60983



C 1.6

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

4181157

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



C 1.6 Serial Number Over 65 000

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360577

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



C 200 40 6 1.6

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)



C+S 200

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)



C1.6

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 6:1

6405703

BKU65703




Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

C200

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete Diameter 9.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece	6175587 / 6405687	BKU65587	
Replacement Middle Gear Complete 1:5	6175603	BKU65603	
Replacement Middle Gear Complete 1:5	6175603	BKU65603WB	
Head Gear Complete Diameter 8.7mm premium head gear assembled with ceramic ball bearings tested in master electric handpiece	6414457 / 6405679	BKU65687	
Head Gear Complete Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece	6414457 / 6405679	BKU65687X	
Replacement O-Ring Special Dental Applications O-Ring for Electric handpiece Special dental applications O-Ring	6328236	BKU68236	
Head Gear Complete Diameter 9.7mm assembled with ceramic ball bearings and tested in master electric handpiece	6175587 / 6405687	BKU65587X	

C200 T1

Position / Description / Application	Manufacturers Part	MBA Part	
Head Beaing High Speed Bearing - Torlon Retainer Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2)		SMR7-032-MC3-GR7-TO-AF12	 Image Coming Soon

Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

C200L 8.7 Diameter

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Premium

6414457 / 6405679

BKU65687M

For Handpiece with Wedge Type Chuck

Diameter 8.7mm assembled with ceramic ball bearings and tested in master electric handpiece



C200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

6049337

BKU69337



C40

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Rear Bearing

Stainless Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)

S683-MC3-JP6-AF12



Replacement Press Button

4 Parts

BKU65083




Replacement O-Ring

Special Dental Application O-Ring

Special dental applications O-Ring

6328236

BKU68236



C40 Mini

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

6049329




BKU69329





Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


C40L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Head Gear 1:1	6000983	BKU60983	
Replacement Head Gear 1:1 Mini	6405695	BKU65695	
Replacement Middle Gear Complete 1:1	6000959	BKU60959	


C40L mini

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	6360593	SMR63-MC3-JP6-AF12	
Bearing for Headgear High Precision Bearing Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)	5887174	SMR74D-MC3-GR5-L23	

C6

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 6:1	6405703	BKU65703	



C6L mini

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)	6360593	SMR63-MC3-JP6-AF12	




Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


EL1

Position / Description / Application	Manufacturers Part	MBA Part	
Xenon Lamp	5945407	BKU65407Xenon	
LED Bulb	5945407	BKU60291LED	



Endo

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 6:1	6000777	BKU60777	
Replacement Press Button	5997528	BKU67528	
Replacement Middle Gear Complete 6:1	5994038	BKU64038	

ENDO 6:1 VDW

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	6000777 + 5994038	BKU67738	

Endo X

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 6:1	H1033C1051015	BKU61015X	
Replacement Middle Gear Complete 6:1	H1033C1051022	BKU61022X	

Key **Usually stocked or Short Delivery**
**** Available while stocks last**

Long Delivery or Large Minimum Quantities
Not in current production

H40 6

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

4181132

S693-ZZS-MC34-TH5-R-L23

High Precision Bearing

Open Abec 5 precision stainless steel bearing with low-torque retainer and light oil lubricant (3x8x3)



H40L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

4181199

S623-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (3x10x4)



Bearing for Shaft

4181108

SMR85-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



Quick Coupling

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb and Cap

6314558

BKU64558LED



Quick Coupling Sirosonic L

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

5940291


BKU60291LED



Xenon Lamp

5940291

BKU60291Xenon



S 12000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Chuck End

SF692-ZZS-MC34-GR5-R-L23

Stainless Flanged High Precision Bearing

Double shielded flanged Abec 5 precision stainless steel bearing with NON-ISO Width and low-torque ribbon retainer (2x6x2.3)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

S 200 40 6 1.6

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Drive

6360601

SMR1035-06-MC3-GR7-TO

High Speed Bearing - Torlon Retainer

Open stainless steel bearing with high speed Torlon retainer. Special size for dental applications. (6x10.35x3)



Siroboost

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

6334598

BKUR220024T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



SL

Position / Description / Application


Manufacturers Part

MBA Part

LED Bulb

5945407

BKU60291LED



Bearing for Rotor

7040801

S684A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x9x4)



Xenon Lamp

5945407

BKU65407Xenon



Bearing for Rotor

7036130

S692-ZZS-MC34-GR5-R-G48

High Precision Bearing Shielded Option

Double shielded Abec 5 precision bearing with Low Torque ribbon retainer and grease lubrication. Non ISO Width. (2x6x2.3)



Smart Plus

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 6:1

H1033C1051015

BKU61015X



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Smart Plus 6:1

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 6:1

H1033C1051022

BKU61022X



T1 Boost Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)



T1 Classic S200

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

5890855

SMR7-032-MC3-GR7-TO-AF12

High Speed Bearing - Torlon Retainer

Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2)



Bearing for Headgear

5890855

SMR7-032-CSN-MC3-GR7-TO-AF12

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (3.18x7x2)



T1 Classic S200 mini

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

5899863

SMR6-3.175-CSN-ZZS-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6x2)



Bearing for Headgear

5899863

SMR6-3.175-CSN-ZZS-MC3-J3-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

3.175x6x2, ceramic-hybrid, Torlon retainer - to check (3.18x6x2)



T1 CLASSIC S200L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6175587 + 6175603

BKU68703

Dental replacement drive kit



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

T1 CLASSIC S200L Mini

Position / Description / Application


Manufacturers Part

MBA Part

Drive Kit

6414457 + 6175603

BKU65703S



Drive Kit

6405679 + 6175603

BKU67903



T1 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part


MBA Part

Replacement O-Ring

7028988

BKU68988

Special Dental Application O-Ring



T1 Control Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

6367648

BKU67648



T1 Control Serial Number Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)



T1 Control till Serial Number 65

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

5828889

BKU68889



T1 Control till Serial Number 65 000

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

5828889

BKU68889



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

T1 Line C1.6L

Position / Description / Application


Manufacturers Part

MBA Part

Drive Kit

6405703 + 6000942

BKU60342



T1 Line C200

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

4 Parts

BKU66967



Bearing for Headgear

5890855

SMR7-032-MC3-GR7-TO-AF12

High Speed Bearing - Torlon Retainer

Open Abec 7 precision stainless steel bearing with high speed Torlon retainer (3.18x7x2)



Bearing for Headgear

5890855

SMR7-032-CSN-MC3-GR7-TO-AF12

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (3.18x7x2)



T1 Line C200 mini

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

5899863

SMR6-3.175-CSN-ZZS-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6x2)



Bearing for Headgear

5899863

SMR6-3.175-CSN-ZZS-MC3-J3-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option

3.175x6x2, ceramic-hybrid, Torlon retainer - to check (3.18x6x2)



T1 Line C200L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6175587 + 6175603

BKU68703

Dental replacement drive kit



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

T1 Line C200L Mini

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6405679 + 6175603

BKU67903



T1 Line C40L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6000983 + 6000959

BKU68359

Dental replacement drive kit



T1 Line C40L Mini

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6405695 + 6000959

BKU69559

Dental replacement drive kit



T1 Line C6L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6405703 + 6000942

BKU60342



T2

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

5868307

BKUR220023T

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



T2 Boost Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)



T2 Boost Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

6367648

BKU67648



Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

T2 Control

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

5868299

BKUR220022T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



Bearing for Rotor Cartridge

5994814

BKU24814

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

3.175x6.35x2.778, ceramic-hybrid, Torlon retainer - Checking Details (3.18x6.35x2.78)



T2 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

5829416

BKU69416

Special Dental Application O-Ring



T2 Control Serial No. Over 602 001

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

6554815

SWR144B-CSN-ZE-MC3-GR7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Double shielded Abec 7 precision ceramic-hybrid bearing with high speed Torlon retainer (3.18x6.35x2.35)



T2 Control Serial Number 602 001X

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

6367648

BKU67648



T2 Line A200L

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6562586 + 6175603

BKU68603



T2 Line A40 Mini

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

6049352 + 6000959

BKU65259



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

T2 Line A40L

Position / Description / Application

Manufacturers Part


MBA Part

Drive Kit

6000983 + 6000959

BKU68359

Dental replacement drive kit



T2 Racer

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

5828897

BKU68897



Bearing for Rotor Cartridge

5994814

BKU24814

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

3.175x6.35x2.778, ceramic-hybrid, Torlon retainer - Checking Details (3.18x6.35x2.78)



T2 Racer Air Driven Handpiece

Position / Description / Application

Manufacturers Part


MBA Part

Replacement O-Ring

5829416

BKU69416

Special Dental Application O-Ring



T2 Revo

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

4 Parts

BKU65083



T2 Revo Mini

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

6049329

BKU69329



T2 Revo R170

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

7044654

SMR55-0235-T4HB-MC3-AF12

High Speed Bearing - Torlon Retainer





Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)




Key **Usually stocked or Short Delivery**
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


T2 Revo R40

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	5943092 + 5943167	BKU69267	
Replacement Head Gear 1:1 Mini	6049352	BKU69352	
Replacement Head Gear 1:1	5943092	BKU63092	
Replacement Middle Gear Complete 1:1	5943167	BKU63167	


T2 Revo R40 mini

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	6405695 + 5943167	BKU69567	

T2 S

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	6562586 + 6175603	BKU68603	

T3 Line E40 Mini

Position / Description / Application	Manufacturers Part	MBA Part	
Drive Kit	6049352 + 6000959	BKU65259	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

T3 Racer

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

5868307

BKUR220023T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



TEH40

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

7025273

SMR105-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x10x3)



TS1 Control Air Driven Handpiece

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

7028988

BKU68988

Special Dental Application O-Ring



VDW

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 6:1

6000777

BKU60777



Replacement Middle Gear Complete 6:1

5994038

BKU64038



VDW 6:1

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Press Button

5997528

BKU67528



X SMART PLUS 6:1

Position / Description / Application

Manufacturers Part

MBA Part

Drive Kit

H1033C1051015 + H1033C1051022

BKU615022



Key

Usually stocked or Short Delivery





Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

SS White®



SS White

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Cap End High Speed Bearing - Polyamide Retainer Precision flanged dental bearing with a high speed polyamide retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-T9H-OPT	
Bearing for Cap End High Speed Bearing - Phenolic Retainer Single shielded flanged precision dental handpiece bearing with phenolic retainer (3.18x6.35x2.78)		SFR144A-D27-ZSO-KC-OPT	
Bearing for Chuck End High Speed Bearing - Phenolic Retainer Narrow precision dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-KC-OPT	 Image Coming Soon
Bearing for Chuck End High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.38)		SR144-D28-CSN-KC-OPT	 Image Coming Soon





Key **Usually stocked or Short Delivery**
 ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


430

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic LUBEFREE	263739	BKURMU290108	 <p style="font-size: 8px; color: red;">Ultra precision sub-assemblies, cartridges</p>
Turbine Sub Assembly Without Bearings For Handpiece with Wedge Type Chuck	263739	BKURM290108	 <p style="font-size: 8px; color: red;">Ultra precision sub-assemblies, cartridges</p>

430 Vista

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	263968	BKUR2901094	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Phenolic retainer - myonic LUBEFREE	263968	BKUR2901083	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	263968	BKUR2910945	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	263968	BKUR290109	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>

430 Vista Lubefree

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)	SG1303X-D	SFR144A-D1303-CSN-ZSO-KC	 <p style="font-weight: bold; color: black;">Image Coming Soon</p>

Key Usually stocked or Short Delivery ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

430 Vista Manual Chuck

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

418WIFU

SFR144A-D37-ZSO-KC

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)



Bearing for Rotor Cartridge

418WIFU-CER

SFR144A-D37-CSN-ZSO-KC

Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)



430SWL

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Washer - Thick

Precision Washer 6.5 x 8.30 x 0.5 Stainless
(6x8x0.5)

W0060-FP-008-0050-S4



Replacement Washer - Thin

Precision Washer 6.5 x 8.36 x 0.2 Stainless
(6x8x0.2)

W0060-FP-008-0020-S4



Replacement O-Ring

Special Dental Application O-Ring
Special sized o-ring used in dental Applications
2 Required

OR-00600-080-N70



Replacement Wave Washer

6.4mm Shaft Size
(6.35x9.32x0.76)

W0064-T-010-008-W3-S2



C6

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

262522

BKU62522Xenon



COLORight System

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

260986

BKU60986Xenon



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Concentrix

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Futura I and II

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Solara

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

263968

BKUR2910945

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic




Replacement Turbine With Bearings

263968

BKUR2901094

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic




Replacement Turbine With Bearings

263968

BKUR2901083

For Handpiece with Wedge Type Chuck

Turbine sub assembly with singular contact ceramic-hybrid bearings with Phenolic retainer - myonic LUBEFREE




Replacement Turbine Excluding Bearings

263968

BKUR290109

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings



Solara Lubefree

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1303X-D

SFR144A-D1303-CSN-ZSO-KC

Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Star 430 Air Driven Handpiece

Position / Description / Application

Manufacturers Part


MBA Part

Replacement O-Ring

O-Ring+Washer+Disk

BKU69997

Special Dental Application O-Ring



StarBright ISO

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

262522

BKU62522Xenon



Starflex Lubefree

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

SG1303X-D

SFR144A-D1303-CSN-ZSO-KC

Lubefree - High Speed Ceramic-Hybrid Bearing - LUBEFREE Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with LUBEFREE retainer (3.18x6.35x2.78)



Ultra

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SG1303X-D

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required



Vista PB

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

418WIFU

SFR144A-D37-ZSO-KC

Daily Lube - High Speed Bearing - Phenolic Retainer

Single shielded flanged stainless steel precision dental bearing with phenolic retainer (3.18x6.35x2.78)



Bearing for Rotor Cartridge

418WIFU-CER

SFR144A-D37-CSN-ZSO-KC

Daily Lube - High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Single shielded flanged stainless steel ceramic-hybrid precision dental bearing with phenolic retainer (3.18x6.35x2.78)



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

Starflite

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR2-5A-D17-ZS-TO

High Speed Bearing - Torlon Retainer

Single shielded Abec 7 precision special dental bearing with high speed Torlon retainer (3.18x7.94x3.57)

2 Required



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

EVO M5 L

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete

97230612

BKU60612

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Head Gear Complete

97230612

BKU60612M

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

Miera

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Thera

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Titan

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Turbina Master

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR2-5A-D17-ZS-KC

High Speed Bearing - Phenolic Retainer

Single shielded Abec 7 precision special dental bearing with high speed phenolic retainer (3.18x7.94x3.57)

2 Required



Turbina Mini Master

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

X 6000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SMR93A-ZZ-MC3-JP6

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3x9x4)

2 Required



Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

115K 14000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

00542400

S625-ZZ-MC3-JP6

Stainless Bearing
Double shielded P6 precision stainless steel bearing (5x16x5)
2 Required



260 Blue and Black

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 7-14-5mm

03202200

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer
Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)



Bearing for Shaft 8-16-6mm

03202300

S688B-ZZ-MC3-GR3-T9H-G36

High Speed Bearing - Polyamide Retainer
Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



300 Series

Position / Description / Application

Manufacturers Part


MBA Part

Turbine Sub Assembly With Bearings

07234102

BKUR2300502

For Handpiece with Wedge Type Chuck
With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic




Replacement Turbine With Bearings

07234102

BKUR2300504

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic




Replacement Turbine Excluding Bearings

07234102

BKUR230050

For Handpiece with Wedge Type Chuck
Turbine sub assembly excluding bearings



400 Series

Position / Description / Application

Manufacturers Part


MBA Part

Replacement Turbine With Bearings

07005900

BKUR230039T

For Handpiece with Wedge Type Chuck
Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic



Key



Usually stocked or Short Delivery

** Available while stocks last





Long Delivery or Large Minimum Quantities

Not in current production

56A

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	BB 5-8-2	SMR85-MC3-JP6-AF12	
Replacement Bearing Stainless Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (4x7x2)	BB 4-7-2	SMR74-MC3-JP6-AF12	

796

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	02765500	BKUR2300255	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	02765500	BKUR2300253	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	02765500	BKUR230025	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02765500	BKUR2300254	 <p style="font-size: 8px; color: red;">Precision sub-assemblies, cartridges with steel inserts & WCC coated actuators</p>





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

798

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	02760000	BKUR230026	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	02760000	BKUR2300265	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR2300263	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR2300264	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR230026T	

896

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	02765500	BKUR2300255	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02765500	BKUR2300254	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	02765500	BKUR2300253	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	02765500	BKUR230025	

Key



Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


898

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	02760000	BKUR2300265	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR230026T	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR2300264	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	02760000	BKUR230026	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	02760000	BKUR2300263	

969

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)	01579600	SMR55-0235-T4HB-MC3-AF12	

99

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	01565800	SMR85-MC3-JP6-AF12	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

999

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01589400

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



999E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

04486000

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



A20

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

01850600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



A25

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

01850600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



ADEC

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01565800

SMR85-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



Airmotor 13000 RPM

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

01850600

SR2-5A-ZZ-MC3-JP6

Stainless Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)

2 Required



AM 25

Position / Description / Application


Manufacturers Part

MBA Part

Xenon Lamp

04300400

BKU60400Xenon



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Contra Angle Handpiece New

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

02649400

SMR74-ZS-T4HB-MC3-GR5-AF12

High Precision Bearing

Single shielded Abec 5 precision stainless steel bearing with high speed Torlon retainer and light oil lubricant. NON-ISO width (4x7x2)

2 Required



Contra Angle Handpiece Old

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

01579300

SMR55-0235-MC3-TO-AF12

High Speed Bearing - Torlon Retainer

Open precision dental handpiece bearing with Torlon retainer (2.35x5.5x2)

2 Required



EB40

Position / Description / Application

Manufacturers Part


MBA Part

Dental Motor Set with 21 parts

1502217-001

BKU62217


Includes parts - 3x 01862300 - 1x 04062500 - 1x 04061800 - 1x 04062300 - 1x 04058800 - 2x 04062100 - 2x 04062200 - 2x 04062400 - 2x 04058800 - 3x 04059200 - 2x 04059300 - 1x 04069600



Replacement Motor Brush

04059300

BKU69300



EB40A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

04059000

SR1-5A-ZZS-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



Bearing for Rotor

04062000

SMR104A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (4x10x4)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

EB40E

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

04062000

SMR104A-ZZ-MC3-JP6

High Precision Bearing Shielded Option
Double shielded P6 precision stainless steel bearing (4x10x4)



Bearing for Rotor

04059000

SR1-5A-ZZS-MC3-JP6

High Precision Bearing Shielded Option
Double shielded P6 precision stainless steel bearing (2.38x7.94x3.57)



HA43A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01672300

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option
Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)



Bearing for Lower Shaft

01672300

SMR85C-ZZS-MC34-GR5-L23

High Precision Bearing Standard Cage Option
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



HD43

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01381900

SMR85-MC3-JP6-AF12

High Precision Bearing
Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)



HD43EM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01672300

SMR85C-ZZS-MC34-GR5-TR-L23

High Precision Bearing Low Torque Option
Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)



Bearing for Lower Shaft

01672300

SMR85C-ZZS-MC34-GR5-L23

High Precision Bearing Standard Cage Option
Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

LA66

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 7-14-5mm

03202200

S687A-ZZS-MC3-GR3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (7x14x5)



Bearing for Shaft 8-16-6mm

03202300

S688B-ZZ-MC3-GR3-T9H-G36

High Speed Bearing - Polyamide Retainer

Double shielded Abec 3 precision stainless steel bearing with a high speed polyamide cage (8x16x6)



LT25

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor

01850600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



Motor MF Perfecta

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft 3-8-4mm

01206200

S693A-ZZS-MC3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)



Bearing for Shaft

02685500

S685A-ZZS-MC34-TH7-R-G48

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)



Motor MF Tectorque

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Shaft

02685500

S685A-ZZS-MC34-TH7-R-G48

High Precision Bearing

Double shielded Abec 7 precision stainless steel bearing with low-torque ribbon retainer (5x11x5)



Bearing for Shaft 3-8-4mm

01206200

S693A-ZZS-MC3-T9H

High Speed Bearing - Polyamide Retainer

Double shielded high speed stainless steel bearing with Polyamide retainer (3x8x4)



Key


Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities




** Available while stocks last

Not in current production




RA 24

Position / Description / Application	Manufacturers Part	MBA Part	
Xenon Lamp	04208000	BKU68000Xenon	


RC95BC

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05421700	BKUR230046	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic	05421700	BKUR2300462	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	05421700	BKUR2300464	

RC95RM

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05421700	BKUR230046	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	05421700	BKUR2300464	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic	05421700	BKUR2300462	


RQ24

Position / Description / Application	Manufacturers Part	MBA Part	
Roto Quick Coupler O-Ring Special dental applications O-Ring	02060200+02060100	BKU60200	




Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


RQ34

Position / Description / Application	Manufacturers Part	MBA Part	
Roto Quick Coupler O-Ring	02060200+02060100	BKU60200	
Special dental applications O-Ring			




S10

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft	01577000	S623-ZZ-MC3-JP6	
High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)			
Bearing for Lower Shaft	01672300	SMR85C-ZZS-MC34-GR5-TR-L23	
High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)			
Bearing for Lower Shaft	01672300	SMR85C-ZZS-MC34-GR5-L23	
High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)			

S11

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft	01577000	S623-ZZ-MC3-JP6	
High Precision Bearing Double shielded P6 precision stainless steel bearing (3x10x4)			



S9

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft	01577000	S623-ZZ-MC3-JP6	
High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)			
Bearing for Lower Shaft	01672300	SMR85C-ZZS-MC34-GR5-TR-L23	
High Precision Bearing Low Torque Option Double shielded Abec 5 precision bearing with low-torque ribbon retainer (5x8x3)			
Bearing for Lower Shaft	01672300	SMR85C-ZZS-MC34-GR5-L23	
High Precision Bearing Standard Cage Option Double shielded Abec 5 precision bearing with standard crown retainer (5x8x3)			


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production




Straight Handpiece 6000 RPM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Centre Stainless Bearing Double shielded P6 precision stainless steel bearing (5x11x4)	02685500 ?	SMR115-ZZ-MC3-JP6	
Bearing for Chuck End Stainless Bearing Open P6 precision bearing with low-torque retainer and light oil lubricant (3.18x6.35x2.38)	00352800	SR144-MC3-JP6-AF12	

Synea

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	02749408	BKU69408	

Synea TA98

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76) 2 Required		W0048-T-007-008-W3-S2	
Replacement O-Ring Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00605-079-VBR75	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required		SR144A-D27-ZS-T9H-OPT	

Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Synea TA98L

Position / Description / Application

Manufacturers Part

MBA Part

Replacement O-Ring

Very High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required

OR-00605-079-VBR75



Replacement Bearing

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required

SR144A-D27-ZS-T9H-OPT



Replacement Wave Washer

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)

2 Required

W0048-T-007-008-W3-S2



Synea Vision

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Spring Washer

06876000

BKU66000



TA96

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine Excluding Bearings

04396700

BKUR230029

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

04396700

BKUR2300295

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

04396700

BKUR2300294

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

04396700

BKUR2300293

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

TA96CM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

04378600

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)



TA96L

Position / Description / Application

Manufacturers Part

MBA Part

Xenon Lamp

02207600

BKU67600Xenon

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)



Bearing for Rotor Cartridge

02770712

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)



TA96LC

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

04378600

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)



TA96LCM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

04378600

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)



Key





Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


** Available while stocks last

Not in current production


TA97C LED

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05307200	BKUR2300305	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05307200	BKUR2300304	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05307200	BKUR2300303	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05307200	BKUR230030	


TA97CLED

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	05288500	SR144A-D27-ZS-T9H-OPT	

TA97CLM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	05288500	SR144A-D27-ZS-T9H-OPT	

TA97LW

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	06307800	SR144A-D27-ZS-T9H-OPT	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TA98

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	04396600	BKUR2300275	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	04396600	BKUR2300274	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	04396600	BKUR230027	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	04396600	BKUR2300273	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	04396600	BKUR230027T	

TA98C

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	04378600	SR144A-D27-ZS-T9H-OPT	

Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

TA98C LED

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05307000	BKUR2300334	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05307000	BKUR2300333	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05307000	BKUR230033	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05307000	BKUR2300335	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05307000	BKUR230033T	




TA98CLED

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	05055800	BKU65800	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TA98L

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	03503000	BKU63000	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	02770712	SR144A-D27-ZS-T9H-OPT	
Xenon Lamp	02207600	BKU67600Xenon	





TA98L New Version

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05306900	BKUR2300333	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05306900	BKUR2300334	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05306900	BKUR2300335	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05306900	BKUR230033	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TC95

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	03305000	BKUR230028	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	03305000	BKUR2300283	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	03305000	BKUR2300284	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	03305000	BKUR2300285	



Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TE95

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05958900	BKUR230036	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05958900	BKUR2300363	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05958900	BKUR2300364	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05958900	BKUR2300365	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)	05541600	SR144A-D27-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer 3.175x6.35x2.778, stainless steel-THB Phenolic - Checking Details (3.18x6.35x2.78)	SG1295X	BKUSG1295XSP	

TE95 production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)	07003500	SFR144A-D74F-CSN-ZSO-KC	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)	07003500	SFR144A-D74F-ZSO-KC	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

TE95BC

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

05541600

SR144A-D27-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)



Bearing for Rotor Cartridge

SG1295X

BKUSG1295XSP

High Speed Bearing - Phenolic Retainer

3.175x6.35x2.778, stainless steel-THB Phenolic - Checking Details (3.18x6.35x2.78)



TE97

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

05959100-05959200

BKUR2300453

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

05959100-05959200

BKUR230045

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

05959100-05959200

BKUR2300454

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

05959100-05959200

BKUR2300455

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



Key

Usually stocked or Short Delivery

** Available while stocks last



Long Delivery or Large Minimum Quantities

Not in current production

TE97 production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07234102	BKUR230050	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings, Deep Groove, Full stainless, Phenolic retainer, myonic	07234102	BKUR230052	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Phenolic retainer - myonic	07234102	BKUR230054	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Single shielded precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)	07003500	SFR144A-D74F-ZSO-KC	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Single shielded ceramic-hybrid precision special dental bearing with two flanges and high speed Phenolic retainer (3.18x6.35x2.77)	07003500	SFR144A-D74F-CSN-ZSO-KC	





TE97BC

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)	05541600	SR144A-D27-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Bearing - Phenolic Retainer 3.175x6.35x2.778, stainless steel-THB Phenolic - Checking Details (3.18x6.35x2.78)	SG1295X	BKUSG1295XSP	




Key Usually stocked or Short Delivery ** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TE97CLED

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05959100-05959200	BKUR2300453	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05959100-05959200	BKUR2300455	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05959100-05959200	BKUR2300454	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05959100-05959200	BKUR230045	

TE97CM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.78)	05958600	SR144A-D27-CSN-ZS-KC-OPT	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Torlon Retainer 3.175x6.35x2.778, ceramic-hybrid, Torlon retainer - Checking Details (3.18x6.35x2.78)	05958600	BKUA23700	
Bearing for Rotor Cartridge High Speed Ceramic-Hybrid Bearing - Phenolic Retainer 3.175x6.35x2.778, ceramic-hybrid, Phenolic retainer - Checking Details (3.18x6.35x2.78)	SG1295X-D	BKUSG1295XCP	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TE97CN

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

05958600

SR144A-D27-CSN-ZS-KC-OPT

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

Precision ceramic-hybrid dental handpiece bearing with high speed phenolic retainer (3.18x6.35x2.78)



Bearing for Rotor Cartridge

05958600

BKUA23700

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

3.175x6.35x2.778, ceramic-hybrid, Torlon retainer - Checking Details (3.18x6.35x2.78)



Bearing for Rotor Cartridge

SG1295X-D

BKUSG1295XCP

High Speed Ceramic-Hybrid Bearing - Phenolic Retainer

3.175x6.35x2.778, ceramic-hybrid, Phenolic retainer - Checking Details (3.18x6.35x2.78)



TE97M

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

05541600

SR144A-D27-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)



Bearing for Rotor Cartridge

SG1295X

BKUSG1295XSP

High Speed Bearing - Phenolic Retainer

3.175x6.35x2.778, stainless steel-THB Phenolic - Checking Details (3.18x6.35x2.78)



Key





Usually stocked or Short Delivery

** Available while stocks last




Long Delivery or Large Minimum Quantities

Not in current production

TE98

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	05958900	BKUR230036	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	05958900	BKUR2300364	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	05958900	BKUR2300365	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	05958900	BKUR2300363	

TE98 300 Series production year 2014

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings	06787501-06787502	BKUR2300494	
Turbine Sub Assembly With Bearings For Handpiece with Wedge Type Chuck With Full stainless bearings Deep Groove	06787501-06787502	BKUR2300492	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	06787501-06787502	BKUR230049	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TE98B

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

05541600

SR144A-D27-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Stainless steel bearing with high speed phenolic retainer for precision dental handpieces (3.18x6.35x2.78)



Bearing for Rotor Cartridge

SG1295X

BKUSG1295XSP

High Speed Bearing - Phenolic Retainer

3.175x6.35x2.778, stainless steel-THB Phenolic - Checking Details (3.18x6.35x2.78)



TE98CLED

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine With Bearings

04396600

BKUR230027T

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

04396600

BKUR2300275

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer




Replacement Turbine With Bearings

04396600

BKUR2300274

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine Excluding Bearings

04396600

BKUR230027

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

04396600

BKUR2300273

For Handpiece with Wedge Type Chuck





Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>


TG97 400 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07006000	BKUR230038	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	07006000	BKUR2300383	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07006000	BKUR2300384	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	07006000	BKUR2300385	

TG97L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	05958600	SR144A-D27-ZS-T9H-OPT	





TG98

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07005900	BKUR230039T	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TG98 400 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07005900	BKUR230039	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	07005900	BKUR2300393	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	07005900	BKUR2300395	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07005900	BKUR2300394	





TG98L

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	05958600	SR144A-D27-ZS-T9H-OPT	





Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TK100 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	07272900	BKUR2300513	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07272900	BKUR2300514	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07272900	BKUR230051	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	07272900	BKUR2300515	

TK94 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07273100	BKUR230052	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	07273100	BKUR2300525	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07273100	BKUR2300524	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	07273100	BKUR2300523	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

TK94L

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

07106200

SR144-D28-CSN-ZE-TO-OPT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)



TK94LM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

07106200

SR144-D28-CSN-ZE-TO-OPT

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Precision ceramic-hybrid single encapsulated shield dental special bearing with high speed Torlon retainer (3.18x6.35x2.38)



TK97 500 Series

Position / Description / Application

Manufacturers Part

MBA Part


Replacement Turbine Excluding Bearings

07272600

BKUR230043

For Handpiece with Wedge Type Chuck

Turbine sub assembly excluding bearings




Replacement Turbine With Bearings

07272600

BKUR2300433

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

07272600

BKUR2300434

For Handpiece with Wedge Type Chuck

Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic




Replacement Turbine With Bearings

07272600

BKUR2300435

For Handpiece with Wedge Type Chuck

Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer



TK97LM

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Rotor Cartridge

06307900

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer





Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)




Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


TK98 500 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer - myonic	07272800	BKUR2300443	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with ceramic-hybrid bearings with Torlon retainer - myonic	07272800	BKUR2300444	
Replacement Turbine With Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly with angular contact ceramic-hybrid bearings with Torlon retainer	07272800	BKUR2300445	
Replacement Turbine Excluding Bearings For Handpiece with Wedge Type Chuck Turbine sub assembly excluding bearings	07272800	BKUR230044	


TK98LM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Rotor Cartridge High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)	06307900	SR144A-D27-ZS-T9H-OPT	

Top Air 100 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required	02722300	SR144-D28-ZS-KC-OPT	

Top Air 200 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Phenolic Retainer Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38) 2 Required	02106302	SR144-D28-ZS-KC-OPT	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Top Air 300 Series

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

02106302

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Top Air 700 Series

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Wave Washer

W0048-T-007-008-W3-S2

4.8mm Shaft Size

Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)



PTFE Washer-Spacer

W0030-F-007-001-PT

0.13mm - half thou - thick

Deleted Item. Available while stocks last. (3x6.56x0.13)

2 Required



Replacement O-Ring

OR-00605-079-VBR75

Very High Temperature

Deleted Item. Available while stocks last. Special o-ring for dental applications

2 Required



Replacement Bearing

02722300

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)





2 Required







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

Top Air 800 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement O-Ring Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00605-079-VBR75	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required	02722300	SR144A-D27-ZS-T9H-OPT	





Top Lite 700 Series

Position / Description / Application	Manufacturers Part	MBA Part	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required	02722300	SR144A-D27-ZS-T9H-OPT	
Replacement O-Ring Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		OR-00605-079-VBR75	

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Top Lite 800 Series

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Bearing High Speed Bearing - Polyamide Retainer Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78) 2 Required	02722300	SR144A-D27-ZS-T9H-OPT	
Replacement O-Ring Very High Temperature Deleted Item. Available while stocks last. Special o-ring for dental applications 2 Required		<i>OR-00605-079-VBR75</i>	
Replacement Wave Washer 4.8mm Shaft Size Additional stock available from USA subject to Prior Sale = 4719 (4.76x6.35x0.76)		W0048-T-007-008-W3-S2	
PTFE Washer-Spacer 0.13mm - half thou - thick Deleted Item. Available while stocks last. (3x6.56x0.13) 2 Required		W0030-F-007-001-PT	

Toplight

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	02749408	BKU69408	

Trend

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Spring Washer	02749408	BKU69408	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

WA56

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	01589600	SMR74-MC34-GR7-R-L23	
Replacement Middle Gear Complete 1:1	03821000	BKU61000	
Head Gear Complete 1:1	03819000	BKU69000	
Replacement Press Button	03819500	BKU69500	
Head Gear Complete 1:1	05431100	BKU61100	
Bearing for Lower Shaft High Precision Bearing Shielded Option Double shielded P6 precision stainless steel bearing (3x10x4)	01577000	S623-ZZ-MC3-JP6	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	

WA56LT NEW

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open Abec 5 precision special dental bearing with low-torque retainer and light oil lubricant (4x7x1.6)	05039400	SMR74D-MC3-GR5-L23	

Key

Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

WA56LT OLD

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

01579300

SMR55-0235-MC3-GR3-AF12

High Precision Bearing

Open Abec 3 precision bearing (2.35x5.5x2)



WA62LTA

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

01579300

SMR55-0235-MC3-GR3-AF12

High Precision Bearing

Open Abec 3 precision bearing (2.35x5.5x2)



WA66

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



Head Gear Complete 1:1

03819000

BKU69000



WA66LTA

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

01579300

SMR55-0235-MC3-GR3-AF12

High Precision Bearing

Open Abec 3 precision bearing (2.35x5.5x2)



WA67LTA

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

01579300

SMR55-0235-MC3-GR3-AF12

High Precision Bearing

Open Abec 3 precision bearing (2.35x5.5x2)



WA86

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

WA956

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



WA966

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



WA969

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Middle Gear

01579600

SMR55-0235-T4HB-MC3-AF12

High Speed Bearing - Torlon Retainer

Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)



WA975

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



WA986

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



WA988

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

01589600

SMR74-MC34-GR7-R-L23

High Precision Bearing



Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)







Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

WA99

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	01565800	SMR85-MC3-JP6-AF12	
Bearing for Middle Gear High Speed Bearing - Torlon Retainer Open precision stainless bearing with high speed Torlon retainer (2.35x5.5x2)	01579600	SMR55-0235-T4HB-MC3-AF12	 Image Coming Soon




WA99A

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer ceramic-hybrid, Torlon retainer - to check	05798700	BKU28700	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (5x8x2)	01565800	SMR85-MC3-JP6-AF12	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	04486000	SMR74-CSN-TH7-TO	 Image Coming Soon
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)	05257200	SMR6.50-3.40-CSN-ZZS-MC3-TO	 Image Coming Soon




Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


WA99LT

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete For Handpiece with Wedge Type Chuck Assembled with ceramic ball bearings, tested in master electric handpiece	05220300	BKU60300	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer ceramic-hybrid, Torlon retainer - to check	05798700	BKU28700	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Shielded Option Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)	05257200	SMR6.50-3.40-CSN-ZZS-MC3-TO	


WD56

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 1:1	03425000	BKU65000	
Replacement Middle Gear Complete 1:1	03403500	BKU63500	
Bearing for Lower Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (3x10x4)	01577000	S623-ZZ-MC3-JP6	

WD66

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	

WD969

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

WD999

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

01589400

S683-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)



WE56

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

05060600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing Shielded Option

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



Bearing for Lower Shaft

05765800

SMR63-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



Bearing for Lower Shaft

06164000

S684-MC3-JP6-AF12

High Precision Bearing


Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)



Head Gear Complete 1:1

05843300

BKU63300



Replacement Middle Gear Complete 1:1

05852500

BKU62500W



Replacement Press Button

05843800

BKU63800



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>

WE56A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

05060600

SR2-5A-ZZ-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (3.18x7.94x3.57)



Head Gear Complete 1:1

05843300

BKU63300



WE56LEDG

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

05844400

BKU64400



WE66

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete 1:1

05843300

BKU63300



Bearing for Lower Shaft

05765800

SMR63-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



WE99A

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Headgear

02649400

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer

Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

WE99LEDG

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

05765800

SMR63-MC3-JP6-AF12

High Precision Bearing
Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



Bearing for Headgear

02649400

SMR74-CSN-TH7-TO

High Speed Ceramic-Hybrid Bearing - Torlon Retainer
Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)



WG56

Position / Description / Application


Manufacturers Part

MBA Part

Head Gear Complete 1:1

07605900


BKU65900WG



Replacement Middle Gear Complete 1:1

05233200


BKU63200



Replacement Middle Gear Complete 1:1

06796000

BKU69600



Replacement Press Button

4 Parts

BKU62800



WG99A

Position / Description / Application

Manufacturers Part

MBA Part

Head Gear Complete - Premium

07473400

BKU63400

For Handpiece with Wedge Type Chuck
Premium head gear assembled with ceramic ball bearings tested in master electric handpiece



Head Gear Complete - Standard

07473400

BKU63400M

For Handpiece with Wedge Type Chuck
Assembled with ceramic ball bearings, tested in master electric handpiece



Key

Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities

** Available while stocks last

Not in current production

WG99LT

Position / Description / Application

Manufacturers Part

MBA Part


Head Gear Complete - Premium

07473400

BKU63400

For Handpiece with Wedge Type Chuck

Premium head gear assembled with ceramic ball bearings tested in master electric handpiece



Head Gear Complete - Standard

07473400

BKU63400M

For Handpiece with Wedge Type Chuck

Assembled with ceramic ball bearings, tested in master electric handpiece



WI75

Position / Description / Application


Manufacturers Part

MBA Part

Replacement Middle Gear Complete 20:1

04667600

BKU67600



Replacement Press Button

04940500

BKU60500



Head Gear Complete 20:1

04939900

BKU69900



WI75 LEDG

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

04033500

SMR63-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (3x6x2)



Key

Usually stocked or Short Delivery

** Available while stocks last

Long Delivery or Large Minimum Quantities

Not in current production

WI75LEDG

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

06164000

S684-MC3-JP6-AF12

High Precision Bearing

Open P6 precision bearing with low-torque cage and light oil lubricant (4x9x2.5)



Bearing for Headgear

04027300

SMR74-MC34-GR7-R-L23

High Precision Bearing

Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)



Bearing for Middle Gear

04733900

SMR55-0235-MC3-GR3-AF12

High Precision Bearing

Open Abec 3 precision bearing (2.35x5.5x2)



WK56

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Middle Gear Complete 1:1

05233200


BKU63200



Replacement Middle Gear Complete 1:1

06796000

BKU69600



Head Gear Complete 1:1

05431100

BKU61100



WK93LT

Position / Description / Application

Manufacturers Part

MBA Part

Bearing for Lower Shaft

06558000

SR166-ZZS-MC3-JP6

High Precision Bearing

Double shielded P6 precision stainless steel bearing (4.76x9.52x3.18)



Key



Usually stocked or Short Delivery

Long Delivery or Large Minimum Quantities


** Available while stocks last

Not in current production





WK99LT

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (4.76x9.52x3.18)	06558000	SR166-ZZS-MC3-JP6	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open ceramic-hybrid special bearing for dental applications with high speed Torlon retainer (3.4x6.5x1.6)	05257200	SMR6.50-3.40-CSN-ZZS-MC3-TO	 Image Coming Soon

WP64

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Lower Shaft High Precision Bearing Open P6 precision bearing with low-torque cage and light oil lubricant (3x7x2)	01589400	S683-MC3-JP6-AF12	




WS75

Position / Description / Application	Manufacturers Part	MBA Part	
Replacement Press Button	04940500	BKU60500	
Bearing for Lower Shaft High Precision Bearing Double shielded P6 precision stainless steel bearing (3x10x4)	01577000	S623-ZZ-MC3-JP6	
Replacement Middle Gear Complete 20:1	06025100	BKU65100	
Head Gear Complete 20:1	04939900	BKU69900	


Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production


WS75 EKM

Position / Description / Application	Manufacturers Part	MBA Part	
Head Gear Complete 20:1	05101300	BKU61300	
Replacement Middle Gear Complete 20:1	04016900	BKU66900	
Bearing for Middle Gear High Precision Bearing Open Abec 3 precision bearing (2.35x5.5x2)	04733900	SMR55-0235-MC3-GR3-AF12	

WS75E KM

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Precision Bearing Open Abec 7 precision special dental bearing with low-torque ribbon retainer and light oil lubricant (4x7x2)	04027300	SMR74-MC34-GR7-R-L23	

WS92E

Position / Description / Application	Manufacturers Part	MBA Part	
Bearing for Headgear High Speed Ceramic-Hybrid Bearing - Torlon Retainer Open Abec 7 precision stainless ceramic-hybrid bearing with high speed Torlon retainer (4x7x2)	04025900	SMR74-CSN-TH7-TO	

Key Usually stocked or Short Delivery
** Available while stocks last

Long Delivery or Large Minimum Quantities
Not in current production

Friction Grip

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144-D28-ZS-KC-OPT

High Speed Bearing - Phenolic Retainer

Precision dental handpiece bearing with narrow design and high speed phenolic retainer (3.18x6.35x2.38)

2 Required



Push Button

Position / Description / Application

Manufacturers Part

MBA Part

Replacement Bearing

SR144A-D27-ZS-T9H-OPT

High Speed Bearing - Polyamide Retainer

Stainless steel bearing with polyamide retainer for precision dental handpieces (3.18x6.35x2.78)

2 Required



Key	Usually stocked or Short Delivery
	** Available while stocks last

Long Delivery or Large Minimum Quantities
<i>Not in current production</i>



Miniature Bearings Australia
<https://minibearings.com.au>
sales@minibearings.com.au

We sell an extensive range of products with varying amounts of technical data. Please assist us to provide you with accurate information by placing your orders and enquiries online or by e-mail.

Dental Handpiece Replacement Parts